

# anchor

## Bath Chair

Available in 2 sizes, the **anchor** bath chair is ideal for growing children with special needs who require assistance while bathing. The **anchor** is designed to provide a range of positioning for optimum fit and comfort in the shower and bath. The aluminum frame is covered with a fast-drying mesh fabric, padded support bars for positioning comfort and non-slip rubber pads for secure placement. The **anchor** folds flat for easy storage.

**Highly adjustable** – the **anchor** features one-handed angle adjustable legs, seat and back. Both the seat and back are angle adjustable at 3 different positions.

**Growth adjustability** – The **anchor** is seat depth adjustable and allows the chair to grow with the child.

### Standard features

Adjustable headrest

Trunk and pelvic positioning strap

Leg straps

*Headrest and positioning straps can be placed anywhere along the length of the chair and are secured when slipping it through the slots of the two vertical strips on the chair.*

### Accessories

Rolling shower base

Lateral support

