


# Rogue Backrest Mount and Hardware



Please note the following: <p>

**\*\*Backrest Bracket Mount Kits (parts 107199 and 107200) and the Backrest Transit Bracket Mount Kit (part 107245) contain complete assemblies for the left and right sides. <p>**

**\*\*This Rogue Backrest Mount and Hardware is effective starting with serial number RG0016612. <p>**

**\*\*If retrofitting to this backrest style, new back frames and backrest tubes are required. <p>**

**\*\*If replacing backrest parts for a Discontinued Rogue with XP option, see Discontinued Rogue Backrest Mount and Hardware page. If you have a Rogue XP, see the Rogue XP Backrest Mount pages.**

Item Number	Part Number	Description	Quantity	Retail
1-4, 5a, 5c, 6-8, 11-21	107199	Backrest Bracket Mount Kit, Black Anodized <i>This kit does not contain a release bar or release bar ends. Select appropriate Backrest Release Bar Kit in addition to this kit for a complete assembly.</i>	1	\$367.41
1-4, 5b, 5d, 6-8, 11-21	107200_ _ _	Backrest Bracket Mount Kit, Color Anodized <i>This kit does not contain a release bar or release bar ends. Select appropriate Backrest Release Bar Kit in addition to this kit for a complete assembly.</i>	1	\$398.31
1-4, 5e, 6-8, 11-21	107245	Backrest Transit Bracket Mount Kit, Black Anodized <i>This kit does not contain a release bar or release bar ends. Select appropriate Backrest Release Bar Kit in addition to this kit for a complete assembly.</i>	1	\$367.41
4, 20	108107	Backrest Pivot Replacement Kit w/Instruction (PAIR)	1	\$116.75

Item Number	Part Number	Description	Quantity	Retail
		<i>This kit contains items 4 (qty 2) and 20 (qty 2), and instructions.</i>		
9a, 10	107190	Backrest Release Bar Kit, 12" Wide, Black Anodized with Backrest Release Bar Ends	1	\$28.68
9b, 10	107191	Backrest Release Bar Kit, 13" Wide, Black Anodized with Backrest Release Bar Ends	1	\$28.68
9c, 10	107192	Backrest Release Bar Kit, 14" Wide, Black Anodized with Backrest Release Bar Ends	1	\$28.68
9d, 10	107193	Backrest Release Bar Kit, 15" Wide, Black Anodized with Backrest Release Bar Ends	1	\$28.68
9e, 10	107194	Backrest Release Bar Kit, 16" Wide, Black Anodized with Backrest Release Bar Ends	1	\$28.68
9f, 10	107195	Backrest Release Bar Kit, 17" Wide, Black Anodized with Backrest Release Bar Ends	1	\$28.68
9g, 10	107196	Backrest Release Bar Kit, 18" Wide, Black Anodized with Backrest Release Bar Ends	1	\$28.68
9h, 10	107197	Backrest Release Bar Kit, 19" Wide, Black Anodized with Backrest Release Bar Ends	1	\$28.68
9i, 10	107198	Backrest Release Bar Kit, 20" Wide, Black Anodized with Backrest Release Bar Ends	1	\$28.68
1	000055	Saddle 1"	4	\$1.58
2	001762	Arc Adjustment Cam	2	\$3.98
3a	003232	M6 x 55 BHCS Blk Zc <i>Used on chairs with 0" or positive center of gravity</i>	4	\$8.94
3b	003241	M6 x 60 BHCS P/THD Custom Blk Zc <i>Used on chairs with negative center of gravity</i>	4	\$11.56
4	003734	Pivot Assembly	2	\$58.37
5a	003711	Back Bracket Right, Black Anodized	1	\$56.80
5b	003711_ _ _	Color Anodized Back Bracket, Right	1	\$69.20
5c	003712	Back Bracket Left, Black Anodized	1	\$56.80
5d	003712_ _ _	Color Anodized Back Bracket, Left	1	\$69.20
5e	004434	Back Bracket with Transit, Black, Rigid	1	\$26.72
6	003713	Strike Hex	4	\$5.25
7	003715	Back Saddle 1 3/8" Tube	2	\$2.43
8	003716	Saddle 1 3/8" Tube	4	\$1.05
9a	003739	Backrest Release Bar 12" Wide, Black Anodized	1	\$21.01
9b	003740	Backrest Release Bar 13" Wide, Black Anodized	1	\$21.01
9c	003741	Backrest Release Bar 14" Wide, Black Anodized	1	\$21.01

Item Number	Part Number	Description	Quantity	Retail
9d	003742	Backrest Release Bar 15" Wide, Black Anodized	1	\$21.01
9e	003743	Backrest Release Bar 16" Wide, Black Anodized	1	\$21.01
9f	003744	Backrest Release Bar 17" Wide, Black Anodized	1	\$21.01
9g	003745	Backrest Release Bar 18" Wide, Black Anodized	1	\$21.01
9h	003746	Backrest Release Bar 19" Wide, Black Anodized	1	\$21.01
9i	003747	Backrest Release Bar 20" Wide, Black Anodized	1	\$21.01
10	003750	Backrest Release Bar End	2	\$2.92
11	103284	M6 x 35 Shoulder Bolt BZ	2	\$4.73
12	100658	M6 Nylock Nut Blk Zn	4	\$1.05
13	101593	M6 x 25 BTN SHCS BLKZC	2	\$1.58
14	101835	M6 Lock Washer Blk Zn	4	\$0.53
15	101853	M6 Washer Blk Zn	2	\$0.53
16	102614	M5 Split Lock Washer BLK ZC	2	\$0.53
17	102621	M5 x 8 BHCS BLK ZC	2	\$0.53
18	103022	M5 Washer, Blk Zn	2	\$0.53
19	103227	M6 x 30 Long Shoulder BHCS BZ	2	\$4.20
20	103279	#8-16 Tri-Lobular Plastic Screw BZ	2	\$0.53
21	103283	M5 Top Lock Nut BZ	2	\$0.53

# Color Anodizing Samples

