


# Rogue XP Backrest Mount - Growth



*Back Release Cable based on seat width The Back Release Cable length will measure 1" less than the specified seat width*

Item Number	Part Number	Description	Quantity
1a, 3	121180	Back Release Cable 8" Wide Black with Lock Out Pins	1
1b, 3	121179	Back Release Cable 9" Wide Black with Lock Out Pins	1
1c, 3	121181	Back Release Cable 10" Wide Black with Lock Out Pins	1
1d, 3	121182	Back Release Cable, 11" Wide, Black with Lock Out Pins	1

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
1e, 3	121183	Back Release Cable, 12" Wide, Black with Lock Out Pins	1
1f, 3	121184	Back Release Cable, 13" Wide, Black with Lock Out Pins	1
1g, 3	121185	Back Release Cable, 14" Wide, Black with Lock Out Pins	1
1h, 3	121186	Back Release Cable, 15" Wide, Black with Lock Out Pins	1
1i, 3	121187	Back Release Cable, 16" Wide, Black with Lock Out Pins	1
1j, 3	121188	BACK RELEASE CABLE, 17" WIDE, BLACK with Lock Out Pins	1
1k, 3	121189	BACK RELEASE CABLE, 18" WIDE, BLACK with Lock Out Pins	1
2	000965	Backrest Cable Release Spring RG	2
3	119660	Lockout Pin	2