



## Arc Front Side Frame - Growth



**Extra Short**

*Is your depth changing? Will your existing side frames achieve your new depth? If no, please select the appropriate side frame. Select color of parts based on Base Frame color of the chair. If quoting a Left Front Frame for Australia or New Zealand, 122921 is required when transit is not present.*

| Item Number | Part Number | Description                                                                    | Quantity |
|-------------|-------------|--------------------------------------------------------------------------------|----------|
| 1a, 2, 3    | 113109_ _ _ | Front Side Frame, Extra Short, Right, Arc<br><i>Accommodates 12-15" Depths</i> | 1        |
| 1b, 2, 3    | 113110_ _ _ | Front Side Frame, Extra Short, Left, Arc<br><i>Accommodates 12-15" Depths</i>  | 1        |
| 1c, 2, 3    | 113111_ _ _ | Front Side Frame, Short, Right, Arc<br><i>Accommodates 14-17" Depths</i>       | 1        |
| 1d, 2, 3    | 113112_ _ _ | Front Side Frame, Short, Left, Arc<br><i>Accommodates 14-17" Depths</i>        | 1        |
| 1e, 2, 3    | 113113_ _ _ | Front Side Frame, Medium, Right, Arc<br><i>Accommodates 16-19" Depths</i>      | 1        |

| Item Number | Part Number | Description                                                              | Quantity |
|-------------|-------------|--------------------------------------------------------------------------|----------|
| 1f, 2, 3    | 113114_ _ _ | Front Side Frame, Medium, Left, Arc<br><i>Accommodates 16-19" Depths</i> | 1        |
| 1g, 2, 3    | 113115_ _ _ | Front Side Frame, Long, Right, Arc<br><i>Accommodates 18-20" Depths</i>  | 1        |
| 1h, 2, 3    | 113116_ _ _ | Front Side Frame, Long, Left, Arc<br><i>Accommodates 18-20" Depths</i>   | 1        |
|             | 122921      | Not AS/NZS 3696.19 Compliant On-Chair Label with Instruction             | 1        |