

# Liberty Wheels - 12" Mag



*If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles and/or axle plates may be needed.*

Item Number	Part Number	Description	Quantity
1a	001424	12" Mag Wheel w/ Pneumatic Assembly Tilt	1
1b	001425	12" Mag Wheel w/ Pneumatic w/ Airless Insert Assembly Tilt	1
2	001426	12" Mag Wheel w/ Full Poly Assembly Tilt	1