


Focus MAXX Mag Wheel / Tire Assemblies



These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles and/or receivers may be needed. For the Catalyst 4, new axle plates may be needed. See Focus Axles page for axle and receiver utilization.

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| 1a | 116870 | 22" MAXX Mag Wheel with Pneumatic w/ Airless Insert Tire, Black | 1 |
| 1b | 116869 | 24" Wheel MAXX Mag 5-Spoke, Pneumatic w/ Airless Insert Tire, Black | 1 |