


# Focus MAXX Mag Wheel / Tire Assemblies



*These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles and/or receivers may be needed. For the Catalyst 4, new axle plates may be needed.*

Item Number	Part Number	Description	Quantity
1a	116870	22" MAXX Mag Wheel with Pneumatic w/ Airless Insert Tire, Black	1
1b	116869	24" MAXX Mag Wheel with Pneumatic w/ Airless Insert Tire, Black	1