


Rogue ALX Seat Upholstery - Growth (Formerly Tsunami)



SCREW (ITEM 3) QTY PER CHAIR	
SEAT DEPTH	SCREW QTY
14" - 15"	8
16" - 18"	10
19" - 21"	12
22"	14

Seat upholstery width measures approximately 1" less than the chair size Seat upholstery depth measures true to the chair size

Item Number	Part Number	Description	Quantity
1a	101000	Seat Upholstery Rigid/Liberty 1414	1
1b	101005	Seat Upholstery Rigid 1415	1

Item Number	Part Number	Description	Quantity
1c	101012	Seat Upholstery Rigid/Liberty 1416	1
1d	101019	Seat Upholstery Rigid 1417	1
1e	101026	Seat Upholstery Rigid/Liberty 1418	1
1f	101033	Seat Upholstery Rigid 1419	1
1g	101074	Seat Upholstery Rigid/Liberty 1420	1
1h	101001	Seat Upholstery Rigid/Liberty 1514	1
1i	101006	Seat Upholstery Rigid 1515	1
1j	101013	Seat Upholstery Rigid/Liberty 1516	1
1k	101020	Seat Upholstery Rigid 1517	1
1l	101027	Seat Upholstery Rigid/Liberty 1518	1
1m	101034	Seat Upholstery Rigid 1519	1
1n	101040	Seat Upholstery Rigid/Liberty 1520	1
1o	101002	Seat Upholstery Rigid/Liberty 1614	1
1p	101007	Seat Upholstery Rigid 1615	1
1q	101014	Seat Upholstery Rigid/Liberty 1616	1
1r	101021	Seat Upholstery Rigid 1617	1
1s	101028	Seat Upholstery Rigid/Liberty 1618	1
1t	101035	Seat Upholstery Rigid 1619	1
1u	101041	Seat Upholstery Rigid/Liberty 1620	1
1v	101542	Seat Upholstery Rigid 1621	1
1w	101543	Seat Upholstery Rigid 1622	1
1x	101003	Seat Upholstery Rigid/Liberty 1714	1
1y	101008	Seat Upholstery Rigid 1715	1
1z	101015	Seat Upholstery Rigid/Liberty 1716	1
1aa	101022	Seat Upholstery Rigid 1717	1
1ab	101029	Seat Upholstery Rigid/Liberty 1718	1
1ac	101036	Seat Upholstery Rigid 1719	1
1ad	101042	Seat Upholstery Rigid/Liberty 1720	1

Item Number	Part Number	Description	Quantity
1ae	101544	Seat Upholstery Rigid 1721	1
1af	101545	Seat Upholstery Rigid 1722	1
1ag	101004	Seat Upholstery Rigid/Liberty 1814	1
1ah	101009	Seat Upholstery Rigid 1815	1
1ai	101016	Seat Upholstery Rigid/Liberty 1816	1
1aj	101023	Seat Upholstery Rigid 1817	1
1ak	101030	Seat Upholstery Rigid/Liberty 1818	1
1al	101037	Seat Upholstery Rigid 1819	1
1am	101043	Seat Upholstery Rigid/Liberty 1820	1
1an	101546	Seat Upholstery Rigid 1821	1
1ao	101547	Seat Upholstery Rigid 1822	1
1ap	101010	Seat Upholstery Rigid 1915	1
1aq	101017	Seat Upholstery Rigid/Liberty 1916	1
1ar	101024	Seat Upholstery Rigid 1917	1
1as	101031	Seat Upholstery Rigid/Liberty 1918	1
1at	101038	Seat Upholstery Rigid 1919	1
1au	101044	Seat Upholstery Rigid/Liberty 1920	1
1av	101548	Seat Upholstery Rigid 1921	1
1aw	101549	Seat Upholstery Rigid 1922	1
1ax	101011	Seat Upholstery Rigid 2015	1
1ay	101018	Seat Upholstery Rigid/Liberty 2016	1
1az	101025	Seat Upholstery Rigid 2017	1
1ba	101032	Seat Upholstery Rigid/Liberty 2018	1
1bb	101039	Seat Upholstery Rigid 2019	1
1bc	101045	Seat Upholstery Rigid/Liberty 2020	1
1bd	101550	Seat Upholstery Rigid 2021	1
1be	101551	Seat Upholstery Rigid 2022	1
1bf	101528	Seat Upholstery Rigid/Liberty 2116	1

Item Number	Part Number	Description	Quantity
1bg	101529	Seat Upholstery Rigid 2117	1
1bh	101530	Seat Upholstery Rigid/Liberty 2118	1
1bi	101531	Seat Upholstery Rigid 2119	1
1bj	101532	Seat Upholstery Rigid/Liberty 2120	1
1bk	101533	Seat Upholstery Rigid 2121	1
1bl	101534	Seat Upholstery Rigid 2122	1
1bm	101535	Seat Upholstery Rigid/Liberty 2216	1
1bn	101536	Seat Upholstery Rigid 2217	1
1bo	101537	Seat Upholstery Rigid/Liberty 2218	1
1bp	101538	Seat Upholstery Rigid 2219	1
1bq	101539	Seat Upholstery Rigid/Liberty 2220	1
1br	101540	Seat Upholstery Rigid 2221	1
1bs	101541	Seat Upholstery Rigid 2222	1
2a	000241	Seat Rail 14" Deep Rigid/Liberty	2
2b	000242	Seat Rail 15" Deep RIGID	2
2c	000243	Seat Rail 16" Deep Rigid/Liberty	2
2d	000244	Seat Rail 17" Deep RIGID	2
2e	000245	Seat Rail 18" Deep Rigid/Liberty	2
2f	000246	Seat Rail 19" Deep RIGID	2
2g	000247	Seat Rail 20" Deep Rigid/Liberty	2
2h	000552	Seat Rail 21" Deep RIGID	2
2i	000553	Seat Rail 22" Deep RIGID	2
3	100734	M5 x 16 PHMS ZC	14