

anchor

bathing chair



718.943.6452



www.circlespecialty.com



info@circlespecialty.com



anchor

bathing chair

The Anchor provides a carefree bathing experience for children and caregivers alike. We make safety and comfort a priority so your child looks forward to bath time.

We designed the Anchor knowing your child's needs are unique and constantly evolving. That's why you can make angle and depth adjustments with one hand on either side of the chair, all while your child is seated. This agile and individualized functionality ensures longevity, comfort, and a hassle free bath time **so every bubble counts.**



Angle Adjustable

back, seat, & calf support



Folds Flat

easy storage & transport

FEATURES

- Depth adjustable for growth
- Angle adjustable back, seat and calf support, configure for individual requirements
- Adjustable for caregiver's height
- Positioning straps and head support can be placed anywhere along length of the chair
- Folds flat for easy storage and transport

ACCESSORIES

01 Rolling Shower Base

AN2100

Four swivel wheels with locking casters at ideal height for the caregiver, can be easily and securely attached to a depth adjustable frame. Accommodates both chair sizes at multiple depth and angle positions.



Product Dimension	47 L x 23.5 W x 15.5 H
Product Weight	9 lbs

02 Lateral Support

AN2110

Wedges and adjustable straps help maintain optimal trunk alignment.




03 Positioning Straps, XL

AN2111

Set of XL trunk & pelvic straps.



ADJUSTMENTS



SPECIFICATIONS

AN1400

Small

Seat Width	14.5"
Seat Depth	10.5" - 13"
Back Height	24"
Back Width	17"
Calf Support Length	11"
Seat Height	5.75"-11"-15.5"-17"
Back Angle	0°, 30°, 60°, 90°
Calf Support Angle	0°, 30°, 60°
Product Weight	14.5 lbs
Weight Capacity	155 lbs

AN1600

Medium

Seat Width	16"
Seat Depth	13"-16"
Back Height	32"
Back Width	18.5"
Calf Support Length	14"
Seat Height	5.75"-13"-18.25"-20.5"
Back Angle	0°, 30°, 60°, 90°
Calf Support Angle	0°, 30°, 60°
Product Weight	18 lbs.
Weight Capacity	155 lbs

NOTE:

We regretfully cannot accept returns for bathing and toileting products due to public hygiene concerns.