





With full adjustability and a minimalist mono-tube frame, the ZRA dishes out style without boundaries. Add the ZRA's TiFit titanium frame,

and you have the complete package. Strength, style, and the ability to set-up and configure your chair just for you.





TECHNICAL SPECIFICATIONS.

Frame Features

Frame Options

Transport Weight* Weight Limit

Total Chair Width* Seat Width Front Seat Width (Seat Taper) Seat Depth Ergonomic Seat Front Seat Height Rear Seat Height Camber Front Frame Angle Footrest Width Options Center of Gravity Adjustment Front Wheels Rear Wheels Handrims Convenience

Aesthetics

Titanium Mono-Tube Rigid TiFit Full Adjustability ZRA Black Edition SuperLite Edition Armed Forces Edition Heavy Duty Edition Power Adaptable Frame Reinforced Frame 10.7 lbs. 265 lbs. Standard 300 lbs. Heavy Duty Edition 23" 10" - 20" 3 Options 10" - 20" 7 Options 16" - 21" 12.5" - 20" 0°, 2°, 4°, 6°, 8°, 12°, 15° or 18° 70°, 75°, 80°, 85° or 90° Standard Taper or V-Front End 6" 23 Options 67 Options 15 Options Transit Tie-Down Option Folding Angle Adjustable Seat Back Depth Adjustable Back/Frame 6 Titanium Frame Finishes 19 Frame Colors 6 Color Anodize Packages 7 Ultrasuede Colors 8 Spinergy Wheel Colors Lifetime on Frame Limited Lifetime on Components

Warranty

*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lbs.



TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit[®] chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. The optimized fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!



TiLite | The Ultimate Ride™ 800.545.2266 | tilite.com | © 2014 TiLite, LLC. TiLite wheelchairs are shown with available options.