

Anti-Rollback System



The Anti-Rollback System is only available as a full kit. For replacement parts please contact Pratiko Medical. If retrofitting, please visit the Back Tubes for Anti-Rollback page for complete assembly. If retrofitting Anti-Rollback System to a Cat4, new axle brackets and quick release axles may be needed see below.

Item Number	Part Number	Description	Quantity
	122959	Anti-Rollback System for 22" Rear Wheels Cat 4	1
	122960	Anti-Rollback System for 24" Rear Wheels Cat 4	1
	122961	Anti-Rollback System for 22" Rear Wheels Cat 5/5Vx	1
	122962	Anti-Rollback System for 24" Rear Wheels Cat 5/5Vx	1
	114491	CAT 4 Axle Bracket, Right <i>Required for Cat 4 Retrofit unless Maxx Mag is present</i>	1
	114492	CAT 4 Axle Bracket, Left <i>Required for Cat 4 Retrofit Required for Cat 4 Retrofit unless Maxx Mag is present</i>	1
	101891	Axle Spacer .15" CAT4 <i>Used as needed</i>	1
	101892	Axle Spacer .35" CAT4 <i>Used as needed</i>	1
	101191	Axle Quick Release w/ Button 1/2" x 110MM <i>Required for Cat 4 Retrofit</i>	1

Catalyst Axle Reference

Quick Release Axles				
Wheel	Receivers			
	000052 / 100756	114491 / 114492	114480	000281
Assembled 6-Spoke Mag	100901	N/A	102371	102432
X Core Mag	101191	N/A	N/A	N/A
Ki Spoke / MAXX Spoke / Spoke / Spinergy	101108	101191	101191	101182
MAXX Mag	N/A	101191	101191	N/A
w/ Quad Quick Release Lever	100755	N/A	N/A	100901

Drum Brake Quick Release Axle	
MAXX Spoke	101191

Fixed Axle	
All Wheels except Assembled 6-Spoke Mag	100899
Assembled 6-Spoke Mag	102086

Axle Part Numbers		
Push Pin Style	Push Button Style	Threaded (Fixed)
100755 (1/2" x 99mm)	101108 (1/2" x 99mm)	100899 (1/2" x 4-3/4")
100901 (1/2" x 104mm)	101182 (1/2" x 104mm)	102086 (1/2-20 x 5-1/2")
102432 (1/2" x 112mm)	101191 (1/2" x 110mm)	
	102432 (1/2" x 112mm)	
	102371 (1/2" x 122mm)	