



ORDER FORM

To see QUICKIE Power Options or Accessories, scan here or visit: https://www.sunrisemedical.com/ power-wheelchair-options-accessories



Q500 with Captain's Seat

EIPW30

Customer	Ship To
Account #	Name
Date	Address
PO#	City
Buyer	State Zip
Mark For	Phone

Client Measurements

d) Elbow to seat

j) Actual sitting depth

m) Foot length


n) Knee to foot

o) Seat width

a) Top of head to seatb) Top of shoulder to seat

in

- c) Underarm to seat
- e) Top of leg to seat
- f) Shoulder width
- g) Chest widthi) Hip width



Notes

Chair C	pecifications
LIAUS	NECHICATIONS
Unun c	poornouciono

Mid-Wheel Drive Base Specifications

		Drive Wheel Position	Mid-wheel drive
HCPC	K0848 / K0849 / K0856 / K0857	User Weight Capacity	300 lbs.
		Base Size	Width 24", Length 33"
Weights starting at: (with Batteries)	NCC: 277 lbs SCC: 301 lbs	Minimum Turning Radius	18.5" (excludes footplates)
		Drive Control	R-net PM 90A / VR-2 90A
Power Seating	50° CG Tilt / 12" seat elevator	Driving Range	Up to 12 miles
		Drive Motors	4-pole 6 mph / 10 kph
Seat Type	Captain's Seat	Drive Wheels	13" w/flat-free inserts (black)
	w/ 25° Adjustable Recline	Casters	7" front / 7" rear (black)
Saat Ciza Panga	Width 16" to 20"	Charger	8A Dual Mode (Compact)
Seat Size Range	Depth 16" to 22"	Battery Size	Group 34 / Group 22
Seat-to-Floor Height Range	17" to 19"	Color	8 shroud colors and 7 color accents
	17 10 13	Ground Clearance	3"

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
Group 3 Stand	ard Duty	/ Mid Wheel Drive Base End Item Number EIPW30 (Captain's Seat)	\$
300PC218	K0849	Quickie Q500 M NCC No Power Captain Seat and Back Included in Base Pricing : Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, 6-wheel Suspension System, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 17".	\$8,541.00
300PC219	K0857	Quickie Q500 M SCC Single Power Option Captains Seat and Back Included in Base Pricing : Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, 6-wheel Suspension System, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 17".	\$8,541.00
300PC220	K0848	Quickie Q500 M NCC No Power Captains Seat Back and Solid Seat Included in Base Pricing : Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, 6-wheel Suspension System, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 17".	\$8,541.00
300PC221	K0856	Quickie Q500 M SCC Single Power Option Captains Seat Back and Solid Seat Included in Base Pricing : Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, 6-wheel Suspension System, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 17".	\$8,541.00
User Weight			\$
300WG23	_	Up to 300 lbs.	NC
Motor Packag			\$
300PF99	e	4-pole 10 kph / 6 mph (Standard for Group 3 K0848, K0849, K0856, K0857)	ہ NC
	_		
	K0100	SureTrac® Module & Motor Encoders (Improves Drive Tracking)	¢1.020.00
300MI01	K0108		\$1,628.00
Drive Wheels	_		\$
300DW69	-	Drive Wheel, 13", Dual Density, Knobby, Black	NC
300DW58	-	Drive Wheel, 14", Dual Density, Knobby, Black	NC
Caster Wheels	5		\$
300CW59	-	Caster Wheels 7", Solid, Black	NC
Shroud Color			\$
300C119	-	Matte Black	NC
300C165	-	Candy Apple Red	NC
300C166	-	Candy Blue	NC
300C167	-	Candy Purple	NC
300C176	-	Green Apple	NC
300C178	-	Pearl Pink	NC
300C182	-	Silver	NC
300C170	-	American Flag	\$602.00

BASE	_					
OPTION	HCPC	DESCRIPTION				RETAIL
Color Accents						\$
	ors and shi	roud colors may not match	l.			1
300EZ20	-	Apple Green				NC
300EZ22	-	Black				NC
300EZ21	-	Blue				NC
300EZ23	-	Orange				NC
300EZ24	-	Red				NC
300EZ25	-	Silver				NC
300EZ19	-	White				NC
SO Transit Tie	Downs					\$
300T050	-	ISO Transit Package				NC
.ED Lights						\$
300LI04	-	LED Lights and Indicator	S			\$1,400.00
300LI08	_	LED Fender Lights (Sing	le Switch Activation)			\$566.00
300LI02	_	Omit LED Lights				NC
Seat Heights						\$
						μ Ψ
		Seat Height	Power Seating with Tilt	Power Seating with Elevate	No Power Seating	
300HT55	-	17"	N/A	N/A	Group 34	NC
300HT104	-	17.5"	N/A	N/A	Group 34	NC
300HT56	-	18"	Group 34	N/A	Group 34	NC
300HT105	-	18.5"	Group 34	Group 34	Group 34	NC
300HT57	-	19"	Group 34	Group 34	Group 34	NC
Batteries and	Chargers	3				\$
iroup 34 Batter	y Options					
300AG27	E2359	GP34 Flat Top Pair w/St	raps MK Gel			\$1,097.00
300AG29	E2359	GP34 Flat Top Pair w/St	raps AGM			\$828.00
Group 22 Batter	у					
300AG31	E2361	GP22 Flat Top Pair w/St	raps			\$613.00
)mit Batteries						
300AG24	-	Battery Harnesses Only	(Omit Batteries)			NC
Chargers			- ()			
300BJ30	-	10A Dual Mode (Compa				\$110.00
300BJ19	-	8A Dual Mode (Compac	L)			NC

Q500 M with Captain's Seat

SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
Power Seating	Options		\$
Seat Elevator			
300P107	E2298	Seat Elevator, Powered, 300 mm / 12"	\$3,791.00
Seat Tilt			
300PA48	E1002	Posterior 50° Tilt	\$8,458.00

Pre Set Seat Fi	rame An	gle (No Power Seating Options, Increases Seat to Floor Height)	\$
300FA10	_	0° (Adds 0" to STFH)	NC
300FA3	-	3° (Adds 1" to STFH)	NC

Backrest Frame	е Туре		\$
300BF08	-	Captain Seat Frame (Includes Back Cushion)	NC

Seat Cushion T	уре				\$
300ST77	-	Captain's Cus	aptain's Cushion (See Section A)		
300ST80	-	JAY® (See Se	ection B)		NC
300ST75	-	Omit			NC
A - Seat Size					\$
A - Seat Size Seat Pan Width			Available Depths		\$
	-	16"	Available Depths 16" to 18"		\$ NC
Seat Pan Width	-	16" 18"			

Seat Pan Depth (Front of Seat Pan to Back Structure. Back Cushions Can Cause 2" Loss of Effective Depth)

300RF51	-	16"	NC
300RF45	-	17"	NC
300RF46	-	18"	NC
300RF47	-	19"	NC
300RF48	-	20"	NC
300RF49	-	21"	NC
300RF50	-	22"	NC

SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
B - JAY [®] Seat	Cushion		¢
D - JAT - Seat	GUSHIOH		φ
300CU59	E2607	JAY® lon® (Sizes 16" - 20" Wide)	\$460.00
300CU61	E2607	JAY® Union® (Sizes 16" - 20" Wide)	\$520.00
300CU64	E2622	JAY® Fusion® (Sizes 16" - 20" Wide)	\$514.00
300CU104	E2622	JAY® Fusion® with Cryo [™] Technology (Sizes 16" - 20" Wide)	\$630.00

OPTION	HCPC	DESCRIPTION	RETAIL
rmrest Struct	ture		
300EG37	-	Flip Back (Includes Waterfall Arm Pads) (Height Adjustable 11" to 14")	N
ixed Centerm	ount Leg	irest	
300H183	-	SEDEO Fixed Centermount Legrest (3" to 17.5" Leg Length)	\$292.0
300H184	-	Fixed Centermount Legrest, Short (Includes Calf Pads and Single Footplate) (11.5" to 14.5" Leg Length)	Ν
300H185	-	Fixed Centermount Legrest, Long (Includes Calf Pads and Single Footplate) (13.5" to 18.5" Leg Length)	Ν
ower Centerr	nount Le	grest	
300H175	E1012	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension (Mechanical Articulation) (5" to 17.5" Leg Length)	\$2,706.
300H109	E1012	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension (Powered Articulation) (3" to 17.5" Leg Length)	\$3,791.
entermount L	egrest A	Accessories	
ootplates	1/0040	Devel 11-1-1-1 MC-10	ф107
300F127	K0040	Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D)	\$197.
300F128	K0040	Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D)	\$197.
ootplate Acces			****
300DI68	K0108	Extra Large Centermount Footplate (Single Footplate Overlay 16"W x 12"D)	\$199.
300D169	K0108	Footplate Connector Link (Attaches footplates together)	\$97.
300DI57	K0108	Gel Footpads (Pair, LH and RH)	\$403.
alf Pad Option 300D134	s _	Calf Pads, Pair, Standard Integral Foam	1
300DI35	_	Calf Pads, Pair, Padded (N/A w/ Std footplate)	\$408.
elvic Position			•••••
300A837		Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	1
300A839	_	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	I
300A840	E0978	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$88.
300A845	K0108	Bodypoint – Evoflex Pelvic Stabilizer Belt	\$203.
300A846	K0108	Bodypoint – Evoflex Pelvic Stabilizer Belt with Extension	\$224.
ccessories			
2 Tank Option	1		
300A718	E2208	02 Tank Holder Mounted to Backrest (D and E Styles)	\$379.
aptain's Acce	ssories		
300A608	E2207	Cane / Crutch Holder	\$303.
300A610	K0108	Walker Holder	\$303.
lisc. Accessor		Deer Manuel Diverse	¢100
3008A52	-	Ram Mount Phone	\$100. \$409
3008A53	-	Ram Mount Tablet	\$408. \$100.
3008A54 3009A10	-	Ram Mount Cup Holder Utility Hooks	\$100. \$60.
ledial Thigh S	-		Φ00.
3008A80	E0957	Medial Thigh Support	\$197
000A00	L0337	inional migh oupport	φ13/

Q500 M with Captain's Seat

JOYSTICK / E	LECTRO	DNICS						
OPTION	HCPC	DESCRIPTION	RETAIL					
VR-2 Joystick			\$					
300J127	_	VR-2 Joystick	NC					
B-net Joystick	R-net Joystick <i>(Remote)</i> \$							
300J121	_	Standard LED Joystick (w/o Screen Display)	NC					
300J116	_	Color Joystick (2.3" Color LCD Screen Display, CJSM1)	NC					
300J118	_	High Resolution Color Joystick, (CJSM2+BT) 3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control	\$1,230.00					
300J119	-	Omit Remote Control	NC					
Joystick Moun	iting		\$					
300JM34	-	Fixed on Armrest	NC					
300JM35	E1028	Parallel Swing-Away	\$345.00					
300JM36	E1028	Link-It Mid Line Mount on Armrest (Swing-Away)	\$991.00					
Joystick Moun	t (Locati	nn)	\$					
300HM16	_	Right Hand Armrest	NC					
300HM17	_	Left Hand Armrest	NC					
Display Langua	age Sele		\$					
300L1	-	English	NC					
300L3	-	French	NC					
300L7	-	Spanish	NC					
Joystick Hand	le (Optio	nal)	\$					
300JH13	E2323	Ball Handle	\$147.00					
300JH15	E2323	T-Handle	\$147.00					
300JH24	E2323	Bodypoint – Rubber Dome Joystick Handle	\$191.00					
300JH25	E2323	Bodypoint – U-Shaped Joystick Handle (4 in)	\$191.00					
300JH26	E2323	Bodypoint – Grooved Mushroom Joystick Handle	\$191.00					
Electronic Prog	grammer	S	\$					
300PR13	-	R-net PC Wireless Programmer	\$613.00					
R-net Power N	lodules		\$					
300ER12	-	Power Module, VR-2 90A, Non-Expandable	NC					
300ER22	_	Power Module, R-net, 90A, Non-Expandable (Only available with single power)	NC					
300ER13	-	Power Module, R-net, 90A, Expandable	\$1,335.00					
300EX11	E2377	Expandable 90A Controller - R-net	\$706.00					
300EX12	E2313	Harness for 90A Expandable Controller	\$629.00					

8

RETAIL

JOYSTICK / ELECTRONICS

OPTION HCPC DESCRIPTION

0.1101			11217112		
Seating Contro	ol		\$		
Powered Seatin	ng Control				
300EO40	-	Omit Thru-Drive (Requires Switch Control)	NC		
300E037	E2310	Thru-drive Control for 1 Actuator	\$2,240.00		
300E038	E2311	Thru-drive Control for 2+ Actuators	\$4,055.00		
Switch Control (Available with Omit Thru-Drive only) (1 actuator for Tilt or Elevate only)					
300E033	-	Dual Toggle for 1 Actuator	NC		
300EO34	-	2.5" Single Button Switch for 1 Actuator (Button switches are not mounted)	NC		
300EO41	-	Ctrl+5 Switch Box for 2 Actuators or More	NC		
Additional Pow	er Seatin	g Control (Available with Thru-Drive only)			
300QG11	-	Ctrl+5 Switch Box with Free Assignable Functions	\$1,459.00		
Dual Toggle Mo	unting Lo	cation			
300TN1	-	Toggle Mount, Right	NC		
300TN2	-	Toggle Mount, Left	NC		
Ctrl+5 Mounting	g Locatior	n (Will swing away if mounted on the same side as the joystick)			
300EY19	-	Right Armrest	NC		
300EY20	-	Left Armrest	NC		
Ctrl+5 Switch T	ype				
300EY22	-	Push Buttons	NC		
300EY23	K0108	Toggle Switches	\$350.00		
R-net Program	mina		\$		
	3		¥		

Assignable Buttons

Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

	300PM12	-	Include Assignable Buttons	NC
Power Centermount Footplate Extension Option				
	300IE03	K0108	Include 7.5" Powered Footplate Extension for Transfer Assist (Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only Available with Powered Articulation	\$1,459.00



Q500 M with Captain's Seat	\$			
Base Price				
Total of additional equipment				
Packaging and transport				
Total Price				



800.333.4000

Sunrise Medical (US) LLC Mount Juliet,, TN 37122 01.2025 MK-129645 Rev. P