


Arc Axles



| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 1 | 101191 | Axle Quick Release w/ Button 1/2" x 110MM <i>Used on 20" Mags wheels</i> | 1 |
| 2 | 101108 | Axle Quick Release w/ Button 1/2" x 99MM <i>Used on 12" and 16" Mags, 20" MAXX Performance Spoke, and 20" MAXX Mag wheels</i> | 1 |
| 3 | 110293 | HHB, 1/2-13 x 4.5, Grade 8, BLZ | 1 |