


# Canada Catalyst 4 Transit



Notice: Chair configuration review required prior to parts being quoted or ordered. <p>

Only available in Canada on Catalyst 4. Not available on Catalyst 4C <p>

Not available with 250 lb or 350 lb weight capacity <p>

Requires Adjustable Axle Plate and Adjustable Axle Plate Rear Frames

Item Number	Part Number	Description	Quantity	Retail
3a	101826	M6 x 50 HHCS Blk Zn <i>Use quantity 2 if only Transit is present; Use quantity 1 if Transit and Anti-Tip bracket are present</i>	4	\$2.21
3b	100676	M6 x 55 CLS10.9 DIN931 HHCS ZC <i>Use quantity 2 if Transit and Anti-Tip bracket are present. Use quantity 1 if Swing Away Arm bracket, Transit, and Anti-Tip bracket are present</i>	4	\$1.66
3c	100677	M6 x 60 CLS10.9 DIN931 HHCS ZC	4	\$1.66

Item Number	Part Number	Description	Quantity	Retail
		<i>Used to attach Front Transit Bracket. Use quantity 2 if Transit, Anti-Tip bracket and Swing Away Arm bracket are present</i>		
4	000056	Saddle 1 1/8"	10	\$1.66
5	100746	M6 Flat Washer Blk Zn	4	\$1.66
6	100558	M6 NYL INS Domed Cap L/N ZC	4	\$1.10
7	100658	M6 Nylock Nut Blk Zn	6	\$1.10
8	100510	Eccentric Washer	4	\$1.73