


# Catalyst 5 Wheels - SpinerGY LX



*Recommended to order the SpinerGY LX Replacement Spoke Tool Kit (110119) when replacing spokes. Replacement SpinerGY Spokes and SpinerGY Wheels with Flex Rims will have an extended lead time. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed.*

Item Number	Part Number	Description	Quantity
1a	100886	22" x 1/2" SpinerGY LX Black Spoke	1
1b	100939	22" x 1/2" SpinerGY LX Red Spoke	1
1c	100940	22" x 1/2" SpinerGY LX White Spoke	1

Item Number	Part Number	Description	Quantity
1d	100938	22" x 1/2" Spinergy LX Blue Spoke	1
1j	100942	24" x 1/2" Spinergy LX Black Spoke	1
1k	100941	24" x 1/2" Spinergy LX Red Spoke	1
1l	100944	24" x 1/2" Spinergy LX White Spoke	1
1m	100943	24" x 1/2" Spinergy LX Blue Spoke	1
1s	118348	24" x 1/2" Spinergy LX Black Spoke w/ Pneumatic w/ Airless Insert Tire (Black)	1
1t	118349	24" x 1/2" Spinergy LX Red Spoke w/ Pneumatic w/ Airless Insert Tire (Black)	1
1u	118350	24" x 1/2" Spinergy LX White Spoke w/ Pneumatic w/ Airless Insert Tire (Black)	1
1v	118351	24" x 1/2" Spinergy LX Blue Spoke w/ Pneumatic w/ Airless Insert Tire (Black)	1
1ab	100994	25" x 1/2" Spinergy LX Black Spoke	1
1ac	100996	25" x 1/2" Spinergy LX Red Spoke	1
1ad	100998	25" x 1/2" Spinergy LX White Spoke	1
1ae	100995	25" x 1/2" Spinergy LX Blue Spoke	1
1ak	100921	26" x 1/2" Spinergy LX Black Spoke	1
1al	100946	26" x 1/2" Spinergy LX Red Spoke	1
1am	100945	26" x 1/2" Spinergy LX White Spoke	1
1an	100947	26" x 1/2" Spinergy LX Blue Spoke	1
2	102854	Spinergy LX Spoke Wrench <i>No charge when ordered with Spinergy wheel or replacement spoke</i>	1
	110119	Spinergy LX Replacement Spoke Tool Kit <i>This tool kit contains: Spinergy LX Spoke Wrench and Spinergy Spoke T-Wrench</i>	1
3a	101929	Replacement Spoke Spinergy LX 22", Black (511291-8.75BLA)	1
3b	101930	Replacement Spoke Spinergy LX 22", Red (511291-8.75RED)	1
3c	101931	Replacement Spoke Spinergy LX 22", White (511291-8.75WHT)	1
3d	101932	Replacement Spoke Spinergy LX 22", Blue (511291-8.75BLU)	1
3i	101937	Replacement Spoke Spinergy LX 24", Black (511291-9.5BLA)	1
3j	101938	Replacement Spoke Spinergy LX 24", Red (511291-9.5RED)	1

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
3k	101939	Replacement Spoke Spinergy LX 24", White (511291-9.5WHT)	1
3l	101940	Replacement Spoke Spinergy LX 24", Blue (511291-9.5BLU)	1
3q	101945	Replacement Spoke Spinergy LX 25", Black (511291-9.84BLA)	1
3r	101946	Replacement Spoke Spinergy LX 25", Red (511291-9.84RED)	1
3s	101947	Replacement Spoke Spinergy LX 25", White (511291-9.84WHT)	1
3t	101948	Replacement Spoke Spinergy LX 25", Blue (511291-9.84BLU)	1
3y	101953	Replacement Spoke Spinergy LX 26", Black (511291-10.47BLA)	1
3z	101954	Replacement Spoke Spinergy LX 26", Red (511291-10.47RED)	1
3aa	101955	Replacement Spoke Spinergy LX 26", White (511291-10.47WHT)	1
3ab	101956	Replacement Spoke Spinergy LX 26", Blue (511291-10.47BLU)	1

### Tire Component Chart

Wheel		1 3/8" Street Pneumatic			1" High Pressure			1 3/8" Euro Trek Knobby			1" Schwalbe				1 3/8" Schwalbe HS-127			2" Mountain		
Wheel	Wheel Part #	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Mara-thon Tire	One Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip
22" Spoke	116730	116735	100588	101908	102961	100609	101908											101475	100588	101908
24" Spoke	116728	116736	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596
18" MAXX Spoke	109285	101587	101588	101137																
20" MAXX Spoke	108244	116843	100912	101137	102957	102844	101137													
22" MAXX Spoke	105135	116735	100588	101908	102961	100609	101908											101475	100588	101908
24" MAXX Spoke	105136	116736	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596
25" MAXX Spoke	107436				102959	101457	101452				101357	103177	101457	101452				100166	100208	101452
26" MAXX Spoke	107437	100636	100590	100608	102960	100582	100597				100865		100582	100597						
20" MAXX Mag	111853	116843	100912	101137																
22" MAXX Mag	111854	116735	100588	101908														101475	100588	101908
24" MAXX Mag	111855	116736	100589	100596				102982	100589	100596					102920	100589	100596	101476	101188	100596
24" Superlight	101159	111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596			
25" Superlight	101158				102959	101457	101452				101357	103177	101457	101452						
20" Spinergy Spox	See Spinergy Spox Page	116843	100912	101137	102957	102844	101137													
22" Spinergy Spox		116735	100588	101908	102961	100609	101908											101475	100588	101908
24" Spinergy Spox		116736	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596
25" Spinergy Spox					102959	101457	101452				101357	103177	101457	101452				100166	100208	10145
26" Spinergy Spox		100636	100590	100608	102960	100582	100597				100865		100582	100597						
22" Spingery LX	See Spinergy LX Page	116735	100588	101908	102961	100609	101908											101475	100588	101908
24" Spingery LX		116736	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596
25" Spingery LX					102959	101457	101452				101357	103177	101457	101452				100166	100208	101452
26" Spingery LX		100636	100590	100608	102960	100582	100597				100865		100582	100597						

See Page 2 for Discontinued Wheels.

**Tire Component Chart (Discontinued Wheels)**

Wheel		1 3/8" Street Pneumatic			1" High Pressure			1 3/8" Euro Trek Knobby			1" Schwalbe				1 3/8" Schwalbe HS-127			2" Mountain			
Wheel	Wheel Part #	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Mara-thon Tire	One Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	
18" Ki Spoke	200529	101587	101588	101137																	
20" Ki Spoke	200530	100911	100912	101137	102957	102844	101137														
22" Ki Spoke	200531	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Ki Spoke	200532	111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Ki Spoke	200533				102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" Ki Spoke	200534	100636	100590	100608	102960	100582	100597				100865		100582	100597							
20" 5-Spoke X Core	101961	100911	100912	101137																	
22" 5-Spoke X Core	101962	100634	100588	101908																	
24" 5-Spoke X Core	100960	111818	100589	100596				102982	100589	100596					102920	100589	100596				

### Tire Compatibility Chart

Wheel		Tire											
Wheel	Wheel Part #	1 3/8" Street Pneumatic	1 3/8" Full Profile Poly	1" Low Profile Poly	1 3/8" Pneumatic w/ Airless Insert (Tire / Insert)	1" Solid Sentinel	1" Solid Shox	1" High Pressure	1 3/8" Euro Trek Knobby	1" Schwalbe Marathon	1" Schwalbe One	1 3/8" Schwalbe HS-127	2" Mountain
22" Spoke	116730	116735	100621	100632	116735 / 100503	103165	102718	102961					101475
24" Spoke	116728	116736	100622	100633	116736 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
18" MAXX Spoke	109285	101587	101992	100918	101587 / 101092								
20" MAXX Spoke	108244	116843	100917	102085	116843 / 100913	102974	102717	102957					
22" MAXX Spoke	105135	116735	100621	100632	116735 / 100503	103165	102718	102961					101475
24" MAXX Spoke	105136	116736	100622	100633	116736 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
25" MAXX Spoke	107436					102977	101453	102959		101357	103177		100166
26" MAXX Spoke	107437	100636			100636 / 100505	102978	100884	102960		100865			
12" Mag <sup>1</sup>	N/A	001424	001426		001425								
16" Mag <sup>1</sup>	N/A	001159	001161	003453	001160								
20" MAXX Mag	111853	116843	100917	102085	116843 / 100913	102974							
22" MAXX Mag	111584	116735	100621	100632	116735 / 100503	103165							101475
24" MAXX Mag	111855	116736	100622	100633	116736 / 100504	102976			102982			102920	101476
24" Superlight	101159	100635					100894	102958	102982	100864	103176	102920	
25" Superlight	101158						101146	102959		101357	103177		
20" Spinerger Spox	See Spinerger Spox Page	116843	100917	102085	116843 / 100913	102974		102957					
22" Spinerger Spox		116735	100621	100632	116735 / 100503	103165		102961					101475
24" Spinerger Spox		116736	100622	100633	116736 / 100504	102976	100894	102958	102982	100864	103176	102920	101476
25" Spinerger Spox						102977	101146	102959		101357	103177		100166
26" Spinerger Spox		100636				100636 / 100505	102978	100895	102960		100865		
22" Spinerger LX	See Spinerger LX Page	116735	100621	100632	116735 / 100503			102961					101475
24" Spinerger LX		116736	100622	100633	116736 / 100504	102976	100894	102958	102982	100864	103176	102920	101476
25" Spinerger LX						102977	101146	102959		101357	103177		100166
26" Spinerger LX		100636				100636 / 100505	102978	100895	102960		100865		

<sup>1</sup> 12" and 16" Mag Wheels are sold as a combination Wheel and Tire Assembly

See Page 2 for Discontinued Wheels

**Tire Compatibility Chart (Discontinued Wheels)**

Wheel		Tire											
Wheel	Wheel Part #	1 3/8" Street Pneumatic	1 3/8" Full Profile Poly	1" Low Profile Poly	1 3/8" Pneumatic w/ Airless Insert (Tire / Insert)	1" Solid Sentinel	1" Solid Shox	1" High Pressure	1 3/8" Euro Trek Knobby	1" Schwalbe Marathon	1" Schwalbe One	1 3/8" Schwalbe HS-127	2" Mountain
18" Ki Spoke	200529	101587	101992	100918	101587 / 101092								
20" Ki Spoke	200530	100911	100917	102085	100911 / 100913	102974	102717	102957					
22" Ki Spoke	200531	100634	100621	100632	100634 / 100503	103165	102718	102961					101475
24" Ki Spoke	200532	111818	100622	100633	111818 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
25" Ki Spoke	200533					102977	101453	102959		101357	103177		100166
26" Ki Spoke	200534	100636			100636 / 100505	102978	100884	102960		100865			
20" 5-Spoke X Core	101961	100911	100917	102085	100911 / 100913		102717						
22" 5-Spoke X Core	101962	100634	100621	100632	100634 / 100503		102718						
24" 5-Spoke X Core	100960	111818	100622	100633	111818 / 100504		100640		102982			102920	

Effective 2/27/2023

### Handrim Hardware Chart

Wheel		Handrim Connection				Handrim										
Wheel	Wheel Part No.	Connection Points	Rivet/Tab	Spacer*	Screw*	Aluminum Anodized	Superlight	Plastic Coated	Projection	Ergonomic Standard	Ergonomic LT	Natural Fit Standard	Natural Fit LT	Flex Rim		
22" Spoke	116730	5 or 6	116732 (Used with 5 Rivet) 100698 (Used with 6 Rivet)	100653 (Not used w/ Superlight Handrim)	Aluminum, Plastic Ctd, Projection Ergonomic: 100654 Natural Fit: 100835	101898 (5 Rivet)		101964 (5 Rivet)	112824 (6 Rivet)	113082 (6 Rivet)	113085 (6 Rivet)	200538 (6 Rivet)	200201 (6 Rivet)			
24" Spoke	116728	5 or 6	Not used w/ Superlight Handrims		Superlight: 100669 Nut: 100657	100975 (5 Rivet)	101161 (6 Tab)	100976 (5 Rivet)	112825 (6 Rivet)	113083 (6 Rivet)	113806 (6 Rivet)	100793 (6 Rivet)	200202 (6 Rivet)			
18" MAXX Spoke	109285	3	100698 (Not used w/ Superlight Handrim)	100653 (Not used w/ Superlight Handrim)	Aluminum, Plastic Ctd, Projection, Ergonomic: 100654 Natural Fit: 100835 Superlight: 100669 Nut: 100657	100206										
20" MAXX Spoke	108244	6				200536		112819								
22" MAXX Spoke	105135					100560		112820	112824	113082	113085	200538	200201			
24" MAXX Spoke	105136					200349	101161	112821	112825	113083	113086	100793	200202			
25" MAXX Spoke	107436					200350		112822	200548	113084	113087	200539	200540			
26" MAXX Spoke	107437					200351		112823	200549			100907	101454			
20" MAXX Mag	111853	6		Aluminum, Plastic Coated, Projection: 100629	Aluminum, Plastic Coated, Projection, Ergonomic: 103545	200536		112819								
22" MAXX Mag	111854			Natural Fit & Ergonomic: 101756	Natural Fit: 101893	100560		112820	112824	113082	113085	200538	200201			
24" MAXX Mag	111855					105137		112821	112825	113083	113086	100793	200202			
24" Superlight	101159	6			100536	100754	101161	100836		113080	113081	100830	100828			
25" Superlight	101160							101160	101091			101464	101460			

\* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 100792 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels, part number 100666 for the screw on 5-Spoke X Core wheels, or part number 100654 for the screw on MAXX Mag wheels.

Effective 11/27/2023



Handrim Hardware Chart (Continued)															
Wheel		Handrim Connection				Handrim									
Wheel	Wheel Part No.	Connection Points	Rivet/Tab	Spacer*	Screw*	Aluminum Anodized	Superlight	Plastic Coated	Projection	Ergonomic Standard	Ergonomic LT	Natural Fit Standard	Natural Fit LT	Flex Rim	
20" Spinerger Spox	See Spinerger Spox Page	6			Ergonomic: 100536	103125		103179							
22" Spinerger Spox						100827		100808				100889	100888		
24" Spinerger Spox						All Others: 100669	100766	101161	100615				100830	100828	
25" Spinerger Spox						Nut: 100657	100767	101160	100765				101464	101460	
26" Spinerger Spox							101477		101148				200200	100950	
22" Spinerger LX	See Spinerger LX Page	6			Ergonomic: 100536	100827		100808				100889	100888		
24" Spinerger LX						100766	101161	100615		113080	113081	100830	100828	See Spinerger LX Page	
25" Spinerger LX						All Others: 100669	100767	101160	100765				101464	101460	
26" Spinerger LX						Nut: 100657	101477		101148				200200	100950	

Handrim Hardware Chart (Discontinued Wheels)													
Wheel		Handrim Connection				Handrim							
Wheel	Wheel Part No.	Connection Points	Rivet/Tab	Spacer*	Screw*	Aluminum Anodized	Superlight	Plastic Coated	Projection	Natural Fit Standard	Natural Fit LT	Flex Rim	
18" Ki Spoke	200529	3	100698	100653	Aluminum, Plastic Ctd, Projection: 100654  Natural Fit: 100835	100206		101106					
20" Ki Spoke	200530	6				200536		200542					
22" Ki Spoke	200531					100560		100576	100569	200538	200201		
24" Ki Spoke	200532					200349		100577	200547	100793	200202		
25" Ki Spoke	200533					200350		101870	200548	200539	200540		
26" Ki Spoke	200534					200351		100578	200549	100907	101454		
20" 5-Spoke X Core	101961	5		Aluminum, Projection, Natural Fit: 100629	100724	101897		101963					
22" 5-Spoke X Core	101962						101898		101964				
24" 5-Spoke X Core	100960						100975		100976	200546			
24" 5-Spoke X Core	100960				Plastic Ctd: 101756	101893				100768	100769		

\* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 100792 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels, part number 100666 for the screw on 5-Spoke X Core wheels, or part number 100654 for the screw on MAXX Mag wheels.

Effective 11/27/2023