


## Focus Wheels - MAXX Spoke for Drum Brake



*If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed.*

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
1a, 2a, 3a	111733	Wheel 24" MAXX Spoke, Drum Brake	1
2a	112197	Replacement Spoke, Drum Brake Wheel, 24" Cross	1
3a	112196	Replacement Spoke, Drum Brake Wheel, 24" Radial	1