


Liberty Tilt Mechanism - Growth

CABLE LENGTHS			
BACK POST TYPE	SEAT DEPTH RANGE		
	SHORT (14-16D)	MEDIUM (16-18D)	LONG (18-20D)
HEIGHT ADJ SHORT	003288	003289	003290
HEIGHT ADJ MEDIUM	003289	003290	003291
HEIGHT ADJ TALL	003290	003291	003292
STROLLER HANDLE SHORT	003289	003290	003291
STROLLER HANDLE TALL	003291	003292	003293



Did your depth change? If yes, evaluate if a different size cable is necessary

Item Number	Part Number	Description	Quantity
1-3, 4a	506801	Cable/Trigger Kit with Short Cable, Liberty	1
1-3, 4b	506802	Cable/Trigger Kit with Medium Cable, Liberty	1
1-3, 4c	506803	Cable/Trigger Kit with Long Cable, Liberty	1
1-3, 4d	506804	Cable/Trigger Kit with XL Cable, Liberty	1
1-3, 4e	506805	Cable/Trigger Kit with XXL Cable, Liberty	1
1-3, 4f	506806	Cable/Trigger Kit with XXXL Cable, Liberty	1
4a	003288	Cable, Tilt, Short	2

Item Number	Part Number	Description	Quantity
4b	003289	Cable, Tilt, Medium	2
4c	003290	Cable, Tilt, Long	2
4d	003291	Cable, Tilt, XL	2
4e	003292	Cable, Tilt, XXL	2
4f	003293	Cable, Tilt, XXXL	2