


Titanium 70 Degree Extension Mount Hanger



Please Note: If chair was built with a fixed front frame or hangers were omitted, then latch blocks need to be ordered with hangers. Please reference the side frame parts manual page.

| Item Number | Part Number | Description | Quantity | Retail |
|-----------------------|-------------|---|----------|----------|
| 1 - 17 | 100381 | Hanger Assembly Titanium 70 Extension Mount Right | 1 | \$235.84 |
| 1 - 17 | 100382 | Hanger Assembly Titanium 70 Extension Mount Left | 1 | \$235.84 |
| 1-7, 9, 10, 12, 14-20 | 504741 | Hanger Assm Titanium 70 Extension Mount 4-Way - In-Line | 1 | \$235.84 |
| 1-7, 9, 10, 12, 14-20 | 504742 | Hanger Assm Titanium 70 Extension Mount 4-Way - Out Left | 1 | \$235.84 |
| 1-7, 9, 10, 12, 14-20 | 504743 | Hanger Assm Titanium 70 Extension Mount 4-Way - Out Right | 1 | \$235.84 |
| 1, 7, 14 | 501019 | Hanger Tube Assembly Titanium 70 Extension Mount | 1 | \$180.64 |
| 8 - 13 | 003417 | Classic Hanger Latch Assembly, Left | 1 | \$63.78 |
| 8 - 13 | 003416 | Classic Hanger Latch Assembly, Right | 1 | \$63.78 |

| Item Number | Part Number | Description | Quantity | Retail |
|------------------|-------------|--|----------|---------|
| 9, 10, 12, 14-20 | 003225 | Assembled 4-Way In Line Hanger Release | 1 | \$36.77 |
| 9, 10, 12, 14-20 | 003225INST | Assembled 4-Way In Line Hanger Release with Instruction | 1 | \$36.77 |
| 9, 10, 12, 14-20 | 003643 | Assembled 4-Way Out Hanger Release, Left | 1 | \$36.77 |
| 9, 10, 12, 14-20 | 003643INST | Assembled 4-Way Out Hanger Release, Left with Instruction | 1 | \$36.77 |
| 9, 10, 12, 14-20 | 003642 | Assembled 4-Way Out Hanger Release, Right | 1 | \$36.77 |
| 9, 10, 12, 14-20 | 003642INST | Assembled 4-Way Out Hanger Release, Right with Instruction | 1 | \$36.77 |
| 2 | 000073 | M5 x 24 Threaded Barrel | 1 | \$3.78 |
| 3 | 100669 | M5 x 10 BTN SHCS BLKZC w/Patch | 2 | \$1.58 |
| 4 | 100673 | M8 x 45 BTN SHCS BLKZC | 1 | \$1.58 |
| 5 | 000032 | Footrest Hanger Pivot Saddle | 1 | \$5.25 |
| 6 | 100659 | M8 Nylock Nut Blk Zn | 1 | \$1.05 |
| 7 | 100804 | Tube End Plug 1" | 1 | \$1.58 |
| 8 | 001158 | Lever Swing Away Short | 1 | \$7.31 |
| 9 | 100657 | M5 Nylock Nut Blk Zn | 1 | \$0.53 |
| 10 | 001374 | M5 x 40 BHCS P/THD BLKZC | 1 | \$2.14 |
| 11 | 000031 | Hanger Bracket | 1 | \$23.73 |
| 12 | 000119 | Hanger Spring | 1 | \$3.05 |
| 14 | 000034 | Sleeve 3/4" x 1 1/2" | 1 | \$3.66 |
| 15 | 000988 | Fastener Retainer | 1 | \$3.94 |
| 16 | 000987 | Threaded Fastener | 1 | \$3.78 |
| 17 | 200427 | M5 x 10 CLS8.8 HHCS ZC w/Patch DIN933 | 1 | \$0.53 |
| 18 | 003228 | Hanger Latch, 4-Way, Lever | 1 | \$3.15 |
| 18 | 003228INST | Hanger Latch, 4-Way, Lever with Instruction | 1 | \$3.15 |
| 19 | 003226 | Hanger Latch, 4-Way, Bracket | 1 | \$5.25 |
| 19 | 003226INST | Hanger Latch, 4-Way, Bracket with Instruction | 1 | \$5.25 |
| 20 | 003227 | Hanger Latch, 4-Way, Pin | 1 | \$19.96 |
| 20 | 003227INST | Hanger Latch, 4-Way, Pin with Instruction | 1 | \$19.96 |