

922



Available in  
Black or White

### Rehab Shower/Commode Chair-Folding

- Packages Now Available
- Stainless Steel Frame
- Several New Seat Options Available
- Easy to Fold and Transport

The unique design of the 922, Rehab Shower/Commode Chair-Folding, was developed to address the wide range of necessities for both user and caregiver.

The padded, ergonomically correct configuration of the seat emphasizes pressure relief and proper sitting form while the Stainless Steel frame provides a safe, non-corrosive, easily maneuverable chair that folds compactly for storage or travel.



**FRAMES**



**P13559**  
18Wx18"D  
Stainless Steel  
Frame

**WHEELS**



**P61095**  
22" Wheels  
(pair)



**P61096**  
24" Wheels  
(pair)



**P60549**  
24" Pneumatic  
Street Tread  
Wheels  
(pair)



**P61059**  
24" Flat Free Street  
Tread Wheels  
(pair)



**P61103**  
22" PVC Handrims  
(replaces standard  
handrims, pair)



**P61104**  
24" PVC Handrims  
(replaces standard  
handrims, pair)

**AXLES**



**P13639**  
Bolted Axles  
(pair)



**P13640**  
Quick Release  
Axles  
(pair)

**WHEEL LOCKS**



**P60532**  
Wheel Locks  
(pair)



**P13135**  
Wheel Locks with  
Extended Handle  
(pair)

**CASTERS**



**60841-2**  
5" Casters  
(pair)



**60800-2**  
5" Dual Locking  
Casters  
(pair)

**SEATS**



**PP60993**  
18"Wx18"D  
Comfortuff 1.5" waterfall  
foam  
Front/Rear Open Seat  
A-5"xB-8" opening



**PP60974**  
18"Wx18"D  
Comfortuff sloped foam  
Front Open Seat  
A-4"xB-6.5" opening



**P61120**  
18"Wx18"D  
Ensolite® 2" dual density  
foam w/coccyx notch  
Front Open Seat  
A-3"xB-5.25" opening



**PP61058**  
18"Wx18"D  
Ensolite® 2" waterfall foam  
Front Open Contoured Seat  
A-5.5"xB-8.25" opening



**PP61086**  
18"Wx18"D  
Ensolite® 2" foam  
Solid Seat



**PP60675**  
18"Wx18"D  
Ensolite® 1" foam  
Cloverleaf Seat  
A-7.5"xB-9" oval opening



**PP60701**  
18"Wx18"D  
Ensolite® 1" foam  
Front/Rear Open Seat  
A-4.55"xB-7.75" opening




**PP60850**  
18"Wx18"D  
Ensolite® 1" waterfall foam  
Front/Rear Open Seat  
A5"xB-8" opening




**PP60913**  
18"Wx18"D  
Ensolite® 2" waterfall foam  
Front/Rear Open Seat  
A5.25"xB-8" opening

**Front/Rear Open Seat  
Opening Measurements**



**Cloverleaf Seat Opening  
Measurements**



**ARMRESTS**



**P13563**  
Left Hand  
Height-Adjustable  
Armrest



**P13564**  
Right Hand  
Height-Adjustable  
Armrest



**P13565**  
Left Hand Flip-Up  
Armrest



**P13566**  
Right Hand Flip-Up  
Armrest

**LEG RESTS**



**PL13569**  
Left Gooseneck  
Leg Rest



**PR13569**  
Right Gooseneck  
Leg Rest



**P13794**  
Left Elevating  
Leg Rest



**P13795**  
Right Elevating  
Leg Rest



**PL13570**  
Left Low-Mount  
Leg Rest



**PR13570**  
Right Low-Mount  
Leg Rest

**FOOT PLATES**



**PL11099**  
Left Foot Plate



**PR11099**  
Right Foot Plate



**PL11100**  
Left Foot Plate with  
Heel Loop



**PR11100**  
Right Foot Plate  
with Heel Loop



**P61137**  
Left Oversized  
Aluminum  
Foot Plate



**P61138**  
Right Oversized  
Aluminum  
Foot Plate



**PL11101**  
Left Padded  
Foot Plate

**BACKS**



**PR11101**  
Right Padded  
Foot Plate



**P11098**  
18"W Solid Back



**P60915**  
18"W Sling Back

**BELTS**



**PW60909**  
Waist Belt  
(48" circumference with  
12" velcro overlap)



**PT60909**  
Trunk Belt  
(48" circumference with  
12" velcro overlap)



**P60879**  
Calf Strap



**P11093**  
Bodypoint®  
Trunk Belt



**P11102**  
Bodypoint®  
Trunk Belt (pair)



**P11111**  
Bodypoint® Trunk  
Belt & Waist Belt



**P11112**  
Bodypoint® XL  
Trunk Belt



**P11113**  
Belt w/Quick  
Release Knobs



**P61081**  
Bodypoint®  
Calf Strap



**P11114**  
Bodypoint®  
Knee Belt

**POSITIONING AIDS**



**PL13574**  
Left Hand  
Swing-Away  
Adjustable Lateral



**PR13574**  
Right Hand  
Swing-Away  
Adjustable Lateral



**P13573**  
18"Wx18"W  
Pan/Hanger

**PANS**

full option details at  
[activeaid.com](http://activeaid.com)



922 with Options



922 with Options

## SPECIFICATIONS

<b>Weight Capacity</b>	<b>350 lbs.</b>
<b>Floor to top of seat height</b>	<b>22.25"</b>
<b>Seat depth</b>	<b>18"</b>
<b>Seat width between arms</b>	<b>18"</b>
<b>Back height from seat</b>	<b>15.5"</b>
<b>Width under seat at wheel lug location</b>	<b>14.25"</b>
<b>Overall width</b>	<b>25"</b>
<b>Overall length</b>	<b>44.5" - 45.5"</b>
<b>Arm height from seat</b>	<b>8.5"</b>
<b>Footrest adjustment</b>	<b>16.5" - 21.25"</b>
<b>Anti-tipper floor clearance</b>	<b>1"</b>
<b>Underseat clearance from floor</b>	<b>18"</b>