



# KidWalk / KidWalk II

## Dynamic Mobility System

Patent Protected 6,832,770 B1



**Order Form  
Price List  
Effective 9-1-18**

Client: \_\_\_\_\_ Age: \_\_\_\_\_  
 Height: \_\_\_\_\_ Weight: \_\_\_\_\_ Dx: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, St, Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Facility: \_\_\_\_\_  
 Therapist: \_\_\_\_\_ Phone: \_\_\_\_\_

PO: \_\_\_\_\_ Date: \_\_\_\_\_  
 Dealer: \_\_\_\_\_  
 Contact: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, St, Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 E-mail: \_\_\_\_\_

### KidWalk / KidWalk II Dynamic Mobility System

**Standard Features Include:** Auto centering pelvic guides; Lateral weight shift; Vertical weight relief with adjustable load compensating spring assembly; Angle adjustable mast; Anti-tippers; Wheel locks; Quick release wheels; Padded seat; Adjustable axle position; Light up rear casters; Two piece frame for easy transport.

**KidWalk**  
**\$3,195.00**  
Dynamic Mobility System  
 Seat height 12" - 22"  
 Client Weight up to 65 pounds

Specifications:  
 Base: 24" wide by 27" long  
 Overall length: 33" - 35"  
 Mast Height: 33.5"  
 Upholstery: Black nylon  
 Wheels: 20" Mag  
 Casters: 4" Lighted



**KidWalk II**  
**\$3,795.00**  
Dynamic Mobility System  
 Seat height 15" - 28"  
 Client Weight up to 90 pounds

Specifications:  
 Base: 27.5" wide by 32" long  
 Overall length: 36" - 38"  
 Mast Height: 37.25"  
 Upholstery: Black nylon  
 Wheels: 24" Mag  
 Casters: 5" Lighted

**Frame Color:**

- Candy Blue  Candy Red  Candy Purple  Neon Pink  Neon Orange  Shocker Green  Razzberry  White  Black

**KidWalk Hip Support Configuration Selection: \***

- Standard Configuration:  Width 5" - 9"  Depth 4" x 7.5"  
 Optional Configurations:  Width 7" - 12"  Depth 5" - 7"  
 Flip Up Pad Brackets:  Straight (std)  10 Degree Angle  
 \*If no selections are marked the standard configuration will be supplied.

**KidWalk II Hip Support Configuration Selection: \***

- Standard Configuration:  Width 7"-12"  Depth 5" x 7"  
 Optional Width Configurations:  5.5" - 9"  8.5" - 15"  
 Optional Depth Configurations:  3.5" - 6" (with non-extended hip pads)  
 Flip Up Pad Brackets:  Straight (std)  10 Degree Angle

**Upper Body Support Options:**

- Lateral Chest System: Includes lateral side supports, swingaway anterior support pads and lumbar support pad .....\$ 495.00
- |   |  |
|---|--|
| <p><b>KidWalk</b></p> <p><input type="checkbox"/> Std 6.5" - 10.5" W x 4" - 6.5" D</p> <p><input type="checkbox"/> Optional Width Adj. 8" - 14.5"</p> <p><input type="checkbox"/> Optional Depth Adj. 5.5" - 8.5"</p> | <p><b>KidWalk II</b></p> <p><input type="checkbox"/> Std 8.5" - 15.5" W x 5.5" - 8.5" D</p> <p><input type="checkbox"/> Optional Width Adj. 6.5" - 10.5"</p> <p><input type="checkbox"/> Optional Depth Adj. 4" - 6.5"</p> |
|---|--|
- Flip Up Chest Pad Brackets:  10 Degree Angle (std)  Straight
- Upper Body Soft Support Wrap: For the client that needs minimal upper body support .....\$ 239.00
- Headrest Options:**
- Headpad: 3" x 4" flat headpad .....\$ 225.00
- Contoured Small Headrest: 4"H x 6"W .....\$ 329.00
- Contoured Medium Headrest: 5"H x 8"W .....\$ 395.00
- 1/2 Square Tube Mount: (Headrest not included) .....\$ 149.00
- 1/2 Flip Down Headrest Mount: .....\$ 339.00
- i2i Flip Down Headrest Mount: .....\$ 339.00
- Lower Extremity Positioning Options:**
- Extended Hip Pads: (Standard on 5" - 7" depth) .....\$ 95.00
- Dynamic Leg Wraps: Assist in supplying feedback to the client while encouraging movement in the lower extremities .....\$ 149.00

**Seat Options:**

- Wide Seat: 3" x 9": Large padded seat .....\$ 195.00
- Extended Seats: Longer sides assist in controlling adduction and scissoring: Specify length of sides below .....\$ 299.00
- Short 4" length
- Medium 6" length
- Long 8" length
- Leg Straps: Restricts movement of the lower extremities (Requires 6" or 8" extended seat).....\$ 195.00
- Additional Standard Seat: (Ordered with unit).....\$ 100.00

**Frame Options:**

- Cruiser Wheels: Allows use as an omni directional walker for clients up to 40" tall (Not available on **KidWalkII**) .....\$ 349.00
- Non-Lighted Caster Wheels: .....\$ no charge
- Sliding Ankle Prompts: Straps hold ankles to help keep the client's feet from crossing .....\$ 289.00
- Directional locks: Locks out caster swivel.....\$ 265.00

**Wheel Options:**

- Anti-Rollback Device:.....\$ 299.00
- 20" Spoke Wheels (**KidWalk** only):.....\$ 425.00
- 24" Spoke Wheels (**KidWalk II** only):.....\$ 425.00
- Spoke Protectors 20":.....\$ 219.00
- Spoke Protectors 24":.....\$ 219.00

**PRIME ENGINEERING**

A Division of Axiom Industries Inc.  
 All specifications subject to change without notice

4202 W. SIERRA MADRE AVE., FRESNO, CA 93722 1-800-827-8263 Fax 800-800-3355  
 E-mail: info@primeengineering.com Web Site: www.primeengineering.com

090118KW\_OF\_PL

