

Item Number	Part Number	Description	Quantity
		<i>Does not include attaching hardware. Quote 502171 if hardware is needed.</i>	
11g	107571	One Arm Drive - Blade Drive Shaft 17.67"	1
11i	108072	One Arm Drive - Blade Drive Shaft 19.67"	1
12	107574	One Arm Drive - Blade Style Axle 5.03"	1
13	100659	M8 Nylock Nut Blk Zn	4
14	101858	Washer Flat M8 Blk Zn	4
17a	101831	M8 x 80 HHCS Blk Zn	4
17b	101832	M8 x 90 HHCS Blk Zn	4
17c	102384	M8 x 100 HHCS Blk Zc Part THD	2
18	101843	M5 x 20 SHCS Blk Zn	8

Focus Blade OAD Drive Shaft Configurations Guide

Focus Seat Width	Part No.	Compressed Drive Shaft Length	Extended Drive Shaft Length
14"	107565	10.12"	11.72"
15"	107566	11.12"	12.72"
16"	107567	12.22"	13.82"
17"	107568	13.07"	14.67"
18"	107569	14.07"	15.67"
19"	107570	15.07"	16.67"
20"	107571	16.07"	17.67"
21"	107572	17.07"	18.67"
22"	108072	18.07"	19.67"