



Quickie Shark S

EISP6, Effective July 2009

ACCOUNT INFORMATION

Account No: _____ Drop Ship/Ship To: _____
 Date: _____ Name: _____
 P.O. #: _____ Address: _____
 Buyer: _____ City: _____ State: _____
 Marked For: _____ Zip Code: _____ Phone Number: _____

Quickie Shark S

■ 157BX9 Quickie Shark S \$4,899.00

FRAMES

Frame Color		
<input type="checkbox"/>	159C91	Amethyst Fantasy \$95.00
<input type="checkbox"/>	159C55	Black NC
<input type="checkbox"/>	159C1	Blue NC
<input type="checkbox"/>	159C13	Candy Teal NC
<input type="checkbox"/>	159C75	Flip Flop Blue / Green \$95.00
<input type="checkbox"/>	159C61	Forrest Green NC
<input type="checkbox"/>	159C85	Graphite Metallic NC
<input type="checkbox"/>	159C92	Illusion Lapis \$95.00
<input type="checkbox"/>	159C93	Mystery \$95.00
<input type="checkbox"/>	159C84	Platinum Metallic NC
<input type="checkbox"/>	159C54	Red Ruby NC
<input type="checkbox"/>	159C90	Silver NC
<input type="checkbox"/>	159C16	Silver Texture NC
<input type="checkbox"/>	159C53	Sun Yellow NC
<input type="checkbox"/>	159C94	Tutti Frutti Orange NC
<input type="checkbox"/>	159C95	White RAL NC
<input type="checkbox"/>	159C34	Omit Color Frame NC

Brilliant Colors:

<input type="checkbox"/>	159C73	Brilliant Lime NC
<input type="checkbox"/>	159C72	Brilliant Marine NC
<input type="checkbox"/>	159C52	Brilliant Purple NC
<input type="checkbox"/>	159C74	Brilliant Red NC
<input type="checkbox"/>	159C76	Brilliant Silver NC

Matte Colors:

<input type="checkbox"/>	159C68	Matte Lime NC
<input type="checkbox"/>	159C96	Matte Candy Teal NC
<input type="checkbox"/>	159C67	Matte Marine NC
<input type="checkbox"/>	159C65	Matte Purple NC
<input type="checkbox"/>	159C63	Matte Red NC
<input type="checkbox"/>	159C62	Matte Silver NC

Frame Length

<input type="checkbox"/>	157FL3	Frame Length Short NC
<input type="checkbox"/>	159FL18	Frame Length Long NC

Frame Length Notes:

Short Frame

Overall with:

20" Wheels 77.56" (197 cm)

26" Wheels 79.53" (202 cm)

without Rear Wheels 55.91" (142 cm)

without steering fork & wheels 41.34" (105 cm)

Transfer Seat Height 8.66" (22 cm)

Long Frame

Overall with:

20" Wheels 79.92" (203 cm)

26" Wheels 81.5" (207 cm)

without Rear Wheels 58.66" (149 cm)

without steering fork & wheels 43.31" (110 cm)

Transfer Seat Height 7.09" (18 cm)

Fork Style

<input type="checkbox"/>	157FK03	Fixed Fork	NC
<input type="checkbox"/>	157FK01	Detachable Fork	\$250.00

Fork Mount


<input type="checkbox"/>	159FJ01	High Mount w/ ball bearing	NC
<input type="checkbox"/>	159FJ02	Low Mount w/fixed bottom bracket height	NC

Position Fork Mount¹


<input type="checkbox"/>	157PZ01	15" Setting	NC
<input type="checkbox"/>	157PZ02	15.75" Setting (Std setting)	NC
<input type="checkbox"/>	157PZ03	16.5" Setting	NC

Fork Notes:

High Mount Fork Preset
@ 15.75" (40 cm)



Low Mount Fork Fixed
height @ 11.48" (30cm)
Not adjustable



SEATING

Seat Width

<input type="checkbox"/>	159S5	14" Seat Width (36 cm)	NC
<input type="checkbox"/>	159S6	15" Seat Width (38 cm)	NC
<input type="checkbox"/>	159S20	15.5" Seat Width (40 cm)	NC
<input type="checkbox"/>	159S21	16.5" Seat Width (42 cm)	NC
<input type="checkbox"/>	159S8	17" Seat Width (44 cm)	NC
<input type="checkbox"/>	159S9	18" Seat Width (46 cm)	NC

BACKREST

Back Type

<input type="checkbox"/>	159BT10	Standard Backrest (58 cm)	NC
<input type="checkbox"/>	159BT1	Folding Backrest ²	\$140.00
<input type="checkbox"/>	159BT58	Short Backrest (16 cm)	NC
<input type="checkbox"/>	159BT85	Medium Backrest ⁴ (30cm)	NC

Backrest Notes:

Standard Backrest (shown on right)

22.83" (58 cm) long

Height and Angle Adjustable
(30 to 90 degree)

Fold Down Backrest



22.83" (58 cm) long

Uses Quick Release Pin and

Retainer Ring for ease of adjustment

Short Backrest (shown on right)

6.3" (16 cm) long



FOOTNOTES:


¹ Not available with low mount fork.

² Not available on short backrest with low fork


³ Only available with Short or Medium Backrest

⁴ Not available in folding

Backrest Position		
<input type="checkbox"/>	159BQ07	Adj (10"-21.5") Preset @ 12" ¹ NC
<input type="checkbox"/>	159BQ08	Adj (14"-26") Preset @ 16" ² NC

<p>Backrest Position Preset @ 12" (30.48 cm)</p> 


Neck Support		
<input type="checkbox"/>	159NT01	Neck Support ³ \$125.00

<p>Neck Support Diagram: Height and Depth Adjustable</p> 

FOOTREST		
Footrest		
<input type="checkbox"/>	159F236	Adj Short (16"-31.5") (40-80 cm) NC
<input type="checkbox"/>	159F237	Adj Long (27.5"-37.5") (70-95 cm) NC

<p>Footrest Notes: Adjustable Short (16"-31.5") Adjustable Long (27.5"-37.5")</p>
--

CONTROL OPTIONS

<p>Control Notes: Standard V Crank 6.89" (17.5 cm) long 19.29" (49 cm) wide (center to center of Handle) Custom Crank (shown to right)</p>  <p>Handle Vertical Small and Std @ 15 Degree Diagonal Brake Parking Brake Std Brake lever can be locked when brake is on</p>
--

Hand Crank		
<input type="checkbox"/>	159CK2	Crank Custom ¹¹ NC
<input type="checkbox"/>	159CK3	Omit Crank -\$150.00

Hand Mount		
<input type="checkbox"/>	159HM5	Frame Mounted Right Side ⁵ NC
<input type="checkbox"/>	159HM6	Grips Mounted on Right Handle NC
<input type="checkbox"/>	159HM7	Steering fork mounted on right handle NC

Hand Crank Width⁶		
<input type="checkbox"/>	159HJ01	15.74" (40 cm) NC
<input type="checkbox"/>	159HJ02	16.53" (42 cm) NC
<input type="checkbox"/>	159HJ03	17.71" (45 cm) NC
<input type="checkbox"/>	159HJ04	18.50" (47 cm) NC
<input type="checkbox"/>	159HJ05	19.29" (49 cm) NC
<input type="checkbox"/>	159HJ06	20.07" (51 cm) NC
<input type="checkbox"/>	159HJ07	21.25" (54 cm) NC
<input type="checkbox"/>	159HJ08	22.44" (57 cm) NC

Hand Crank Length⁶		
<input type="checkbox"/>	159HK01	6.88" (17.5 cm) ¹³ NC
<input type="checkbox"/>	159HK02	7.28" (18.5 cm) ¹⁰ NC
<input type="checkbox"/>	159HK03	7.67" (19.5 cm) NC
<input type="checkbox"/>	159HK04	8.07" (20.5 cm) NC
<input type="checkbox"/>	159HK05	8.46" (21.5 cm) NC

Gear Type (Choose 1)		
<input type="checkbox"/>	159GT02	Shimano Ultegra ⁹ NC

Chain Plate Kit		
<input type="checkbox"/>	159FZ01	4 Arm Chain Plate (26/36/48) ⁷ NC
<input type="checkbox"/>	159FZ02	5 Arm Chain Plate (30/42/52) ⁷ \$60.00

WHEELS

Drive Wheel		
<input type="checkbox"/>	159DW34	26 X 1" Hi Perf NC
<input type="checkbox"/>	159DW35	Omit Drive Wheel -\$50.00

Rear Wheel Size		
<input type="checkbox"/>	159WS4	20" Wheel NC
<input type="checkbox"/>	159WS18	26" Wheel Dt Rader NC
<input type="checkbox"/>	159WS20	Omit Rear Wheels ⁸ -\$55.00

Rear Wheel Tire		
<input type="checkbox"/>		1.35" Marathon Tire NC
<input type="checkbox"/>	159RT20	1" Hi Performance Reifen NC

Axle		
<input type="checkbox"/>	159AX17	Threaded Axle NC
<input type="checkbox"/>	159AX01	Quick Release Axle NC

ACCESSORIES

Accessories Sports		
<input type="checkbox"/>	159AZ01	Front Fender (26") \$50.00
<input type="checkbox"/>	159AZ02	Rear Fender (20") \$75.00
<input type="checkbox"/>		Rear Fender (26") \$75.00
<input type="checkbox"/>	159AZ03	Front Headlight \$50.00
<input type="checkbox"/>	159AZ04	Rear Light \$50.00
<input type="checkbox"/>	159AZ05	Cushion W/Composite Foam \$125.00
<input type="checkbox"/>	159AZ06	Waist Belt \$125.00
<input type="checkbox"/>	159AZ07	Signal Flag NC
<input type="checkbox"/>	159AZ08	Rear Bumper 20" Rear Wheels \$150.00
<input type="checkbox"/>	159AZ09	Rear Bumper 26" Rear Wheels \$150.00

HANDCYCLE TOTAL _____

Handling Charge¹²		
<input type="checkbox"/>		Standard Handling NC
<input type="checkbox"/>		Expedited Handling \$895.00

Footnotes: **TOTAL** _____

- ¹ Only available with High Mount Fork
- ² Only available with Low Mount Fork
- ³ Not available with Short Backrest option and 36 cm Seat Width.
- ⁴ Not available with Low Mount Fork.
- ⁵ Only available with Custom Crank option and fixed fork style.
- ⁶ Select option if Custom Crank option is selected.
- ⁷ Only available with Hand Crank Custom
- ⁸ Must select Axle Type.
- ⁹ Only available with Steering Fork right side.
- ¹⁰ Only available with Hand Crank width of 42, 45, 47, or 49 cm.
- ¹¹ Not available with 40 mm vertical handle.
- ¹² Standard Handling is a combination of sea & ground freight
Expedited Handling is a combination of air & ground freight
Please contact Customer Service for manufacturing lead times.
- ¹³ NA with hand crank width 51,54,57.