

Little Wave Handrims - Natural Fit



For more information see Handrim Hardware Chart. Quantity in Handrim Hardware Kits are for one wheel.<p>

Tab mounted handrims must be installed prior to the tire being installed.

Item Number	Part Number	Description	Quantity	Retail
1a	200538	Handrim Natural Fit Standard 6 Rivet Nut 22" <i>Used on 22" Ki Spoke and 22" MAXX Spoke</i>	1	\$213.41
1a	107168	Handrim Natural Fit Standard 6 Rivet Nut 22" w/ Micro Grip <i>Used on 22" Ki Spoke and 22" MAXX Spoke</i>	1	\$236.43
1b	100793	Handrim Natural Fit Standard 6 Rivet Nut 24" <i>Used on 24" Ki Spoke and 24" MAXX Spoke</i>	1	\$213.41
1b	107169	Handrim Natural Fit Standard 6 Rivet Nut 24" w/ Micro Grip <i>Used on 24" Ki Spoke and 24" MAXX Spoke</i>	1	\$251.11
1c	100768	Handrim Natural Fit Standard 5 Rivet Nut 24" <i>Used on 24" 5-Spoke Mag</i>	1	\$213.41
1c	107172	Handrim Natural Fit Standard 5 Rivet Nut 24" w/ Micro Grip <i>Used on 24" 5-Spoke Mag</i>	1	\$239.15
1d	100889	Handrim Natural Fit Standard 6 Tab 22" <i>Used on 22" Spinergy</i>	1	\$213.41
1d	107173	Handrim Natural Fit Standard 6 Tab 22" w/ Micro Grip <i>Used on 22" Spinergy</i>	1	\$236.43
1e	100830	Handrim Natural Fit Standard 6 Tab 24" <i>Used on 24" Spinergy</i>	1	\$213.41
1e	107174	Handrim Natural Fit Standard 6 Tab 24" w/ Micro Grip <i>Used on 24" Spinergy</i>	1	\$239.15

Item Number	Part Number	Description	Quantity	Retail
2a	200201	Handrim Natural Fit LT 6 Rivet Nut 22" <i>Used on 22" Ki Spoke and 22" MAXX Spoke</i>	1	\$234.83
2a	107177	Handrim Natural Fit LT 6 Rivet Nut 22" w/ Micro Grip <i>Used on 22" Ki Spoke and 22" MAXX Spoke</i>	1	\$256.83
2b	200202	Handrim Natural Fit LT 6 Rivet Nut 24" <i>Used on 24" Ki Spoke and 24" MAXX Spoke</i>	1	\$234.83
2b	107178	Handrim Natural Fit LT 6 Rivet Nut 24" w/ Micro Grip <i>Used on 24" Ki Spoke and 24" MAXX Spoke</i>	1	\$260.64
2c	100769	Handrim Natural Fit LT 5 Rivet Nut 24" <i>Used on 24" 5-Spoke Mag</i>	1	\$234.83
2c	107181	Handrim Natural Fit LT 5 Rivet Nut 24" w/ Micro Grip <i>Used on 24" 5-Spoke Mag</i>	1	\$259.55
2d	100888	Handrim Natural Fit LT 6 Tab 22" <i>Used on 22" Spinergy</i>	1	\$234.83
2d	107182	Handrim Natural Fit LT 6 Tab 22" w/ Micro Grip <i>Used on 22" Spinergy</i>	1	\$269.67
2e	100828	Handrim Natural Fit LT 6 Tab 24" <i>Used on 24" Spinergy</i>	1	\$234.83
2e	107183	Handrim Natural Fit LT 6 Tab 24" w/ Micro Grip <i>Used on 24" Spinergy</i>	1	\$272.53
3a	200204	Natural Fit Thumb Grip Regular 22"	1	\$144.34
3b	100795	Natural Fit Thumb Grip Regular 24"	1	\$144.34
4a	100887	Natural Fit Thumb Grip Super 22"	1	\$189.55
4b	100829	Natural Fit Thumb Grip Super 24"	1	\$189.55
5-7	502169	Handrim Hardware Kit, 6 Rivet Nut, Natural Fit <i>Used with Ki Spoke and MAXX Spoke Wheels</i>	1	\$38.61
8-9	502172	Handrim Hardware Kit, 5 Rivet Nut, Natural Fit <i>Used with 5-Spoke X Core Wheels</i>	1	\$33.10
10-11	502174	Handrim Hardware Kit, 6 Tab, EACH <i>Used with Spinergy</i>	1	\$16.55
5	100698	Handrim Tab <i>Used with Ki Spoke and MAXX Spoke wheels. Tabs are included / installed with MAXX Spoke wheels.</i>	6	\$4.17
6	100835	10/32 x 1 3/4" Pan Head Screw <i>Used with Ki Spoke and MAXX Spoke Wheels</i>	6	\$1.66

Item Number	Part Number	Description	Quantity	Retail
7	100653	Handrim Spacer Spoke <i>Used with Ki Spoke and MAXX Spoke Wheels</i>	6	\$1.10
8	100629	Handrim Spacer Alum Ano Riv <i>Used with 5-Spoke X Core Wheels</i>	5	\$1.16
9	101893	10-32 x 2 PHTH ZN <i>Used with 5-Spoke X Core Wheels</i>	5	\$4.41
10	100669	M5 x 10 BTN SHCS BLKZC w/Patch <i>Used with Spinergy</i>	6	\$1.66
11	100657	M5 Nylock Nut Blk Zn <i>Used with Spinergy</i>	6	\$0.56

Handrim Hardware Chart

Wheel		Handrim Connection				Handrim							
Wheel	Wheel Part No.	Connection Points	Tab	Spacer*	Screw*	Aluminum Anodized	Superlight	Plastic Coated	Projection	Natural Fit Standard	Natural Fit LT	Flex Rim	
18" Ki Spoke	200529	3	100698	100653	Aluminum, Plastic Ctd, Projection: 100654 Natural Fit: 100835	100206		101106					
20" Ki Spoke	200530					200536		200542					
22" Ki Spoke	200531					100560		100576	100569	200538	200201		
24" Ki Spoke	200532					200349	101161	100577	200547	100793	200202		
25" Ki Spoke	200533					200350		101870	200548	200539	200540		
26" Ki Spoke	200534					200351		100578	200549	100907	101454		
24" Superlight	101159	6			100536	100754	101161	100836		100830	100828		
25" Superlight	101160						101160	101091		101464	101460		
20" 5-Spoke X Core	101961	5		Aluminum, Projection, Natural Fit: 100629 Plastic Ctd: 101756	100724	101897		101963					
22" 5-Spoke X Core	101962						101898		101964				
24" 5-Spoke X Core	100960						100975		100976	200546			
24" 5-Spoke X Core	100960						101893				100768	100769	
20" Spinergy Spox	See Spinergy Spox Page	6			Screw: 100669 Nut: 100657	103125		103179					
22" Spinergy Spox						100827		100808		100889	100888		
24" Spinergy Spox						100766	101161	100615		100830	100828		
25" Spinergy Spox						100767	101160	100765		101464	101460		
26" Spinergy Spox						101477		101148		200200	100950		
22" Spinergy LX	See Spinergy LX Page	6			Screw: 100669 Nut: 100657	100827		100808		100889	100888		
24" Spinergy LX						100766	101161	100615		100830	100828	200213	
25" Spinergy LX						100767	101160	100765		101464	101460	200208	
26" Spinergy LX						101477		101148		200200	100950		
18" MAXX Spoke	109285	6	100698 (Not used w/ Superlight Handrim)	100653 (Not used w/ Superlight Handrim)	Aluminum, Plastic Ctd, Projection: 100654 Natural Fit: 100835 Superlight: Screw: 100669 Nut: 100657	100206		101106					
20" MAXX Spoke	108244					200536		200542					
22" MAXX Spoke	105135					100560		100576	100569	200538	200201		
24" MAXX Spoke	105136					200349	101161	100577	200547	100793	200202		
25" MAXX Spoke	107436					200350		101870	200548	200539	200540		
26" MAXX Spoke	107437					200351		100578	200549	100907	101454		

* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 000240 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels or part number 100666 for the screw on 5-Spoke X Core wheels.