


## **Jay Soft Combi P Cushion Measuring Guide**



PART NUMBER	CHAIR SIZE	ACTUAL WIDTH	ACTUAL LENGTH	OVERALL HEIGHT
		A	B	C
2200P	18 x 16	18.00	16.00	3.13
2201P	18 x 18	18.00	18.00	3.00
2202P	20 x 18	20.00	18.00	3.25
2203P	18 x 20	18.00	20.00	3.00
2205P	16 x 16	15.50	16.00	3.25
2206P	16 x 18	15.50	18.00	3.00
2207P	16 x 20	16.00	20.00	3.13
2210P	14 x 14	14.00	14.00	3.25
2220P	22 x 20	21.50	19.81	3.07
2228P	22 x 18	21.50	17.75	3.07
2240P	24 x 20	23.50	19.75	3.07
2248P	24 x 18	23.50	17.69	3.07

All listed measurements are in inches