




Discontinued Focus CR Caster Forks and Stems



| | |
|--------|------|
| Std/HD | 3.5" |
| Tall | 4.7" |



| Fork | Height | Width |
|---------|--------|-------|
| 5" Std | 4.7" | 2.1" |
| 6" Std | 5.7" | 2.1" |
| 6" Wide | 5.8" | 2.9" |
| 7" Std | 7.0" | 2.1" |
| 7" Wide | 7.0" | 2.9" |



| Item Number | Part Number | Description | Quantity | Retail |
|-------------|-------------|----------------------------------------------------------------------------------------------------------|----------|----------|
| 1a, 2, 5, 6 | 503816 | Fork Assembly 5" Std Stem, Tilt | 1 | \$58.65 |
| 1b, 2, 5, 6 | 503814 | Fork Assembly 6" Std Stem, Tilt | 1 | \$58.65 |
| 1c, 2, 5, 6 | 503817 | Fork Assembly 6" Wide Std Stem, Tilt <i>Only used with 2" wide casters</i> | 1 | \$82.14 |
| 1d, 2, 5, 6 | 503815 | Fork Assembly 7" Std Stem, Tilt | 1 | \$58.65 |
| 1e, 2, 5, 6 | 503818 | Fork Assembly 7" Wide Std Stem, Tilt <i>Only used with 2" wide casters</i> | 1 | \$82.14 |
| 1e, 3, 5-7 | 503122 | Fork Assembly 7" Wide Heavy Duty Std Stem Focus <i>Only used with 2" wide casters</i> | 1 | \$213.00 |
| 8a, 9, 10a | 502653 | Caster Attachment Hardware Kit CAT/FCS | 1 | \$7.88 |
| 8b, 9, 10b | 502654 | Caster Attachment Hardware Kit for 6" and 7" Wide Fork, CAT/FCS <i>Only used with 2" wide casters</i> | 1 | \$17.29 |
| 1a, 5 | 503823 | Caster Fork 5" Focus | 1 | \$44.95 |
| 1b, 5 | 503824 | Caster Fork 6" Focus | 1 | \$44.95 |
| 1c, 5 | 503825 | Caster Fork 6" Wide Focus <i>Only used with 2" wide casters</i> | 1 | \$65.12 |
| 1d, 5 | 503826 | Caster Fork 7" Focus | 1 | \$44.95 |
| 1e, 5 | 503827 | Caster Fork 7" Wide Focus <i>Only used with 2" wide casters</i> | 1 | \$65.12 |
| 2, 5 | 503819 | Caster Stem Standard, Focus | 1 | \$17.82 |
| 3, 5 | 503820 | Caster Stem Standard HD, Focus | 1 | \$108.41 |
| 5 | 100529 | M12 x 1.25mm Lock Nut Black | 2 | \$1.34 |
| 6 | 100682 | Washer 1/2" .518 x .875 x .047 F/W Black Zinc | 3 | \$1.05 |
| 7 | 002593 | Caster Stem Washer Heavy Duty <i>Used with Heavy Duty</i> | 2 | \$19.27 |
| 8a | 101829 | M8 x 65 HHCS Blk Zn | 1 | \$3.15 |
| 8b | 101831 | M8 x 80 HHCS Blk Zn <i>Used with wide forks</i> | 1 | \$3.07 |
| 9 | 100659 | M8 Nylock Nut Blk Zn | 1 | \$1.05 |
| 10a | 100792 | Caster Spacer D8 x D12 x 10mm <i>Used with all casters except 6" x 2" Poly</i> | 2 | \$2.14 |
| 10b | 001246 | Caster Spacer Black 6 x 2 Poly .472ODX.319IDX.370LG | 2 | \$5.84 |

| Item Number | Part Number | Description | Quantity | Retail |
|-------------|-------------|---------------------------------------|----------|--------|
| | | <i>Only used with 6" x 2" casters</i> | | |