

# Anterior Trunk Supports

Anterior Trunk Supports encourage an upright posture and support the torso for improved respiratory function, head positioning, and upper extremity control. Options include six support styles and either structured or dynamic construction.

## CLASSIC STYLE



	DYNAMIC	STRUCTURED	A	B
EI	ATS8107	–	6.5"	9.75"
XS	ATS8100	ATS8300	7.5"	12.5"
S	ATS8101	ATS8301	8.5"	15"
M	ATS8102	ATS8302	9.5"	19"
L	ATS8103	ATS8303	10.5"	22.5"
XL	ATS8104	ATS8304	11.5"	25.5"

The Classic style is designed for optimal distribution of support and maximum stability.