


# Catalyst Straight Backrest



*If Caregiver Push Handles are present, please visit either the Catalyst 5 – 5Ti Stright No Push Handles and Caregiver Height Adjustable Stroller Handles page or Catalyst 4 Straight No Push Handles and Caregiver Height Adjustable Stroller Handles page If Anti-Rollback is present, see Back Tubes for Anti-Rollback System*

Item Number	Part Number	Description	Quantity
1a, 2	100718	Back Post Assm Short Straight CAT	1
1b, 2	100719	Back Post Assm Medium Straight CAT	1

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
1c, 2	100720	Back Post Assm Tall Straight CAT	1
2	100803	Tube End Plug 7/8"	1
3	100681	M6 x 40 CLS8.8 DIN931 HHCS ZC	1
4	100746	M6 Flat Washer Blk Zn	1
5	100658	M6 Nylock Nut Blk Zn	1