


Ethos Padded Velcro Adjustable Calf Strap



Sizes based on seat width

Item Number	Part Number	Description	Quantity
1	004462	Calf Strap Clamp A	2
2	004463	Calf Strap Clamp B	2
3	004464	Calf Strap Rod	2
4	103589	6-13 x 3/4 Phillips Pan HD, Plastite 45, BLK ZN Screw	4
1-4	108224	Calf Strap Clamp Assembly, Ethos	1
5a	101142	Calf Strap Padded, Short 10-14	1
5b	101143	Calf Strap Padded, Medium 15-17	1
5c	101144	Calf Strap Padded, Long 18-24	1
1-4, 5a	108228	Padded Velcro Adjustable Calf Strap with Clamp Assembly, Short 10-14	1
1-4, 5b	108229	Padded Velcro Adjustable Calf Strap with Clamp Assembly, Medium 15-17	1
1-4, 5c	108230	Padded Velcro Adjustable Calf Strap with Clamp Assembly, Long 18-22	1