

☰ **Comfort, convenience and value in an accommodating skin protection and positioning cushion.**



Reaction Force Reduction



Reduction of reaction forces at the maximum points of compression

Tension Reduction



Materials and design that promote greater immersion

Pre-contouring



Pre-contoured cushions help redistribute loads

Foam Density vs Foam Firmness



Dense foams have greater longevity

JAY® ION™ CUSHION

≡ The JAY Ion™ cushion is a comfortable and convenient foam cushion well suited for the client requiring moderate skin protection and positioning.

Comfortable

- Dual layered contoured foam base
- Reticulated foam in the outer cover promotes air movement
- Soft, stretchable outer cover

Convenient

Outer cover

- Machine washable
- Quick to dry


Inner cover:

- Easy to wipe off
- Aquaguard™ zipper resists moisture
- Anti-wicking seam thread protects foam base


Anti-Microbial

Silver impregnated X-static® stretch fabric

- Naturally inhibits bacteria growth
- Safe non-toxic material



Dual layered foam base



Outer and inner covers



Anti-microbial silver thread

SPECIFICATIONS

JAY®

JAY ION Cushion

HCPCS Code	Designed to meet specifications of E2607 & E2608
Weight	2.1 lbs. (12" x 12"), 3.1 lbs. (18" x 18")
Product Width	14" to 24"
Product Depth	14" to 24"
Product Height	4"
Weight Capacity	300 lbs. 14" to 21" widths, 500 lbs. 22" to 24"
Outer Cover	X-static®