

SIGVARIS

LIFE FOR LEGS

For 230 Cotton, 360 Cushioned Cotton, 770 Truly Transparent, and 860 Select Comfort
 (For 500 Natural Rubber, see 500 Sizing Chart)

	Small				Medium				Large				X-Large (Calf only)*			
Ankle Circumference (b)	7" - 8 1/2" (18 - 21.5cm)				8 1/2" - 10" (21.5 - 25.5cm)				10" - 11 1/2" (25.5 - 29cm)				11 1/2" - 13" (29 - 33cm)			
Calf Circumference (c)	Average		Full		Average		Full		Average		Full		Average		Full	
	11" - 13" (28 - 33cm)		13" - 15" (33 - 38cm)		13" - 15 1/2" (33 - 39.5cm)		15 1/2" - 18" (39.5 - 45.5cm)		15" - 18" (38 - 45.5cm)		18" - 21" (45.5 - 53.5cm)		17" - 20" (43 - 51cm)		20" - 23"	
A-D Calf Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	Up to 16" (41cm)	Over 16" (41cm)	Up to 16" (41cm)	Over 16" (41cm)	Up to 16" (41cm)	Over 16" (41cm)	Up to 16" (41cm)	Over 16" (41cm)	Up to 16" (41cm)	Over 16" (41cm)	Up to 16" (41cm)	Over 16" (41cm)	Up to 16" (41cm)	Over 16" (41cm)	Up to 16" (41cm)	Over 16" (41cm)
Thigh Circumference (g)	Average		Full		Average		Full		Average		Full		<i>*X-Large not available in 770 Truly Transparent Series</i>			
	16" - 19 1/2" (40.5 - 49.5cm)		19 1/2" - 23" (49.5 - 58.5cm)		19" - 23" (48.5 - 58.5cm)		23" - 27" (58.5 - 68.5cm)		23" - 27 1/2" (58.5 - 70cm)		27 1/2" - 32" (70 - 81.5cm)					
A-G Thigh Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	Up to 28" (71cm)	Over 28" (71cm)	Up to 28" (71cm)	Over 28" (71cm)	Up to 28" (71cm)	Over 28" (71cm)	Up to 28" (71cm)	Over 28" (71cm)	Up to 28" (71cm)	Over 28" (71cm)	Up to 28" (71cm)	Over 28" (71cm)	Up to 28" (71cm)	Over 28" (71cm)	Up to 28" (71cm)	Over 28" (71cm)
Guidelines for Hip Circumference	S1	S2	S3	S4	M1	M2	M3	M4	L1	L2	L3	L4				
	Up to 42" (106.5cm)		Up to 48" (122cm)		Up to 48" (122cm)		Up to 54" (137cm)		Up to 55" (138cm)		Up to 61" (155cm)					

Guideline for Shoe Measurement for the SIGVARIS system

	Women		Men	
	Short Leg	Long Leg	Short Leg	Long Leg
Shoe Size	5 1/2 - 9 1/2	8 - 12	7 - 11	9 1/2 - 14