


# Catalyst Standard Back Upholstery



Standard Back Upholstery Compatibility		
Back Post	Height Range	Back Height
Standard	Short	8.5" - 11.5"
	Medium	12.5" - 15.5"
	Tall	16.5" - 19.5"
Angle Adjustable	Short	14" - 16"
	Medium	16" - 18"
	Tall	18" - 20"

The Catalyst 4 only utilizes the "Tall" back upholstery options. The Catalyst 4C utilizes all back upholstery options <p>Attaching hardware is not included with the back upholstery

Item Number	Part Number	Description	Quantity	Retail
3, 4, 5	503727	Zip Tie Kit <i>This kit contains one zip tie, one bolt and one nut.</i>	1	\$4.73
3, 4, 5	503745	Zip Tie Kit PAIR	1	\$9.46
1a	100145	Back Upholstery Short 14" Wide CAT	1	\$115.57
1b	100148	Back Upholstery Short 15" Wide CAT	1	\$115.57
1c	100151	Back Upholstery Short 16" Wide CAT	1	\$115.57
1d	100154	Back Upholstery Short 17" Wide CAT	1	\$115.57
1e	100157	Back Upholstery Short 18" Wide CAT	1	\$115.57

Item Number	Part Number	Description	Quantity	Retail
1f	100160	Back Upholstery Short 19" Wide CAT	1	\$126.07
1g	100163	Back Upholstery Short 20" Wide CAT	1	\$126.07
1h	100146	Back Upholstery Medium 14" Wide CAT	1	\$115.57
1i	100149	Back Upholstery Medium 15" Wide CAT	1	\$115.57
1j	100152	Back Upholstery Medium 16" Wide CAT	1	\$115.57
1k	100155	Back Upholstery Medium 17" Wide CAT	1	\$115.57
1l	100158	Back Upholstery Medium 18" Wide CAT	1	\$115.57
1m	100161	Back Upholstery Medium 19" Wide CAT	1	\$126.07
1n	100164	Back Upholstery Medium 20" Wide CAT	1	\$147.08
1o	100147	Back Upholstery Tall 14" Wide CAT	1	\$115.57
1p	100150	Back Upholstery Tall 15" Wide CAT	1	\$115.57
1q	100153	Back Upholstery Tall 16" Wide CAT	1	\$115.57
1r	100156	Back Upholstery Tall 17" Wide CAT	1	\$115.57
1s	100159	Back Upholstery Tall 18" Wide CAT	1	\$115.57
1t	100162	Back Upholstery Tall 19" Wide CAT	1	\$126.07
1u	100165	Back Upholstery Tall 20" Wide CAT	1	\$157.59
1v	100874	Back Upholstery Tall 21" Wide CAT	1	\$157.59
1w	100246	Back Upholstery Tall 22" Wide CAT	1	\$157.59
2	100734	M5 x 16 PHMS ZC	2	\$1.58
3	100685	Zip Tie, Black	2	\$1.58
4	100666	M5 x 45 BTN SHCS BLKZC	2	\$1.58
5	100657	M5 Nylock Nut Blk Zn	2	\$0.53