



Züm-AP / SP

Folding Mobile Shower
Commode Chairs

The Züm's unique crossbrace design provides better toilet alignment and frame rigidity. Two models available: Züm-AP (attendant-propel) and Züm-SP (self-propel) configurations.



FEATURES

- **Unique Folding Frame** – Low-profile cross brace is contoured to fit around the toilet, allowing for better rearward positioning of seat aperture over the bowl

- **Lightweight aluminum frame** facilitates transportability

- **Removable back assembly** to facilitate transportation

- 17" (43 cm) back height with molded, vinyl backrest

- Compatible with all standard-width Raz seats, which include the exclusive IPAS¹

- Flip-up arms with optional locks



Züm-SP:

- 22" (56 cm) or 24" (61 cm) wheels
- 3" (7.5 cm) fore/aft rear wheel center-of-mass adjustment

SPECIFICATIONS	Züm-AP	Züm-SP
WEIGHT CAPACITY	350 lbs (160 kg)	350 lbs (160 kg)
OVERALL WIDTH – Frame Open	23.25" (59 cm)	25.25" (64 cm)
– Frame Folded	11" (28 cm)	13" (33 cm)
– Frame Folded	7.5" (19 cm) – (back assembly removed)	7.5" (19 cm) (back assembly and wheels removed)
OVERALL LENGTH	34" (86 cm)	40" (102 cm) with 22" (56 cm) wheels
- with standard footrests		41" (104 cm) with 24" (61 cm) wheels
- for SP, rear wheels in center position		
SEAT WIDTH	18" (46 cm)	18" (46 cm)
CLEARANCE HEIGHT UNDER SEAT		
– with pan brackets	17.25" (43 cm)	17.25" (43 cm)
– without pan brackets	17.75" (45 cm)	17.75" (45 cm)
ARMREST HEIGHT ABOVE SEAT	8" (20 cm) ²	8" (20 cm) ²
SEAT-TO-FLOOR HEIGHT	21.5" (55 cm) ²	21.5" (55 cm) ²
PRODUCT WEIGHT	44 lbs (20 kg)	57 lbs (26 kg)
FRAME MATERIALS		
– base frame	Aluminum ³	Aluminum ³
– back frame and footrests	Stainless Steel ³	Stainless Steel ³
– hardware	Stainless Steel ³	Stainless Steel ³

¹ Ischial Pelvic Alignment System (IPAS) allows 2" (5 cm) of fore/aft seat depth adjustment

² Measurements taken with standard molded seat and standard arm pads

³ Lifetime warranty against corrosion

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications and options are subject to change without notice.