

# pango

## Activity Chair

*Seating to boost positioning and posture*

Available in 2 sizes, the **pango** is an attractive activity chair designed for kids with mild to moderate positioning needs. It is designed to promote better posture and improves the stability and comfort of the user.

In addition to standard key features of a contoured seat and back, angle adjustable back, height adjustable seat, and pelvic belt, the **pango** has the unique feature of *seat angle adjustability*.

**Below is a list of our range of accessories which are catered to the individual's needs and circumstances.**

*Headrest*

*Adjustable laterals*

*Foot ankle positioner*

*Stabilizer (can also be used as rocker)*

Trunk harness

Hip guides

Abductor

Tray

Mobility wheels

Adjustable footrest

Guide bar

