


Rogue2 Frog Legs Suspension Fork (Ultra Sport Fork)



Item Number	Part Number	Description	Quantity
1, 2, 4-6	114194	Frog Legs Ultra Sport Fork with 1.25" Stem, Rogue2	1
4-6	103201	Frog Legs Suspension Axle w/Hdwre PAIR, RG/Clik/Ethos <i>This hardware may incur a longer lead time</i>	1
3a	105053	Frog Legs Polymer < 80lbs	1
3b	105054	Frog Legs Polymer 80 - 120lbs	1
3c	105055	Frog Legs Polymer 120 - 160lbs	1
3d	105056	Frog Legs Polymer 160 - 220lbs	1
3e	105057	Frog Legs Polymer 220 - 260lbs	1
3f	105058	Frog Legs Polymer 260 - 300lbs	1