



Phone (877) 227-0097 | Fax (704) 882-0751

Manatee

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/20/2025

Manatee

	0	1	2	3
Effective seat width	330 mm (13")	330 mm (13")	330 mm (13")	400 mm (15¾")
Effective seat depth	220-250 mm (8¾-9¾")	300-370 mm (11¾-14½")	390-460 mm (15¼-18")	400-470 mm (15¾-18½")
Backrest height	430 mm (17")	430 mm (17")	580 mm (22¾")	640 mm (25¼")
Lower leg length	280 mm (11")	280 mm (11")	350 mm (13¾")	350 mm (13¾")
Seat surface height at front edge	140-430 mm (5½-17")	140-430 mm (5½-17")	140-430 mm (5½-17")	140-430 mm (5½-17")
Overall width	420 mm (16½")	420 mm (16½")	420 mm (16½")	490 mm (19¼")
Folded length	1030-1100 mm (40½-43¾")	1100-1170 mm (43-46")	1360-1430 mm (53½-56¾")	1440-1510 mm (56¾-59½")
Folded height	130 mm (5")	130 mm (5")	130 mm (5")	130 mm (5")
Mass device	5.5 kg (12.1 lb)	5.5 kg (12.1 lb)	6 kg (13.2 lb)	7 kg (15.4 lb)
Maximum user mass	50 kg (110.2 lb)	50 kg (110.2 lb)	75 kg (165.3 lb)	75 kg (165.3 lb)
Maximum rated load	57.5 kg (126.7 lb)	57.5 kg (126.7 lb)	83 kg (182.9 lb)	84 kg (185.1 lb)

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1259.1 Pricing and specifications subject to change without notice.



Manatee Seat and Frame Options

Size and Color

Shower seat only



Blue

Size 0, Blue (882005)	\$765
Size 1, Blue (882015)	\$765
Size 2, Blue (882025)	\$765
Size 3, Blue (882035)	\$835

Red

Size 1, Red (882014)	\$765
Size 2, Red (882024)	\$765
Size 3, Red (882034)	\$835

Manatee Frame Options

Additional frame used for mounting bath chair onto additional frame



Rolling Shower Frame

Rolling Shower Frame (882100)	\$2499
-------------------------------	--------



Tub Slider Frame

Tub Slider without Manatee W/ Casters, Size 1 and 2 (MC6000PED)	\$3563
Tub Slider without Manatee W/ Casters, Size 3 (MC6200PED)	\$4343



Basic Manatee Frame

Basic Manatee Frame (882101)	\$382
------------------------------	-------

Accessories and Straps

Accessories




Leg Straps

Leg Straps, Size 0,1 and 2, Blue (883055)	\$84
Leg Straps, Size 0,1 and 2, Red (883054)	\$84
Leg Straps, Size 3, Blue (883105)	\$107
Leg Straps, Size 3, Red (883104)	\$107



Head/ Lateral Support

Head Lateral Support, Size 0,1 and 2, Blue (883015)	\$121
Head Lateral Support, Size 0,1 and 2, Red (883014)	\$121
Head Lateral Support, Size 3, Blue (883065)	\$121
Head Lateral Support, Size 3, Red (883064)	\$121



Pommel

Pommel, Size 0,1,2, Blue (883045)	\$107
Pommel, Size 0,1,2, Red (883044)	\$107
Pommel, Size 3, Blue (883095)	\$107
Pommel, Size 3, Red (883094)	\$107



Torso Belt

Torso Belt, Size 0,1 and 2, Blue (883025)	\$84
Torso Belt, Size 0,1 and 2, Red (883024)	\$84
Torso Belt, Size 3, Blue (883075)	\$84
Torso Belt, Size 3, Red (883074)	\$84