


Focus CR Seat Frame - Growth



Offset seat tube (1d-1i)

Has your seat width and/or depth changed? Will your current base strut and seat tube work? If not, please select accordingly.

Item Number	Part Number	Description	Quantity	Retail
1a	001038	Seat Tube Weld Front Short TILT <i>Accommodates 14" - 18" Depths</i>	2	\$89.67
1b	001039	Seat Tube Weld Front Medium TILT <i>Accommodates 16" - 20" Depths</i>	2	\$89.67
1c	001040	Seat Tube Weld Front Long HD <i>Accommodates 18" - 22" Depths</i>	2	\$151.99
1d	001046	Seat Tube Weld Front Offset Right Short TILT <i>Accommodates 14" - 18" Depths with Offset Hangers</i>	1	\$89.67

Item Number	Part Number	Description	Quantity	Retail
1e	001047	Seat Tube Weld Front Offset Right Medium TILT <i>Accommodates 16" - 20" Depths with Offset Hangers</i>	1	\$89.67
1f	001048	Seat Tube Weld Front Offset Right Long HD TILT <i>Accommodates 18" - 22" Depths with Offset Hangers</i>	1	\$151.99
1g	001049	Seat Tube Weld Front Offset Left Short TILT <i>Accommodates 14" - 18" Depths with Offset Hangers</i>	1	\$89.67
1h	001050	Seat Tube Weld Front Offset Left Medium TILT <i>Accommodates 16" - 20" Depths with Offset Hangers</i>	1	\$89.67
1i	001051	Seat Tube Weld Front Offset Left Long HD TILT <i>Accommodates 18" - 22" Depths with Offset Hangers</i>	1	\$151.99
2a	001011	Growth Strut, Base, Small Tilt <i>Accommodates 14" - 18" Widths</i>	2	\$18.91
2b	001012	Growth Strut, Base, Medium Tilt <i>Accommodates 16" - 20" Widths</i>	2	\$18.91
2c	001013	Growth Strut, Base, Large Tilt <i>Accommodates 18" - 22" Widths</i>	2	\$18.91
2d	112005	Growth Strut, Base, X-Large Tilt <i>Accommodates 20" - 24" Widths</i>	2	\$21.00
2a, 3, 4, 5	501126	Seat Frame Strut Assm Small Focus <i>Accommodates 14" - 18" Widths</i>	1	\$47.44
2b, 3, 4, 5	501127	Seat Frame Strut Assm Medium Focus <i>Accommodates 16" - 20" Widths</i>	1	\$47.44
2c, 3, 4, 5	501128	Seat Frame Strut Assm Large Focus <i>Accommodates 18" - 22" Widths</i>	1	\$47.44
2d, 3, 4, 5	114276	Seat Frame Strut Assm X-Large Focus <i>Accommodates 20" - 24" Widths</i>	1	\$63.04