



COMPARE & SAVE!

A less expensive alternative!



Otter Bathing System

- OT 1000 Small, 1/cs
- OT 2000 Medium, 1/cs
- OT 3000 Large, 1/cs



FEATURES:

- Plastic frame.
- Fabric can be removed and machine-washed.
- Seat fabric can be slipped off to bathe a child sitting on the floor of the tub, supported by the bath chair back.
- Seat and back are angle adjustable and have five angle adjustments at 0, 22.2, 45, 67.5 & 90.
- Adjustable slip-resistant legs raise chair up to 7". Height can be adjusted with child in chair.
- **NEW!...** Uni-Bars allow seat and back angle adjustments to be made with one hand.
- Leg straps control abduction and adduction.
- Two positioning straps secure child.
- Head cushions are height and width adjustable and provide support and comfort. They can also be used as lateral supports.
- Folds flat for easy storage.
- Two year limited warranty.



DIMENSIONS	OT 1000	OT 2000	OT 3000
SEAT DEPTH	13"	13"	16.5"
SEAT WIDTH	14.5"	14.5"	14.5"
BACK HEIGHT	18"	25"	32"
SEAT HEIGHT	2" - 7"	2" - 7"	2" - 7"
OVERALL OUTSIDE WIDTH	17"	17"	17"
OVERALL OUTSIDE LENGTH	30"	30"	30"
USER HEIGHT (APPROX)	Up to 36"	32" - 50"	46" - 68"
WEIGHT LIMIT	60 lbs.	120 lbs.	160 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.