



Zippie TS

March 2023

EIZ1A

Date: _____	SUBMITTING FOR: ☐ Quote ☐ Order
Dealer Acct #: _____	
Dealer: _____	PO#: _____
Dealer Contact: _____	Ship To: _____
Dealer Address: _____	Attention: _____
Dealer City: _____ ST: _____ ZIP: _____	Address: _____
Dealer Phone: () _____ Fax: () _____	Address: _____
Confirm Via: ☐ Fax ☐ Email	Ship To City: _____ ST: _____ ZIP: _____
Mark For: _____	Ship To Phone: () _____ Fax: () _____

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
ZIPPIE TS									

1. FRAME

All Zippie TS models grow in width and depth.

Base Models			
☐	187WM63	E1234	Folding Growing Zippie TS Without Seating Transit Included (includes True Fit) \$3,347
☐	187WM65	E1234	Folding Growing Zippie TS Without Seating Transit Optional (includes True Fit) \$3,185
☐	187WM64	E1233	Rigid Growing Zippie TS Without Seating Transit Included (includes True Fit) \$3,214
☐	187WM66	E1233	Rigid Growing Zippie TS Without Seating Transit Optional (includes True Fit) \$3,046
☐	187F02		Additional Rigid Strut Tube NC

Note: HCPCS Codes shown as submitted. Some sizes fit other HCPCS codes i.e. 15"-18" widths = E1161


Transit System (For use when Transit Optional is selected)			
Tie Down Requires: Solid Seat and Solid Back or Growing Seat Pan and JAY Fit or Conventional Seat and Back or Growing Seat Pan and Solid Back or JAY Fit and/or Conventional Seat & Back			
☐	187TO12	K0108	Transit Optional Seating (Req. Selection of Seat & Back) \$328
☐	187TO11	K0108	Transit Optional w/o Seating \$328
☐	187TO20		Transit Omit N/C

Tilt Options			
☐	187TT12		Dual Trigger Cable Release NC
☐	187VZ01		Crisscross Cables (cables are longer individually) NC
☐	187VZ02		Non-crisscross Std
☐	187TT8		Cable Foot Release (may limit folding width) NC

2. SIZE/SEATING

Frame Width			
☐	187W25		Small Cross Brace NC
☐	187S1	10"	NC
☐	187S2	11"	NC
☐	187S3	12"	NC
☐	187W26		Medium Cross Brace NC
☐	187S3	12"	NC
☐	187S4	13"	NC
☐	187S5	14"	NC
☐	187W27		Large Cross Brace NC
☐	187S5	14"	NC
☐	187S6	15"	NC
☐	187S7	16"	NC
☐	187W28		X-Large Cross Brace NC
☐	187S7	16"	NC
☐	187S8	17"	NC
☐	187S9	18"	NC

Frame Depth and Back Post Location			
☐	187FD4		13" Frame Depth NC
☐	187FD5		14" Frame Depth NC
☐	187FD6		15" Frame Depth NC
☐	187FD7		16" Frame Depth NC
☐	187FD8		17" Frame Depth NC
☐	187FD9		18" Frame Depth NC
☐	187FD10		19" Frame Depth (uses rear frame extension plug) NC
☐	187FD36		20" Frame Depth (uses rear frame extension plug) NC



Units are in inches.

2. SIZE/SEATING CONT.

Seats			
☐	187ST55	E2231 Solid Seat (Requires Solid Back, N/A w/ 10"-13" Wide x 19"-20" Frame Depth or 17"-18" Wide x 13"-17" Frame Depth)	\$224
☐	187ST14	E2231 Growing Seat Pan Depth (Fig. 7.1)	\$224
☐	187ST5	Omit Seat	Std

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list. Cushion name or part number must be filled out with desired width and depth.

☐	187JC20	JAY Bundled Cushion	Price List
		Name / Part #	
		Width (in.): Depth (in.):	

REFERENCE ONLY

Seat Depth Pan Size	Depth Range	
	Minimum	Maximum
Medium 13"-14"	11-3/4"	16-1/4"
Long 15"-17"	13-3/4"	18-1/4"
X-Long 18"-20"	15-3/4"	20-1/4"




Figure 7.1 - Frame Depth Growth Range Matrix

3. CASTERS

Caster Wheel

Some interference may occur when choosing 8" Caster w/ Z-Finity

☐	187CW30	3" Micro Light Up	\$72
☐	187CW31	4" Micro Light Up	\$78
☐	187CW51	5" Micro Light Up	\$96
☐	187CW5	5" Low-Profile Polyurethane	STD
☐	187CW7	E2219 5" x 1.5" Semi-Pneumatic	\$96
☐	187CW33	E2214 E2215 6" Pneumatic	\$84
☐	187CW34	6" Polyurethane	NC
☐	187CW35	E2219 6" x 1.5" Semi-Pneumatic	\$96
☐	187CW21	E2214 E2215 8" Pneumatic	\$96
☐	187CW32	8" Polyurethane	NC
☐	187CW39	E2219 8" x 1.5" Semi-Pneumatic	\$96
☐	187CW36	E2214 E2215 8" x 2" Pneumatic	\$162
☐	187CW37	E2214 K0108 8" x 2" Pneumatic w/ Airless Insert	\$252

Caster Options

☐	187CS3	K0073 Caster Pin Locks (N/A with 8" x 2" and 8" x 1.5" casters)	\$218
--------------------------	--------	--	-------

CASTER SIZE

Rear Wheel Size	3"	4"	5"	6"	8"x1"	8"x1.5"	8"x2"
☐ 12"	☐ 17.0"	☐ 16.5"	☐ 16.5"	☐ 17.0"	N/A	N/A	N/A
	☐ 17.5"	☐ 17.0"	☐ 17.0"	☐ 17.5"	N/A	N/A	N/A
	N/A	☐ 17.5"	☐ 17.5"	N/A	☐ 18.5"	☐ 18.0"	☐ 18.5"
☐ 12" Low Poly	☐ 16.0"	☐ 16.0"	☐ 16.0"	☐ 16.5"	N/A	N/A	N/A
	☐ 16.5"	☐ 16.5"	☐ 16.5"	N/A	N/A	N/A	N/A
	☐ 15.5"	☐ 15.5"	☐ 15.5"	N/A	N/A	N/A	N/A
☐ 16"	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.5"	☐ 18.0"	☐ 18.5"
	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 19.0"	☐ 18.5"	N/A
	☐ 15.5"	☐ 15.5"	☐ 15.5"	N/A	N/A	N/A	N/A
☐ 20"	☐ 16.0"	☐ 16.0"	☐ 16.0"	N/A	N/A	N/A	N/A
	☐ 16.5"	☐ 16.5"	☐ 16.5"	☐ 16.5"	N/A	N/A	N/A
	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 18.5"	☐ 18.0"	☐ 18.5"
☐ 22"	☐ 17.5"	☐ 17.5"	☐ 17.5"	☐ 17.5"	N/A	N/A	N/A
	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.5"	☐ 18.0"	☐ 18.5"
	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 19.0"	☐ 18.5"	☐ 19.0"
☐ 24"	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 19.0"	☐ 19.0"	N/A
	N/A	☐ 19.0"	☐ 19.0"	☐ 19.0"	☐ 19.5"	N/A	N/A
	N/A	N/A	☐ 19.5"	☐ 19.5"	☐ 19.5"	☐ 19.5"	N/A

Chart A - Seat to Floor Height Matrix

Note: Rear seat-to-floor ht. may be up to 1.5" lower than front seat-to-floor ht.

4. REAR WHEELS / SEAT HEIGHT

Rear Wheel Size			
☐	187RW13	12" Mag (N/A with handrims)	STD
☐	187RW14	16" Mag (N/A with handrims)	NC
☐	187RW3	20" Mag	NC
☐	187RW8	22" Mag	NC
☐	187RW1	24" Mag	NC
☐	187RW90	20" Lite Spoke	\$224
☐	187RW91	22" Lite Spoke	\$224
☐	187RW54	24" Lite Spoke	\$224

Rear Wheel Tire			
☐	187RT1	E2211 E2212	Pneumatic \$96
☐	187RT2	E2211 E2213	Pneumatic Airless Insert \$96
☐	187RT3		Full Polyurethane (N/A w/ 12" and 16" rear wheels) NC
☐	187RT4		Low Polyurethane (Only available with 12" or 24" wheels) NC
☐	187RT6	E2211 E2212	High Pressure Clincher (Only available with 24" spoke wheels) \$106

Handrim			
☐	187HR1	Aluminum Anodized	NC
☐	187HR4	Plastic Coated	\$123
☐	187HR7	Omit (N/A with wheel/axle credit)	NC

Projections (N/A on Plastic Coated)			
*Only available with 24" wheels & N/A with spoke guards.			
☐	187HR8	8 Oblique	\$363

Axle			
☐	187AX26	One-Arm Drive Axle (Only avail. w/ One-Arm Drive Option)	STD
☐	187AX1	Quick Release	NC
☐	187AX11	Quad Quick Release (N/A with Lite Spokes)	\$141

Axle Plate			
☐	187AL18	12" Rear Wheel Axle Plate (Only avail. w/ 12" or 16" rear wheels)	STD
☐	187AL34	Short Adjustable (N/A with 12" rear wheels)	NC
☐	187AL35	Long Adjustable (N/A with 12" rear wheels)	NC
☐	187AL17	E0958 One Arm Drive (Avail. w/ 20", 22", 24" Mag Rear Wheels & 13"-18" width)	\$1,006
☐	187OA2	Standard Handrim	NC
☐	187HM1	Right	NC
☐	187HM2	Left	NC
☐	187OA1	Small Handrim (N/A w/ 20" wheels)	NC
☐	187HM1	Right	NC
☐	187HM2	Left	NC

5. WHEEL LOCK / HUB LOCK

Wheel Lock			
Wheel Locks can be ordered in combination with Hub Lock option			
☐	187WL1	High Mount Push	NC
☐	187WL2	High Mount Pull	NC
☐	187WL25	Foot Lock (may limit folding width)	NC
☐	187WL10	Omit	NC

Extension Handle			
☐	187EH1	E0961 6" Pair	\$84

Hub Lock			
*Only available with 12" or 16" wheels. Not avail. w/ One-Arm Drive.			
☐	187LK01	12" Hub Lock	\$111
☐	187LK02	16" Hub Lock	\$111

6. FOOTRESTS

Hanger Mounted Footrest (Footplates Mount to Hanger Tubing)

Hanger Type

☐	187H140		70° Swing In/Out Front Mount	NC
☐	187H142		80° Swing In/Out Front Mount	NC
☐	187H158		90° Swing In/Out Front Mount	NC
☐	187H9		70° HD Front Mount Lift-Off	NC
☐	187H10		80° HD Front Mount Lift-Off	NC
☐	187H11		90° HD Front Mount Lift-Off	NC
☐	187H161	E0990	90° Elevating (N/A w/ 10"-11" Wide or Z-Finity)	\$352

Extensions (may interfere with some seat heights & casters)

☐	187E23		2" (not mounted)	\$61
☐	187E24		4" (not mounted)	\$61

Footplate

☐	187F9	K0040	Kids Angle Adjustable (N/A w/ 17"-18" widths)	\$207
☐	187F8	K0040	Adult Angle Adj. (Available with 10"-18" widths)	\$207
☐	187F62	K0040	Footboard (Available with 10"-16" widths)	\$212

Ext. Tube Mounted Footrests (Footplates Mount to Ext. Tube)

Hanger Type

☐	187H145		60° Swing In/Out	NC
☐	187H139		70° Swing In/Out	STD
☐	187H141		80° Swing In/Out	STD
☐	187H159	K0053	Articulating (N/A w/ 10"-13" composite footplates)	\$413
☐	187H160	E0990	Elevating (N/A w/ composite or platform footplates, or <16" width)	\$352

Extensions

☐	187E1		Short (11.5" - 14")	NC
☐	187E2		Med (14" - 17") - (N/A w/ ELR)	NC
☐	187E5		Long (17" - 19")	NC
☐	187E16		High Mount	\$84

Footplate

☐	187F1		Composite (N/A In 10", 11" or 12" width)	STD
☐	187F3	K0040	Aluminum Angle Adjustable	\$207
☐	187F4	K0040	Kids Angle Adjustable (N/A w/ 17"-18" width)	\$207
☐	187F5		Aluminum (Only available with elevating leg rests 14"-18")	NC
☐	187F48	K0040	Platform Flip-Up (Available on 60°, 70°, and 80° ext. mount hangers only)	\$212

Articulating / Elevating Hanger Cover

☐	187DI47		Neoprene Cover	\$380
☐	187DI48		Gel Pad Cover	\$447

Z-Finity Hangers/Footrest System

Some interference may exist in certain configurations.

☐	187H162		Z-Finity Hanger	STD
--------------------------	---------	--	-----------------	-----

Z-Finity Footrest Mount (Must Choose see Illustration in appendix)				
☐	187DI22		Direct Hanger Mt/Fixed Frame - High	NC
☐	187DI23		Direct Hanger Mt/Fixed Frame - Low	NC
☐	187DI24		Extension Tube Mount	NC

Z-Finity Footrest Pivot Arm

(Must Choose, use matrix in appendix for leg length)

☐	187DI20		Short 2"	NC
☐	187DI21		Long 3-1/4"	NC

Z-Finity Footplate


Some interference may exist in certain configurations.

☐	187F320	K0040	Multi Angle Adjustable - Pair (N/A w/ 8"-9" widths)	\$207
☐	187F321	K0040	Angle Adjustable Platform Footboard	\$168
☐	187F6		Omit Footplate	NC

Z-Finity Extension Tube Length (When Choosing Ext. Tube Mount)

☐	187E62		Hemi	NC
☐	187E1		Short	NC
☐	187E2		Medium	NC
☐	187E5		Long	NC

Direct Hanger Mount



Extension Tube Mount



Direct Hanger Mounted

Footrest Adaptor	Hanger Mounting Position	Min	Max
Short	High Mount	-1.3"	3.0"
	Low Mount	2.0"	6.3"
Long	High Mount	-2.5"	4.3"
	Low Mount	0.8"	7.5"

Extension Tube Mounted

Footrest Adaptor	Extension Tube Length	Min	Max
Short	Hemi Short	4.3"	12.0"
	Short	5.3"	13.0"
	Medium	7.0"	15.0"
	Long	9.0"	17.0"
Long	Hemi Short	3.0"	13.3"
	Short	4.5"	14.3"
	Medium	6.0"	16.3"
	Long	8.0"	18.3"

Footrest Hanger Options

☐	187IG2		Neoprene Covers (Avail. On 60° / 70°/80° ext. mount hangers only and Z-Finity except in highest mounted position).	\$101
☐	187HL1	E0951	Heel Loops (N/A w/ Z-Finity widths below 14")	\$45
☐	187A110	E0952	Toe Loop (Only avail. w/ adult angle adjust. Footplates N/A w/ Z-Finity)	\$45
☐	187HL5	K0038	Leg Strap	\$45
☐	187A55		Extra Leg Strap	\$40
☐	187H32		Omit Hanger	-\$15
☐	187F7		Omit Footplate & Ext. Tube	-\$15
☐	187HL4		Omit Heel Loop	NC
☐	187HL6		Omit Single Leg Strap	NC

7. BACKREST / SEATING

Backs

(Solid Back Depth = 2" Less Than Frame Depth)

☐	187PU4	Short Solid Back 12" (Requires solid seat)	NC
☐	187PU5	Med. Solid Back 14"	NC
☐	187PU6	Tall Solid Back 16"	NC
☐	187PU9	Omit Backrest	NC

JAY Backs (Optional)

Ordering JAY backs on-chair, See the JAY Price List. Requires Part #

☐	187JB02	JAY Bundled Backrest	Price List
		Backrest Name or Part #	
		Width (in.): Height (in.):	

Stroller Handle Back Posts - Backrests

Fixed Height

☐	187BP9	Low (19" Height)	NC
☐	187BP10	Medium (22" Height)	NC
☐	187BP11	Tall (25" Height)	NC
☐	187BP13	Removable Extension	\$96
☐	187RG1	Folding Adjustable Rigidizer	\$196


Height Adjustable

Rigid Frame w/ Angle Adjust Push Handle

☐	187B202	Std. (16"-22" Height)	NC
☐	187PH8	Std. Angle Adj. Push Handle Length	NC
☐	187PH9	Long Angle Adj. Push Handle Length	NC
☐	187B203	Tall (19"-25" Height)	NC
☐	187PH8	Std. Angle Adj. Push Handle Length	NC
☐	187PH9	Long Angle Adj. Push Handle Length	NC

Folding Frame with Push Handle

☐	187B204	Std. (17"-23" Height)	NC
☐	187BP47	Swing Away Handle (N/A with dual triggers)	\$61
☐	187B205	Tall (20"-26" Height)	NC
☐	187BP47	Swing Away Handle (N/A with dual triggers)	\$61



8. ARMRESTS / SIDEGUARDS

Armrest and Armpad - Height Adjustable

☐	187AA12	E0973	T-Post Regular w/ Foam	\$296
☐	187AA13	E0973	T-Post Tall w/ Foam	\$296
☐	187AA3		Flip-Up	NC
☐	187AA1	E0973	Adjustable Locking Flip-Up	\$296
☐	187AA7	E0973	Length Adjust. Locking Flip-Up	\$308
☐	187AM23		Desk Length - 10"	NC
☐	187AM24		Full Length 14"	NC
☐	187AM9		Foam Short (N/A w/ T-Post)	NC
☐	187AM10		Foam Long (N/A w/ T-Post)	NC

Armrest Receiver for Height Adjustable Armrest

☐	187AR17		Low (N/A w/ 20", 22" & 24" wheels)	NC
☐	187AR18		Standard	NC

Side Guard

(N/A with height adjustable armrests or T-Post armrests)

☐	187SG33	K0108	Plastic Std	\$157
☐	187SG29	K0108	Plastic Low	\$157


Armrest - Height Adjustable

☐	187AA11	E0973	Single-Post Height Adj. Std.	\$319
☐	187AM23		Desk Classic	NC
☐	187AM21		Desk Contour	NC
☐	187AM24		Full Classic	NC
☐	187AM22		Full Contour	NC
☐	187AA29	E0973	Single-Post Height Adj. Low	\$319
☐	187AM23		Desk Classic	NC
☐	187AM21		Desk Contour	NC
☐	187AM24		Full Classic	NC
☐	187AM22		Full Contour	NC

CLASSIC PAD						
Armrest Size	Transfer Tube Length (D)	Receiver		Height Range (H)	Pad Length (L)	
					Desk	Full
STD	8"	STD		10-3/4" to 15-3/4"	10-1/4"	14"
		LOW		8-3/4" to 13-3/4"		
LOW	5-1/4"	STD		9" to 13"		
		LOW		7-1/4" to 11-1/4"		

CONTOUR PAD						
Armrest Size	Transfer Tube Length (D)	Receiver		Height Range (H)	Pad Length (L)	
					Desk	Full
STD	8"	STD		10" to 15"	10-1/2"	14-1/4"
		LOW		8-1/4" to 12-1/2"		
LOW	5-1/4"	STD		8-1/2" to 12-1/2"		
		LOW		6-3/4" to 10-3/4"		

The diagram illustrates the physical assembly of the armrest. Dimension 'L' represents the horizontal length of the top pad. Dimension 'H' represents the vertical height from the desk surface to the top of the pad. Dimension 'D' represents the horizontal distance from the vertical support post to the center of the pad's front edge.



Armrest Matrix for Single-Post

9. KEY OPTIONS

Back Pack (Kids Backpack included)

☐	187A1	Utility Backpack - Black	\$84
☐	187A2	Additional Kids Black Backpack	\$67

Seat Pouch

☐	187A8	Black	\$56
--------------------------	-------	-------	------

Miscellaneous

☐	187A83	Tool Kit	\$61
--------------------------	--------	----------	------

Spoke Guard

(N/A with 12" and 16" rear wheels & Projections Handrims)

☐	187A23	K0065 Black	\$128
☐	187A24	K0065 Clear Zippie	\$128

Anti Tips

☐	187AO1	E0971 Front Anti-Tip (N/A with caster pin locks) (N/A with transit option)	\$168
☐	187AT1	E0971 Rear Anti-Tip	\$168

Positioning Belt (w/ clamps)

☐	187A112	E0978 1" Dual Pull Side Squeeze Padded	\$96
☐	187A117	E0978 Belt 2" Aircraft Buckle	\$123
☐	187A118	E0978 Belt 2" Aircraft Buckle Padded	\$133
☐	187A141	E0978 Belt 1.5" Auto Buckle w/ clamps	\$78
☐	187A814	E0978 1" Comfort Fit APB5116	\$133

Canopy (SunShade)

(Need 2" back post clearance to mount canopy hardware, may require next height based on seating; please note some interference may occur with some headrest configurations)

☐	1878A06	K0108 Canopy (Sunshade) 10"-14" Width	\$308
☐	1878A07	K0108 Canopy (Sunshade) 15"-20" Width	\$308

Upper Extremity Support Tray

When ordering a tray with full armrest or desk armrest pads, add 1" to seat width to determine best tray size.

☐	187A29	E0950 X-Sm (Suggested 10"-12" Width)	\$285
☐	187A30	E0950 Sm (Suggested 13"-14" Width)	\$285
☐	187A31	E0950 Med (Suggested 15"-17" Width)	\$285
☐	187A32	E0950 Lrg (Suggested 18"-20" Width)	\$285

10. FRAME COLOR

Standard Colors

☐	187C140	Aztec Gold	NC
☐	187C2	Black (N/A as base w/ Zebra)	NC
☐	187C141	Black Cherry	NC
☐	187C10	Black Opal	NC
☐	187C38	Blue Opal	NC
☐	187C36	Candy Blue	NC
☐	187C14	Candy Purple	NC
☐	187C35	Candy Red	NC
☐	187C142	Electric Blue	NC
☐	187C143	Evergreen	NC
☐	187C126	Glow	NC
☐	187C144	Green Apple	NC
☐	187C147	Hot Sparkle Pink	NC
☐	187C148	Mauve Pink	NC
☐	187C103	Rootbeer	NC
☐	187C114	Sparkle Silver	NC
☐	187C145	Sunrise Orange	NC
☐	187C113	Titanium Color Paint	NC
☐	187C8	Yellow	NC

Standard Matte Colors

☐	187C116	Matte Black	NC
☐	187C151	Matte Black Cherry	NC
☐	187C152	Matte Purple	NC
☐	187C150	Matte Electric Blue	NC
☐	187C149	Matte Evergreen	NC

Personalized Colors

☐	187C43	Mossy Oak Camouflage	\$441
☐	187C42	Stars and Stripes	\$441
☐	187C108	Zebra (Choose base color above: White Std. N/A black, black opal, matte black or titanium)	\$441

Additional Options Through Built-4-Me

To inquire about any unique personalization for Zippie products, contact: Phone: (800) 333-4000, Fax: (800) 333-9011



Commonly Requested Options Include:

- ✓ Dynamic Back
- ✓ Reclining Backrest
- ✓ Personalized Anodizing
- ✓ Personalized Paint
- ✓ O2 Holder with Transit
- ✓ IV Pole with Transit



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA