



# Raz-SP600

## RazBari Self Propel

MSRP – USD – USA  
Effective Jan 03, 2024

Date: \_\_\_\_\_

Dealer / Location: \_\_\_\_\_ T: \_\_\_\_\_ F: \_\_\_\_\_

PO# / Client: \_\_\_\_\_ Contact: \_\_\_\_\_ Email: \_\_\_\_\_

**Raz-SP600** Heavy Duty Mobile Shower Commode Chair with stainless steel frame, Fixed Back, fore / aft-adjustable 24" rear wheels, flip-up arm supports, heavy duty foot supports, front dual-locking 5" casters and commode pan. Also included is a 20"H Multi-Position Back with 2" fore/aft adjustment range, two back angles (90°/97°) and tension-adjustable vinyl upholstery. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 2" of fore/aft seat adjustment. **User weight capacity 605 lbs.** User weight over 450 lbs. requires 6" casters (ZCUSP6). For help with justification, please see our [Evaluation and Medical Justification Form](#).

**Need help with this order form? Call Raz customer service at (877) 720-5678**

### 1. Chair Codes by size – **Must Choose One**

Chair Code	Seat Frame Width	MSRP
Z26022	22"W Frame	\$3,495
Z26024	24"W Frame	\$3,795
Z26026	26"W Frame	\$4,195
Z26028	28"W Frame	\$4,395
Z26030	30"W Frame – Requires custom seat	\$4,995

### 2. Seats – **Must Choose One**

#### Molded and Visco Foam Interface Seats

ZMS23FRO	Molded 23"W x 18"D – Front Access	Fits 22" and 24"W <sup>1</sup> frames only	N/C
ZMS26FRO	Molded 26"W x 18"D – Front Access	Fits 26" and 28"W <sup>1</sup> frames only	N/C
ZVFI23FRO	Visco Foam Interface 23"W x 18"D – Front Access	Fits 22" and 24"W <sup>1</sup> frames only	\$625
ZVFI26FRO	Visco Foam Interface 26"W x 18"D – Front Access	Fits 26" and 28"W <sup>1</sup> frames only	\$745
ZSSKASHD	Seat Support Kit for Raz-AP600 – for seats that are 2" narrower than seat frames		\$195

<sup>1</sup> Required for seats that are narrower than seat frames

#### Custom Seats

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. [See Custom Seat Assessment Form](#) – Resources tab at [razdesigninc.com](http://razdesigninc.com)

ZCSEATHD	Heavy Duty Custom Seat (enter specs below) <b>Select if both width and depth are ≤ 28"</b>	\$1,095
ZCSEATXHD	Heavy Duty Custom Seat for widths 28" and 30" (enter specs below) <b>Select if width and/or depth is ≥ 28"</b>	\$1,295

Seat Frame Width:    \_\_ 22" / 24" / 26" / 28" / 30"  
Custom Seat Width:    \_\_ 22" / 24" / 26" / 28" / 30" / 32" (cannot be more than 2" wider than frame width)  
Seat Depth:            \_\_ 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" [≥ 21" requires custom longer frame]  
Access:                \_\_ FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged  
Aperture Width:       \_\_ 6.5" recommended

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom seat widths can be up to 2" wider than the frame width to provide additional thigh support without increasing overall width of the chair.
- For Custom Seats, the Ischial Pelvic Alignment System (IPAS) allows +/- 1" fore / aft adjustment from specified custom seat depth.
- Access for Rear-Open seats is limited due to rear seat support tube.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, pelvic obliquity accommodation, etc.

**3. Symphony Back** (Recessed solid back that replaces the standard Tension-Adjustable Vinyl Upholstery [recommended for infection-control environments]). All Symphony Back panels are 13"H.

ZSB22	Symphony Back – 22"W	\$465
ZSB24	Symphony Back – 24"W	\$485
ZSB26	Symphony Back – 26"W	\$505
ZSB28	Symphony Back – 28"W	\$575
ZSB30	Symphony Back – 30"W	\$675

#### 4. Rear Wheel, Caster and Anti-Tipper Options

ZRW24HD	24" HD Wheels w/ Ergonomic Handrims	Std	ZATSPRR	Rear Anti-Tippers – telescoping	\$115
ZWLE	Wheel Lock with Extensions (pr)	\$75	ZATAS	Front Anti-Tippers (pr)	\$245
ZVH24	Vinyl-Coated Handrims (pr)	\$425	ZCUSP6	6" Caster Upgrade (2x)	\$175
ZPH24	Vertical Projection Handrims (pr)	\$655			
ZCSLSP	User-Actuated Front Directional Caster Locks w/ extension handles (pr)	\$395			

#### 5. Arm Support Accessories

ZALKHD	Arm Support Lock Kit – HD (pr)	\$95	ZPAML	Pivoting Arm Mount – Left (for ZFAPLL / ZTRL)	\$165
ZFAPLL	Long Molded Flat Armpad <sup>1</sup> – Left (14"L x 3.25"W)	\$90	ZPAMR	Pivoting Arm Mount – Right (for ZFAPLR / ZTRR)	\$165
ZFAPLR	Long Molded Flat Armpad <sup>1</sup> – Right (14"L x 3.25"W)	\$90	ZAPSB_ _ *	Anterior Posture Support Bar * Enter frame width to complete code (22" / 24")	\$245
ZTRL	Molded Arm Trough <sup>1</sup> – Left (18.5"L x 4.5"W)	\$110	ZAPSB_ _ *	Anterior Posture Support Bar * Enter frame width to complete code (26" / 28" / 30")	\$275
ZTRR	Molded Arm Trough <sup>1</sup> – Right (18.5"L x 4.5"W)	\$110			

<sup>1</sup> May require ZPAM for arm support to flip up

#### 6. Foot and Leg Supports

ZFSDH	Adjustable Footrests HD (pr) – 16"-20.5"	Std	ZELRL	Elevating Legrest <sup>1</sup> – Left (incl ZADML)	\$495
ZFET2	Footrest Extension Tube (pr) – 2" longer	\$50	ZELRR	Elevating Legrest <sup>1</sup> – Right (incl ZADMR)	\$495
ZADM	Angle/Depth Adj Footplates (pr)	\$120	ZFOR	Forward Offset Foot Support Receiver - Required when seats extend beyond front edge of frame [provides 2.5" offset] Specify 1 or 2 _ _	\$95
ZFET2A	Footrest Extension Tube (pr) for ZADM – 2" longer	\$50	ZRLSL	Universal Residual Limb Support – Left	\$195
ZFBAS	Flip-Back Foot Support – 1-piece 16"-20.5" length adjustment	\$295	ZRLSR	Universal Residual Limb Support – Right	\$195

<sup>1</sup> ZELR compatible with adjusted seat depths ≤ 20"D

**7. Foot and Leg Support Accessories**

ZFSLK	Foot Support Lock (pr)	\$95	ZNC	Neoprene Covers for Standard Footplates (pr)	\$90
ZCS__*	Snap-On Calf Strap – Bodypoint Aeromesh * Enter frame width to complete code (22" / 24")	\$105	ZNCA	Neoprene Covers for Angle/Depth Adjustable Footplates	\$90
ZCS__*	Snap-On Calf Strap – Bodypoint Aeromesh * Enter frame width to complete code (26" / 28" / 30")	\$245	ZNCFB	Neoprene Cover for Flip-Back Foot Support	\$90
ZCSIC__*	Snap-On Infection Control Calf Strap * Enter frame width to complete code (22" / 24" / 26" / 28" / 30")	\$125	ZFOR	Forward Offset Foot Support Receiver - Required when seats extend beyond front edge of frame [provides 2.5" offset] Specify 1 or 2 ____	\$95
ZHS2224	H-Strap – for 22"/24"W frames – hook and loop attachment	\$110	ZLMOR	Lateral / Medial Offset Foot Support Receiver – moves foot support 2" laterally or medially. Specify 1 or 2 ____	\$95
ZAHL	Adjustable Heel Loops (pr)	\$60			

**8. Positioning Pads** – Pad Sizes: Small 5.5" x 3.5" / Large – 6.75" x 4.3".

<b>Head Supports</b>					
ZHSAS1	Head Support with small pad				\$495
ZHSAS2	Head Support with large pad				\$545
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment				\$75
<b>Lateral Thoracic Supports</b> – Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For up to 4.5" medial offset, specify ZLAEP extension plate.					
ZLTS2LHD	Lateral Thoracic Support w/ Large Pad for HD chairs – Left				\$225
ZLTS2RHD	Lateral Thoracic Support w/ Large Pad for HD chairs – Right				\$225
<b>Leg Adductors</b>					
ZLGAS2L	Leg Adductor w/ Large Pad – Left				\$225
ZLGAS2R	Leg Adductor w/ Large Pad – Right				\$225
<b>Extension Plates for Lateral Thoracic Supports and Hip / Leg Adductors</b>					
These plates allow for an additional 1.5" medially and can be used on Lateral Thoracic Support, Hip and Leg Adductors					
ZLAEP	Lateral/Adductor Extension Plate (ea)				\$28
Lateral Thoracic Support	Left	Right			
Leg Adductor	Left	Right			
Hip Adductor	Left	Right			
Total Quantity ____					
<b>Medial Thigh Supports (Pommel)</b>					
ZMTSAS2	Medial Thigh Support w/ Large Pads				\$595

## 9. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts			Chest Belts		
ZPBXLHDAS	Pelvic Belt – XL (seat frame widths 24"-26", 2-piece w/ side-release buckle)	\$75	ZCBXL	Chest Belt – XL <sup>2</sup> (back frame widths 24"-26", 1-piece w/ hook and loop closure)	\$55
ZPBXXLHDAS	Pelvic Belt – XXL (seat frame widths 28"-30", 2-piece w/ side-release buckle)	\$110	ZCBXXL	Chest Belt – XXL <sup>2</sup> (back frame widths 26"-30", 1-piece w/ hook and loop closure)	\$65
ZPBICXLHDAS	Infection Control Pelvic Belt – XL <sup>1</sup> (seat frame widths 24"-26", 2-piece w/ side-release buckle)	\$125	ZCBICXL	Infection Control Chest Belt – XL <sup>1,2</sup> (back frame widths 24"-26", 1-piece w/ side-release buckle)	\$95
ZPBICXXLHDAS	Infection Control Pelvic Belt – XXL <sup>1</sup> (seat frame widths 26"-30", 2-piece w/ side-release buckle)	\$135	ZCBICXXL	Infection Control Chest Belt – XXL <sup>1,2</sup> (back frame widths 26"-30", 1-piece w/ side-release buckle)	\$105
ZPBBPXLHDAS	Pelvic Belt – Bodypoint Aeromesh – XLarge	\$235	ZCBBPXL	Chest Belt – Bodypoint Aeromesh – XLarge	\$195
Recommended sizes for Bodypoint belts		Chest belt	ZCBBPXL	22"-30"W back frame	
		Pelvic belt	ZPBBPXLHDT	22"-30"W base frame	
<sup>1</sup> Non-absorbent belts – recommended for infection-control environments. <sup>2</sup> Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCB). When the chair is supplied with a Symphony Back, Chest Belt loops are provided.					

## 10. Urine Management

ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$67
ZSS	SplaShield – Molded urine deflector (must be held in place during use)	\$125
ZSC23HD	Splash Curtain <sup>1</sup> for 23"W HD Molded and VFI Seats (surrounds entire aperture, raises seat 1/2").	\$345
ZSC26HD	Splash Curtain <sup>1</sup> for 26"W HD Molded and VFI Seats (surrounds entire aperture, raises seat 1/2").	\$375
ZSCCHD	Custom Splash Curtain for Custom Seats – Raises seat 1/2"	\$365
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$345
ZBUDBH	Bard urine drainage bag holder for 22" / 24" / 26" / 28" / 30"W frames	\$155
<sup>1</sup> Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for pricing.		

## 11. Other

ZPSK1	Commode Pan Spacer Kit <sup>1</sup> – 1" lower	\$200
ZPSK2	Commode Pan Spacer Kit <sup>1</sup> – 2" lower	\$215
ZPSK3	Commode Pan Spacer Kit <sup>1</sup> – 3" lower	\$230
ZPSK4	Commode Pan Spacer Kit <sup>1</sup> – 4" lower	\$245
<sup>1</sup> Commode Pan Spacer Kits provide extra clearance for perianal care with pan in place. Increases height required to clear toilet.		

**12. Credits and Omissions**

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z907HD	Omit HD Foot Supports – if no other Foot / Leg Support Option is selected	-\$45
------	--	-------	--------	---	-------

Some assembly of Raz products is required.

**13. Notes:****Specifications****Raz-SP600 (Self Propel)**

Weight Capacity	605 lbs	
Overall Width	Seat frame width plus 7.75" (with Standard Arm Support Pads)	
Overall Length (with foot supports in middle position)	44" (with rear wheel in middle position)	
Seat Depth	18"-20"	
Clearance Height (under seat)	with 5" casters	with 6" casters
- with pan brackets	14.35"-18.35"	15.35"-19.35"
- without pan brackets	15.12-19.12"	16.12"-20.12"
Arm Support Height - above molded seat	10"	
Seat-to-Floor Height (with molded seat)	with 5" casters	with 6" casters
	20" and 21"	20", 21", 24"
Back Height (Molded Seat to top of upholstery)	20"	
Product Weight (with 22"W seat frame)	68 lbs	
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)	