


Ethos User Activated Anti-Tip



Straight Anti-Tip	
Wheel Size	Rear Seat Height
22"	16.5" - 20"
24" - 25"	17" - 20"
26"	18" - 20"

Bent Anti-Tip	
Wheel Size	Rear Seat Height
22"	13.5" - 16"
24" - 25"	14.5" - 16.5"
26"	15.5" - 17.5"

A pair of User Activated Anti-Tips is not available if a tipping lever is present

Item Number	Part Number	Description	Quantity
1a, 2-5	108098	User Activated Anti-Tip Assembly Straight Single Ethos	1
1a, 2-5	108099	User Activated Anti-Tip Assembly Straight PAIR Ethos	1
1b, 2-5	108100	User Activated Anti-Tip Assembly Bent Single Ethos	1
1b, 2-5	108101	User Activated Anti-Tip Assembly Bent PAIR Ethos	1
1a	004350	Anti-Tip Upper Assembly Straight User Activated Ethos	1
1b	004351	Anti-Tip Upper Assembly Bent User Activated Ethos	1
1, 2, 10, 12	119668	Anti-Tip Lower Assembly	1
2	004306	Anti-Tip Anti Rattle Bushing	1
4	101835	M6 Lock Washer Blk Zn	2
5	101721	M6 X 16 Btn Shcs W/ Patch	2