

Raz-SP Self Propel

MSRP – USA – USD Effective Nov 13, 2023

RAZ	Self Propei		Date:	
Dealer / Location:		T:	F:	
PO# / Client:	Contact:	Email:		

Raz-SP: Self Propel Mobile Shower Commode Chair with height-adjustable stainless-steel frame, Fixed Back, flip-up padded arm supports, foot supports, 22" or 24" rear wheels, front dual-locking 5" (125mm) casters and commode pan. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 2" of fore/aft seat adjustment. User weight capacity 407 lbs.

For help with justification, please see our Evaluation and Medical Justification Form.

Need help with this order form? Call Raz customer service at (877) 720-5678

1. Standard	Raz-SP Frame (18"W) – If you require a special width, depth and / or height see [Sec. 2]	
Z20018	Raz-SP w/ 18"W base frame – Standard base frame width, will fit over most toilets	\$2,350

2. Custom / Sp Z20016	Z20016 Raz-SP w/ 16"W frame width ¹ (base chair price – custom upcharge below will be added)							
Z2SW16	Upcharge for custom 16"W frame – will not fit over most toilets	\$495						
Z20020	Raz-SP w/ 20"W frame width ² (base chair price – custom upcharge below will be added)	\$2,350						
Z2SW20	Upcharge for custom 20"W frame	\$495						
Z20022	Raz-SP w/ 22"W frame width (base chair price – custom upcharge below will be added)	\$2,350						
Z2SW22	Upcharge for custom 22"W frame	\$595						
Custom Frame Co	onfigurations							
ZLNGF3	Upcharge for custom 3" Longer Frame	\$395						
Z2LOWF4	Upcharge for custom Low Frame (16"-18" STF, 4" lower than standard) -22 " rear wheels only	\$495						
¹ 18"W Visco Foan	n Interface seats can be modified to fit.							

²18"W seats can be selected but must also select Seat Support Kit (ZSSKAS) [Sec.3]

3. Seats – Standard – Please note that the Molded and VFI seats, formerly designated 16"D and 19"D, are now designated 17"D and 20"D, which is the center of the IPAS range and the position in which they are shipped from the factory.

	J-				
Molded Seats	(thickness ~2.5"-3")		Visco-Foam Int	erface Seats (thickness ~3.5")	
ZMS17FRO	18"W x 17"D – Front Access	Standard	ZVFI17BRI	18"W x 17"D – Front Bridged	\$425
ZMS17BRI	18"W x 17"D – Front Bridged	N/C	ZVFI17FRO	18"W x 17"D – Front Access	\$425
ZMS20FRO	18"W x 20"D – Front Access	\$135	ZVFI17LTO	18"W x 17"D – Left Access	\$425
ZMS20BRI	18"W x 20"D – Front Bridged	\$135	ZVFI17RTO	18"W x 17"D – Right Access	\$425
Specialty Seat	S		ZVFI17RRO	18"W x 17"D – Rear Access	\$425
ZSOM	Shower-Only Seat – Molded	\$425	ZVFI20BRI	18"W x 20"D – Front Bridged	\$495
	18"W x 16"-21"D		ZVFI20FRO	18"W x 20"D – Front Access	\$495
ZEJRM	E&J Replica Seat – Molded	\$195	ZVFI20LTO	18"W x 20"D – Left Access	\$495
75101/	– 18"W x 17"D	4005	ZVFI20RTO	18"W x 20"D – Right Access	\$495
ZEJRV	E&J Replica Seat – Visco Foam – 18"W x 17"D	\$895	ZVFI20RRO	18"W x 20"D – Rear Access	\$495
Seat Options					
ZSSKASSeat Support Kit for Raz-SP- These are used to support the sides of seats that are narrower than seatSframes on which they are mounted. The reason for a narrower seat on a wider frame is typically to select a Molded Seat and avoid the upcharge of a custom seat. A typical example would be when a 20"W chair is needed but an 18"W standard seat (Molded, VFI, EJR, or Shower Only Seat) is selected.					
ZSS05	Spacers for Front Seat Brackets – 1/	2" (pr) – for ad	ding dump angle	(based on 18"D seat)	\$85
ZSS10	Spacers for Front Seat Brackets – 1"	(pr) – for addi	ng dump angle (b	based on 18"D seat)	\$95

4. Seats – Custom

4. Seals - Cusion								
	ecify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Reappening is not required. See <u>Custom Seat Assessment Form</u> – Resources tab at razdesigninc.com							
ZCSEAT Custom Dual-Visco Foam Seat (enter specs below) Select if both width and depth are ≤ 21" \$799								
ZCEJR	Custom E&J Replica Seat (enter specs below) Select if both width and depth are $\leq 21''$	\$995						
ZCSOV	Custom Dual-Visco Foam Shower-Only Seat (no aperture, no drain hole)	\$595						
ZCSEATHD	Custom Dual-Visco Foam Seat – HD (enter specs below) Select if width and/or depth is $\geq 22''$	\$1,095						
Custom Specifications	(for above custom seats)							
Base Frame Width: $16'' / 18'' / 20'' / 22''$ Custom Seat Width: $16'' / 18'' / 20'' / 22'' / 24'' (cannot be more than 2'' wider than frame width)Custom Seat Depth:15'' / 16'' / 17'' / 18'' / 20'' / 21'' / 22'' / 23'' / 24'' [\geq 21'' requires Long Frame (ZLNGF3)]Access Opening:FR - Front / LT - Left / RT - Right / RR - Rear / BR - BridgedAccess Opening Width:2.5'' / 3.0'' / 3.5'' / 4.0'' / 4.5'' / 5.0'' (3.5''W is standard)Aperture Width:3.0''-8.0'' (6.5'' recommended) - enter "NA" if aperture is not required$								
 Unless otherwise requested 	, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seat	s.						
• Custom Seat width should not be narrower than frame width.								
 Custom Seat can be up to 2" 	wider than seat frame width to provide additional thigh support without ordering wider chair frame.							
 For Custom Seats, the Ischia 	I Pelvic Alignment System (IPAS) allows +/- 1" fore / aft adjustment from specified custom seat depth.							

• Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, pelvic obliquity accomodation, etc.

5. Fixed Back Frames (18"W only)

This is called a "	'Fixed" Back because the arm supports are not adjustable in height. Includes molded vinyl back s	upport, fixed-height					
flip-up padded a	flip-up padded arm supports, back angle adjustment (90°/97°) with cane-style posts. Back height is measured from top of molded						
seat to top of ba	ack support. Hip Adductors [ZHPG] are not compatible with Fixed Back.						
ZFB1818	Fixed Back – 18"W x 18"H (Armrest Lock Kit (ZALK) [Sec.11] is not included.)	Standard					
ZFBU1818	Tension-Adjustable Back Upholstery for Fixed Back	\$185					

6. Adjusta-Back Frames – Back frame widths will be the same as the seat frame widths unless otherwise selected in Sec 8.

Includes height- and angle-adjustable arms with Armrest Lock Kit (ZALK), tension-adjustable back upholstery, flip-up padded arm supports, back angle adjustment (90°/97°). Back heights of 15"/17"/20" include cane-style posts, on which back upholstery can be positioned 2.5" higher than nominal height indicated below. e.g. 20"H can be moved up to 22.5"H. Adjusta-Back – 24"H includes straight posts, which do not allow the back upholstery to be raised.

ZAB 15*	Adjusta-Back ¹ – 15"H * insert back frame width to complete the code (14" / 16" / 18" / 20" / 22")	\$225					
ZAB17*	Adjusta-Back ¹ – 17" H $*$ insert back frame width to complete the code (14" / 16" / 18" / 20" / 22")	\$225					
ZAB 20*	Adjusta-Back ¹ – 20"H * insert back frame width to complete the code (14" / 16" / 18" / 20" / 22")	\$295					
ZAB 24*	Adjusta-Back ¹ – 24"H $*$ insert back frame width to complete the code (16" / 18" / 20" / 22")	\$345					
-	Reclining Backs – 30° Reclining back frame width must be the same as the seat frame width. 2-Point Arms (Z2PLAS / Z2PSAS) are included with Reclining Backs.						
ZRBS 20*	Reclining Adjusta-Back with 20"H back (anti-tippers included) * insert back frame width to complete the code (16" / 18" / 20" / 22")	\$945					
ZRBS 24*	Reclining Adjusta-Back with 24"H back (anti-tippers included) * insert back frame width to complete the code (16" / 18" / 20" / 22")	\$995					
ZDTL	Dual Triggers – both on left side	N/C					
ZDTR	ZDTR Dual Triggers – both on right side						
¹ Dual Tilt Triggers of	option positions both triggers on one side, allowing for single-handed activation.						

7. Custom Back Frame Width Narrower Than Chair Frame Width – Must choose back frame [Sec.7] (configurations below are not available for recliner backs)

(*****								
NOTE: For back frame width variances on recliner frames, please request custom quotation.								
ZNB*	Custom Back Frame 2" Narrower	\$345	Base Frame Width	16″	18″	20″	22″	
	than base frame		Back Frame Width	14" ¹	16" ²	18″ ³	20″4	
	ack frame widths to complete code. fies a 20"W chair frame with ne				1		1	
ZNB*	Custom Back Frame 4" Narrower	\$355	Base Frame Width	18″	20″	22″	24"	0 0
	than base frame – requires custom seat		Back Frame Width	14" ¹	16" ²	18″ ³	20" ⁴	
 * Enter seat and back frame widths to complete code. eg. ZNB2016 specifies a 20"W seat frame with a 16"W back frame 								
¹ Requires Custom Se	eat							•

² Compatible with all 18"W seats but with reduced seat depth range. Also compatible with Custom or modified Visco-Foam Interface seats.

³ Custom Seat recommended. Also compatible with all 18"W seats with Seat Support Kit ZSSKAS [Sec.3].

⁴ Custom Seat required.

8. Back Support Options – These back supports replace the molded vinyl upholstery on a Fixed Back or tension-adjustable upholstery on an Adjusta-Back. Symphony and Harmony Backs are recommended for infection-control environments.

ZSB14	14"W Symphony Back ¹ – recessed solid back support	\$295
ZSB16	16"W Symphony Back ¹ – recessed solid back support	\$295
ZSB18	18"W Symphony Back ¹ – recessed solid back support	\$295
ZSB20	20"W Symphony Back ¹ – recessed solid back support	\$395
ZSB22	22"W Symphony Back ¹ – recessed solid back support	\$465
ZHB	Harmony Back ^{1,2} (contoured solid back – 18"W x 13"H) Fits 18"W back frames only Fits 17" / 20" back heights only	\$295
¹ Does not	fit 15"H backs	

² Not recommended for reclining backs.

9. Rear Wheels, Anti-Tippers and Caster Options – Rear Anti-Tippers Are Strongly Recommended

		· · · · · · · · · · · · · · · · · · ·				1
	ZRW22	22" Diameter Wheels w/ Textured Handrims	Std	ZCSLSP	User-Actuated Caster Swivel Lock w/	\$395
	ZRW24	24" Diameter Wheels w/ Textured Handrims	N/C		extension handles (pr)	
	ZQRSP	Quick-Release Axles (pr)	\$45	ZATSPRR	Anti-Tippers – Rear (pr) standard and	\$115
	ZWLE	Wheel Lock Extension Kit ¹ (pr)	\$75		long frames ²	
	ZVH22	Vinyl-Coated Handrims for 22" (pr)	\$395	ZATSPLWRR	Anti-Tippers – (pr) Rear for low frame ²	\$115
	ZVH24	Vinyl-Coated Handrims for 24" (pr)	\$425	ZATAS	Anti-Tippers – Front (pr) ²	\$245
	ZPH22	Vertical Projection Handrims – 22" (pr)	\$595	ZCUSP4	4" Caster Upgrade (2x)	\$60
	ZPH24	Vertical Projection Handrims – 24" (pr)	\$655	ZCUSP6	6" Caster Upgrade (2x)	\$175
				² Front Anti-Tippe	rs are provideded according to caster diameter	
1	¹ Compatible with old-style aluminum Wheel Locks only			ordered.		

10. Arm Supports and Accessories

	ports and Accessories				
ZARML	Long Arm Supports (pr) – 13"L pad	Std	ZPAML	Pivoting Arm Mount – Left	\$165
ZARMS	Short Arm Supports ¹ (pr) – 10"L pad	N/C		(for ZFPLL/ZTRL)	
ZALK	Arm Support Locks (pr) for Fixed Back (ZFB1818)	\$75	ZPAMR	Pivoting Arm Mount – Right (for ZFPLR / ZTRR)	\$165
Z2PLAS	2-Point Arm Supports ² (pr) with 13"L pads (std with reclining backs)	\$345	ZASK05	Arm Support Spacer Kit – 1/2" per side for all Adjusta-Backs (pr)	\$60
Z2PSAS	2-Point Arm Supports ² (pr)	\$345	ZASK10	Arm Support Spacer Kit – 1" per side for all Adjusta-Backs (pr)	\$115
	with 10"L pads (std with reclining backs)	¢00	ZPHGL	Pivoting Hand Grips (pr) for 13"L arms ⁴	\$395
ZFAPLL	Long Molded Flat Armpad ³ – Left (14"L x 3.25"W – fits ZARM13)	\$90	ZPHGS	Pivoting Hand Grips (pr) for 10"L arms ⁴	\$395
ZFAPLR	Long Molded Flat Armpad ³ – Right (14″L x 3.25″W – fits ZARM13)	\$90	ZAPSB*	Anterior Posture Support Bar ⁵ * Enter frame width to complete code	\$225
ZFAPSL	Short Molded Flat Armpad ³ – Left (10″L x 3″W – fits ZARM10)	\$90	ZAPSB22	(14" / 16" / 18" / 20") Anterior Posture Support Bar ⁵	\$245
ZFAPSR	Short Molded Flat Armpad ³ – Right (10″L x 3″W – fits ZARM10)	\$90		for a 22"W Back Frame	
ZTRL	Molded Arm Trough ³ – Left (18.5″L x 4.5″W)	\$110			
ZTRR	Molded Arm Trough ³ – Right (18.5″L x 4.5″W)	\$110			
¹ Not compatik	ole with ZFAPL / ZFAPS / ZTR / ZPAM				
	h same back width and frame width only			combined with 2 Point Arms requires a custom	
³ May require 2	ZPAM / ZASK05 / ZASK10 to flip up, depending on		upcharge.		
medial / lateral	position of the Flat Arm Pad or Arm Trough		⁵ Not compatible	e with ZFAPL / ZFAPS / ZTR	

11. Foot and Leg Supports

ZFS	Foot Supports 16"-20.5" (pr)	Std	ZELRL	Elevating Leg Support ³ – Left	
ZFS16	Foot Supports for 16"W frames (pr) (16"-20.5")	N/C	ZELRR	(includes ZADML)	\$495
ZFET2	Footrest Extension Tubes ¹ (pr) – 2" longer	\$50	ZELNN	Elevating Leg Support ³ – Right (includes ZADMR)	Ş495
ZMFX8	MFX8 Foot Supports 8"-12" (pr)	\$40	ZFOR	Forward Offset Foot Support Receiver –	\$95
ZMFX12	MFX12 Foot Supports 12"-16.5" (pr)	\$70		(ea) Moves foot support forward 2.5" to	
ZADM	Adj Angle/Depth Molded Footplates ² (pr)	\$120		accomodate deeper seats and reduced knee flexion range. Specify 1 or 2	
ZFET2A	Footrest Extension Tubes (pr)	\$50	ZRLSL	Universal Residual Limb Support – Left	\$195
	for ZADM – 2" longer		ZRLSR	Universal Residual Limb Support – Right	\$195
ZADABS16	Angle/Depth Adj ABS Footplates For 16"W frames (pr)	\$195			
ZVS*	V-Style Foot Support (1-piece) 16"-20.5" * Enter frame width to complete code (16" / 18")	\$295			
ZVS*	V-Style Foot Support (1-piece) 16"-20.5" * Enter frame width to complete code (20" / 22")	\$345			
ZFBAS	Flip-Back Foot Support (1-piece) 16"-20.5"	\$295			
^I Not compatible	Not compatible with ZMFX12 / ZADM / ZVS / ZFBAS			patible with frame widths $\geq 18''$ and adjusted seat	ŀ
² Not compatible	Not compatible with Z10016			or $\leq 21^{"}$ w/ 3" longer seat frames (ZLNGF3)	-

12. Foot and Leg Support Accessories

ZFSLK	Foot Support Lock (pr)	\$95	ZNC	Neoprene Cover for Standard	\$90
ZCS*	Calf Strap, Snap-On – Bodypoint	\$105		Footplate (pr)	
	Aeromesh * Enter frame width to complete code (16" / 18" /20" / 22")		ZNCA	Neoprene Cover for Angle/Depth Adjustable Footplates	\$90
ZCSIC*	* Calf Strap, Snap-On – Infection Control \$95 ZNC16 * Enter frame width to complete code		Neoprene Covers for Narrow Molded Footplates – 16"W frame (pr)	\$90	
	(16" / 18" /20")		ZNCABS16	Neoprene Cover for Narrow Molded ABS	\$90
ZCSIC22	Calf Strap, Snap-On – Infection Control	\$125		Angle/Depth Adjustable Footplates	
	For a 22"W frames		ZNCVS	Neoprene Cover for V-Style Foot Support	\$90
ZCPL	Calf Panel – Bodypoint Aeromesh ¹ Fits frame widths (16" / 18" / 20")	\$225	ZNCFB	Neoprene Cover for Flip-Back Foot Support	\$90
ZHS1820	H-Strap – for 18" / 20"W frames – hook and loop attachment	\$90	ZLMOR	Lateral / Medial Offset Foot Support Receiver (ea) – moves foot support 2"	\$95
ZHS2224	H-Strap – for 22"W frames	\$110		laterally or medially. Specify 1 or 2	
	 hook and loop attachment 		ZFOR	Forward Offset Foot Support Receiver –	\$95
ZAHL	Adjustable Heel Loops (pr) (not available for ZVS / ZFBAS)	\$60		Moves foot support forward 2.5" to accomodate reduced knee flexion. Specify 1 or 2	
Not compatible	e with ZMFX12 / ZFBAS / ZVS				

13. Positioning Pads – Pad Sizes: Small 5.5" x 3.5" / Large – 6.75" x 4.3".

10.1 051001					
Head Suppo	orts				
ZHSAS1	Head Support with small pad	\$495			
ZHSAS2	Head Support with large pad	\$545			
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment				
	racic Supports – Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For up to 4 ZLAEP extension plate.	.5" medial			
ZLTS1L	Lateral Thoracic Support w/ Small Pad – Left	\$225			
ZLTS1R	Lateral Thoracic Support w/ Small Pad – Right	\$225			
ZLTS2L	Lateral Thoracic Support w/ Large Pad – Left	\$225			
ZLTS2R	Lateral Thoracic Support w/ Large Pad – Right	\$225			
Leg Adduct	ors				
ZLGAS1L	Leg Adductor w/ Small Pad – Left	\$225			
ZLGAS1R	Leg Adductor w/ Small Pad – Right	\$225			
ZLGAS2L	Leg Adductor w/ Large Pad – Left	\$225			
ZLGAS2R	Leg Adductor w/ Large Pad – Right	\$225			
Hip Adducto	ors				
ZHP1L	Hip Adductor w/ Small Pad – Left	\$225			
ZHP1R	Hip Adductor w/ Small Pad – Right	\$225			
ZHP2L	Hip Adductor w/ Large Pad – Left	\$225			
ZHP2R	Hip Adductor w/ Large Pad – Right				
Hip Adductor	s For Reclining Backs				
ZHPR1L	Hip Adductor for Reclining Back w/ Small Pad – Left	\$225			
ZHPR1R	Hip Adductor for Reclining Back w/ Small Pad – Right	\$225			
ZHPR2L	Hip Adductor for Reclining Back w/ Large Pad – Left	\$225			
ZHPR2R	Hip Adductor for Reclining Back w/ Large Pad – Right	\$225			
	tes for Lateral / Medial Supports Illow for an additional 1.5" medially and can be used on Lateral Thoracic Support, Hip and Leg Adductors				
ZLAEP	Lateral/Adductor Extension Plate (ea)	\$28			
Lateral Thorac	cic Support Left Right				
Leg Adductor	Left Right				
Hip Adductor	Left Right				
Total Quantity	/ <u> </u>				
Medial Thigh	Supports (Pommel)				
ZMTSAS1	ZMTSAS1 Medial Thigh Support w/ Small Pads				
ZMTSAS2	2 Medial Thigh Support w/ Large Pads				

14. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts						Chest Belts					
ZPBS	Pelvic Belt 2-piece w/	•	ame widths 14 ise buckle)	-16",	\$50	ZCBS		Belt ² – Small (bac ce w/ hook and lo		,	\$45
ZPBM			i (frame widths release buckle		\$50	ZCBM		Belt ² – Medium (-piece w/ hook ar			\$45
ZPBL	Pelvic Belt 2-piece w/	0 (ame widths 18 Ise buckle)	<i>"-</i> 22 <i>"</i> ,	\$50	ZCBL		Belt ² – Large (bac ce w/ hook and lo			\$45
ZPBXL		0,	seat frame wid ide-release buo		\$65	ZCBXL	Chest	Belt – XLarge ² (b 24", 1-piece w/ h	ack frame v	vidths	\$55
ZPBICS			on Control, Sma	ll -	\$75		closure)				
	(frame wid w/ side-rele	ease buck	le)			ZCBICS		Belt ^{1,2} – Infectior ce w/ side-release			\$75
ZPBICM	Pelvic Belt ¹ – Infection Control, Medium (frame widths 16"-18", 2-piece		dium	\$75		widths 14-16")					
	(frame wid w/ side-rele		· ·			ZCBICM		Belt ^{1,2} – Infectior um 1-piece w/ sid	,	uckle	\$75
ZPBICL	Pelvic Belt ¹ – Infection Control, Large		\$85		(back widths 16"-18")						
	(frame wid w/ side-rele		•			ZCBICL	1-pied	Belt ^{1,2} – Infectior ce w/ side-release		0	\$85
ZPBICXL					\$95			s 18"-22")	•		
	(seat frame w/ side-rele		2"-24", 2-piece le)			ZCBICXL	1-piece w/ side-re				\$95
ZPBBPS	Pelvic Belt	– Bodypoi	nt Aeromesh, S	Small	\$245	700000		e widths 22"/ 24")		c	4075
ZPBBPM	Pelvic Belt	– Bodypoi	nt Aeromesh, N	Ved	\$215	ZCBBPS		Belt – Bodypoint			\$275
ZPBBPL	Pelvic Belt	– Bodypoi	nt Aeromesh, L	arge	\$215	ZCBBPM		Belt – Bodypoint			\$245
ZPBBPXL	Pelvic Belt	– Bodypoi	nt Aeromesh, >	KLarge	\$175	ZCBBPL		Belt – Bodypoint		-	\$245
						ZCBBPXL	Chest	Belt – Bodypoint	Aeromesh,	XLarge	\$195
Recommended sizes for	Chest belt ZCBBPS 14" ZCBBPM 16"			W e frame	ZPBBPL	18"/20"W base frame	ZPBBPXL	22"/24"\ base frar			
Bodypoint belts			W k frame	ZCBBPL	18"/20"W back frame	ZCBBPXL	22"/24"\ back frar				

² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM / ZCBL / ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.

15. Urine Management

ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$67
ZSS	SplaShield – Molded urine deflector (must be held in place during use)	\$125
ZSC17	Splash Curtain ¹ for 17"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	\$325
ZSC20	Splash Curtain ¹ for 20"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	\$345
ZSCFR	Splash Curtain for 17" / 20"D (at front of aperture only, not for Front Open Seats)	\$345
ZSCC	Custom Splash Curtain for Custom Seats – Raises seat 1/2"	\$345
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$345
ZBUDBH	Bard urine drainage bag holder (available for all frame widths)	\$155
¹ Bridged Front S	eats only. For Left, Right, Front, and Rear Open, call for pricing.	

16. Other

ZPSK1	Commode Pan Spacer Kit ¹ – 1" lower	\$200			
ZPSK2	Commode Pan Spacer Kit ¹ – 2" lower	\$215			
ZPSK3	Commode Pan Spacer Kit ¹ – 3" lower	\$230			
ZPSK4	Commode Pan Spacer Kit ¹ – 4" lower	\$245			
¹ Commode P	¹ Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Increases height required to clear toilet.				

17. Credits and Omissions

Z905	Omit Commode Pan and	-\$20	Z907	Omit Foot Supports - if no other Foot /	-\$45
	Mounting Brackets			Leg Support Option is selected	

Some assembly of Raz products is required.

18. Notes:

Specifications

Raz-SP (Self Propel)

	22" wheels	24" wheels			
Weight Capacity	407 lbs				
Overall Width (with Standard Arm Support Pads)	Seat width plus 7.75"				
Overall Length (foot supports in middle position)	41"	42"			
Seat Depth	16" or 21" depending of	on seat/cushion length			
Clearance Height (under seat) - with pan brackets - without pan brackets	15.25" / 16.25" / 17.25 / 19.25" 16" / 17" / 18" / 20"	15.25" / 16.25" / 17.25 / 18.25" 16" / 17" / 18" / 19"			
Clearance Width (to fit around toilet) – 18"W frame	16"				
Arm Support Height - above molded seat (1/2" incre	ments)				
 - 18"H Fixed Back - standard - 15"H Adjusta-Back - 17"H Adjusta-Back - 20"H Adjusta-Back - 24"H Adjusta-Back - 20"H Reclining Adjusta-Back - 24"H Reclining Adjusta-Back 	8.4 5.5"- 6.5"- 7.5"- 7.5"- 7.5"- 7.5"-	.8.5" 9.5" 10.5" 10.5" 10.5"			
Seat-to-Floor Height (with molded seat)	20" / 21" / 22" / 24"	20" / 21" / 22" / 23"			
Product Weight (with 18"W seat)	47 lbs	49 lbs			
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)				