

A vibrant playroom with colorful toys and furniture. In the background, there's a wooden cabinet with colorful drawers (red, green, yellow) and a shelf with various toys like blocks, a doll, and a stacking ring. A table with colorful chairs is also visible. The floor is light-colored. The word "pango" is written in a large, dark blue font on a white background in the foreground.

pango



718.943.6452



www.circlespecialty.com



info@circlespecialty.com



pango

for Daily Living

The Pango is our new activity classroom chair for children with mild to moderate positioning needs. Its attractive and intelligent design features a cushioned seat and back, angle adjustable back, height and depth adjustable seat, and pelvic belt. In this iteration of Circle Specialty equipment, we included the unique feature of seat to floor angle adjustability.

We understand how posture and stability increase your child's confidence and willingness to engage in daily living activities. That is why we included a variety of accessories and customizations **to make this chair truly yours.**



OPTIONAL ACCESSORIES

01 Stabilizer & Rocker

PA2301, PA2302

Stabilizer anchors the chair, preventing it from moving. Easily converts to Rocker, encouraging dynamic movement.



02 Wheels Kit for Wheeled Mobility Base

PA2310

Four casters convert Stabilizer to a Wheeled Mobility Base.



03 Headrest

PA2100

Adjustable height, depth, and angle for posterior head and neck support.



04 Trunk Harness

PA2201, PA2202

Provides child with anterior trunk support.



05 Laterals

PA2110

Height & width adjustable to stabilize and align child's trunk.



06 Abductor

PA2212

Limits forward sliding, prevents legs from scissoring, and is depth adjustable.



07 Hip Guides

PA2211

Provides hip and thigh support.



OPTIONAL ACCESSORIES

08 Footrest

PA2401, PA2402

Height & angle adjustable.



09 Foot & Ankle Positioners

PA2410

Adjustable heel & toe loops.



10 Tray

PA2213

Tray fits both size chairs.



11 Guide Bar

PA2501, PA2502

Bar is height adjustable.



SPECIFICATIONS

PA1200



Small

Seat Width	12"
Seat Depth	9" - 12"
Back Height	12"
Armrest Height	1"-5"
Back Angle	90°, 95°, 100°
Seat Angle	85°, 90°, 95°
Back Width	11"
Seat to Floor	8", 9", 10" 11"

PA1400



Medium

Seat Width	14"
Seat Depth	12.5" - 16.5"
Back Height	16"
Armrest Height	1"-7"
Back Angle	90°, 95°, 100°
Seat Angle	85°, 90°, 95°
Back Width	13"
Seat to Floor	13", 14", 15" 16"