


Focus Spoke Wheel / Tire / Handrim Kits



Not all components of the Wheel, Tire and Handrim kits will be listed. These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed. For Focus chairs, new wheel locks might be needed when changing wheel size or style.

Item Number	Part Number	Description	Quantity
	118974	24" Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, 6 Rivet Nut, Standard Mount	1
	118975	24" Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, 6 Rivet Nut, Close Mount	1
	118756	24" Spoke, Pneumatic Tire (Black), Plastic Coated Handrim Kit, 6 Rivet Nut, Standard Mount	1
	118757	24" Spoke, Pneumatic Tire (Black), Plastic Coated Handrim Kit, 6 Rivet Nut, Close Mount	1
	118760	24" Spoke, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount	1
	118761	24" Spoke, Pneumatic Tire (Black), Projection Handrim Kit, Close Mount	1
	118978	24" Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118979	24" Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118980	24" Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit w/ Micro Grip, Standard Mount	1
	118981	24" Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit w/ Micro Grip, Close Mount	1
	118808	24" Spoke, Full Poly Tire, Plastic Coated Handrim Kit, 6 Rivet Nut, Standard Mount	1
	118809	24" Spoke, Full Poly Tire, Plastic Coated Handrim Kit, 6 Rivet Nut, Close Mount	1
	118812	24" Spoke, Full Poly Tire, Projection Handrim Kit, Standard Mount	1
	118813	24" Spoke, Full Poly Tire, Projection Handrim Kit, Close Mount	1