

Invacare® XPO₂™ Portable Oxygen Concentrator

The new Invacare® XPO₂™ portable concentrator is the latest addition to Invacare's proven line of home oxygen technologies. Now providers and patients have a choice of ambulatory systems from Invacare, the world leader in oxygen. Named for its extreme portability, the XPO₂ is designed to be lightweight, clinically robust, and easy to operate. The supplemental battery and AC/DC power options make it easy for patients to go just about anywhere without the fear of running out of oxygen.



The patented, integrated electronic conserver provides a substantial bolus of oxygen at each setting and breath rate.

Setting	Average Pulse Volume @20 bpm
1	15 ml
2	23 ml
3	31 ml
4	37 ml
5	42 ml

Small, lightweight design to enhance patient convenience.



Features

- Pulse dose settings 1-5 meet the needs of a broad range of patients
- Battery, AC, or DC operation for patient convenience
- No separate battery recharger. Charge internal and external batteries with either the AC or DC adapter
- Supplemental battery pack allows for extended battery duration

Provider Benefits

- Durable, reliable product backed by Invacare's proven history of oxygen technologies
- Clinically robust design to keep patients oxygenated at all levels of activity
- A great alternative or adjunct to HomeFill® to meet patients' ambulatory needs
- Supports a non-delivery model for lower operating costs

Patient Benefits

- Small, lightweight and inconspicuous
- Designed for ease of use
- No more fear about running out of oxygen
- Encourages mobility, travel and independence
- Easy to understand control panel for simple operation

Invacare® XPO₂™ Portable Oxygen Concentrator



Convenient DC adapter and battery pack allow patients the freedom to travel whenever and wherever they want.



The Invacare® XPO₂™ Portable Concentrator includes the components shown. Other accessories are available separately.

Model Selection Guide

Model number	Description
XPO100B	Concentrator, carrying case, supplemental battery, AC adapter, DC adapter
XPO100	Concentrator, carrying case, AC adapter, and DC adapter
XPO110	Supplemental battery pack
XPO120	Wheeled cart
XPO130	AC power adapter
XPO140	DC power adapter
XPO150	Concentrator carrying case
XPO160	Accessory bag



Wheeled cart and accessory bag.

Product Specifications

Oxygen flow	Pulse Dose Delivery, Settings 1-5
Oxygen concentration	90% +/- 3% at all settings
Weight	6 pound concentrator, 1.3 pounds supplemental battery
Dimensions	10" H x 7" W x 4" D
Battery duration	2.5 hours at setting 2, internal battery 2.5 hours at setting 2, supplemental battery
Battery recharge time	3 hours*
Power consumption	35 watts
AC power adapter	100-240 VAC 50/60 Hz, 1.5 amps at 120 VAC
DC power adapter	12-16 VDC, 3.3 amps at 18 VDC
Operating altitude	Up to 10,000 ft (3046 m) above sea level
Oxygen capacity	900 ml per minute
Concentrator limited warranty	5 years
Battery limited warranty	1 year

* recharge time increases if unit is running while charging