



# Raz-CAT

## Compact Attendant Tilt

MSRP – USA – USD  
Effective May 1, 2025

Date: \_\_\_\_\_

Dealer / Location: \_\_\_\_\_ T: \_\_\_\_\_ F: \_\_\_\_\_

PO# / Client: \_\_\_\_\_ Contact: \_\_\_\_\_ Email: \_\_\_\_\_

**Raz-CAT –Compact Attendant Tilt** Mobile Shower Commode Chair with 40° tilt range, head support fixture with small pad, height adjustable stainless steel frame, Fixed Back, flip-up padded arm supports, foot supports, four dual-locking 5" (125mm) casters and commode pan. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 2" of fore/aft seat adjustment. **User weight capacity 300 lbs.**

For help with justification, please see our [Evaluation and Medical Justification Form](#).

**Need help with this order form? Call Raz customer service at (877) 720-5678**

### 1. Standard Raz-CAT Frame (18"W) – If you require a special width, depth and / or height see [Sec. 2]


Z30318	Raz-CAT w/ 18"W base frame – Standard base frame width, will fit over most toilets	\$3,395
--------	--	---------

### 2. Custom / Special Frame Widths, Depths and Heights – Adjusta-Back must be selected for all custom widths see [Sec.7]

Z30316	Raz-CAT w/ 16"W frame width <sup>1</sup> (base chair price – custom upcharge below will be added)	\$3,395
Z3CSW16	Upcharge for custom 16"W seat and base frame – will not fit over most toilets	\$595
Z30320	Raz-CAT w/ 20"W frame width <sup>2</sup> (base chair price – custom upcharge below will be added)	\$3,395
Z3CSW20	Upcharge for custom 20"W seat and base frame	\$595
<b>Custom Frame Configuration for 18"W frame. (Call Raz for codes and pricing for 16" / 20" widths)</b>		
Z3CLOWF5	Upcharge for custom Low Frame (16"-20" STF)	\$1,595

<sup>1</sup> 18"W Visco-Foam Interface seats can be modified to fit.  
<sup>2</sup> 18"W seats can be selected but must also select Seat Support Kit (Z4SSKCT) [Sec.5]

### 3. Seats – Standard – Please note that the Molded and VFI seats, formerly designated 16"D, are now designated 17"D, which is the center of the IPAS range and the position in which they are shipped from the factory.

Standard CAT Seat – Visco-Foam w/ Tapered Oblong Aperture (thickness ~2.5")							
ZCSC1816	18"W x 11"-16"D (adjustable)	standard w/ 18"W frame				N/C	
ZCSC1616	16"W x 11"-16"D (adjustable)	standard w/ 16"W frame				N/C	
ZCSC2016	20"W x 11"-16"D (adjustable)	standard w/ 20"W frame				N/C	
Set seat depth at: 11" 12" 13" 14" 15" 16"							
Standard aperture width: 3.5" at front, 4.5" at rear							
Other Seats Available – for seat depth adjustment up to 17" – Modified to fit CAT seat frame							
Molded Seats (thickness ~2.5"-3")				Visco Foam Interface Seats (thickness ~3.5")			
ZMSC17FRO	18"W x 17"D – Front Access		N/C	ZVFIC17FRO	18"W x 17"D – Front Access	\$425	
ZMSC17BRI	18"W x 17"D – Front Bridged		N/C	ZVFIC17BRI	18"W x 17"D – Front Bridged	\$425	
				ZVFIC17LTO	18"W x 17"D – Left Access	\$425	
				ZVFIC17RTO	18"W x 17"D – Right Access	\$425	
				ZVFIC17RRO	18"W x 17"D – Rear Access	\$425	
				ZVFIC17FAR	18"W x 17"D – Front and Rear Access	\$425	

#### 4. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. 18" is the standard width for CAT seat frames.

See [Custom Seat Assessment Form](#) – Resources tab at [razdesigninc.com](http://razdesigninc.com)

ZCSEAT	Custom Dual-Visco Foam Seat (enter specs below)	\$795
Seat Frame Width:	__ 16" / <b>18"</b> / 20"	
Custom Seat Width:	__ 16" / 17" / 18" / 19" / 20"	
Custom Seat Depth:	__ 11" / 12" / 13" / 14" / 15" / 16" / 17"	
Access Opening:	__ FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged	
Aperture:	Tapered Oblong: __"W (at front) __"W (at rear) [standard is 3.5"W at front and 4.5"W at rear] Keyhole: __"W (3.0"-8.0"W – enter "NA" if aperture is not required)	
<ul style="list-style-type: none"> <li>• Oblong-shaped apertures can be specified in conjunction with any of the access openings listed.</li> <li>• For Custom Seats, the Ischial Pelvic Alignment System (IPAS) allows +/- 1" fore / aft adjustment from specified custom seat depth.</li> <li>• Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, pelvic obliquity accommodation, etc.</li> </ul>		

#### 5. Back Frames – for standard 18"W seat frames (back heights measured from top of seat to top of upholstery)

**Multi-Position Adjusta-Back:** The "Multi-Position" feature allows for fore/aft adjustment of the backrest to provide seat depth adjustability from 11"-16". Includes tension-adjustable vinyl back support, flip-up padded arm supports, back angle adjustment (90° / 97°) with cane-style back posts. Back height is measured from top of molded seat to top of back support. Hip Adductors [ZHPG] are not compatible with Fixed Back.

Back frames for standard 18"W seat frame (Z30318)			Reclining Backs for 18"W seat frame		
ZCAB845	14" Back Width / 15" Back Height	N/C	ZRBC848	Recliner Back for CAT – 14"W x 18"H <sup>2,3</sup>	\$1,095
ZCAB848	14" Back Width / 18" Back Height	\$95	ZRBC868	Recliner Back for CAT – 16"W x 18"H <sup>2,3</sup>	\$1,095
ZCAB865	16" Back Width / 15" Back Height	N/C	ZRBCR888	Recliner Back for CAT – 18"W x 18"H <sup>2,3</sup>	\$1,095
ZCAB868	16" Back Width / 18" Back Height	\$95			
ZCAB885	18" Back Width / 15" Back Height	N/C			
ZCAB888	18" Back Width / 18" Back Height	\$95			
ZDTL	Dual Triggers – left side <sup>1</sup>	N/C			
ZDTR	Dual Triggers – right side <sup>1</sup>	N/C			
<sup>1</sup> Standard trigger configuration for tilt-only is one on each side. Dual Triggers option positions both triggers on one side, allowing for single-handed operation (tilt only).			<sup>2</sup> Standard trigger locations for Raz-CAT with recline are: tilt triggers on the right and recline triggers on the left.		
			<sup>3</sup> Tilt is limited to 25°		

## 6. Back Frames for 16"W and 20"W Seat Frames (back heights measured from top of seat to top of upholstery)

<b>Multi-Position Adjusta-Back:</b> The "Multi-Position" feature allows for fore/aft adjustment of the backrest to provide seat depth adjustability from 11"-16". Includes tension-adjustable vinyl back support, back angle adjustment (90° / 97°) and flip-up padded arm supports. Back upholstery can be positioned 2.5" higher than factory settings stated below.			
<b>16"W Seat Frames</b>			<b>Reclining Backs for 16"W seat frame</b>
ZCAB645	14" Back Width / 15" Back Height	N/C	2-Point Arms (Z2PLT / Z2PST) are included with Reclining Backs.
ZCAB648	14" Back Width / 18" Back Height	\$95	ZRBC648 Recliner Back for CAT – 14"W x 18"H <sup>2,3</sup> \$1,095
ZCAB665	16" Back Width / 15" Back Height	N/C	ZRBC668 Recliner Back for CAT – 16"W x 18"H <sup>2,3</sup> \$1,095
ZCAB668	16" Back Width / 18" Back Height	\$95	
<b>20"W Seat Frame</b>			<b>Reclining Backs for 20"W seat frame</b>
ZCAB065	16" Back Width / 15" Back Height	N/C	2-Point Arms (Z2PLT / Z2PST) are included with Reclining Backs.
ZCAB068	16" Back Width / 18" Back Height	\$95	ZRBC088 Recliner Back for CAT – 18"W x 18"H <sup>2,3</sup> \$1,095
ZCAB085	18" Back Width / 15" Back Height	N/C	
ZCAB088	18" Back Width / 18" Back Height	\$95	
ZCAB005	20" Back Width / 15" Back Height	N/C	
ZCAB008	20" Back Width / 18" Back Height	\$95	
ZDTL	Dual Triggers – left side <sup>1</sup>	N/C	
ZDTR	Dual Triggers – right side <sup>1</sup>	N/C	
<sup>1</sup> Standard trigger configuration for tilt-only is one on each side. Dual Triggers option positions both triggers on one side, allowing for single-handed operation (tilt only).			<sup>2</sup> Standard trigger locations for Raz-CAT with recline are tilt triggers on the right and recline triggers on left. <sup>3</sup> Tilt is limited to 25°

## 7. Back Support Options – These back supports replace tension-adjustable upholstery on an Adjusta-Back. Symphony Backs are recommended for infection-control environments (for 18"H backs only).

ZSB14	Symphony Back – recessed solid back support – 14" back frame width	\$295
ZSB16	Symphony Back – recessed solid back support – 16" back frame width	\$295
ZSB18	Symphony Back – recessed solid back support – 18" back frame width	\$295
ZSB20	Symphony Back – recessed solid back support – 20" back frame width	\$395

## 8. Arm Supports and Accessories

ZARMTS	Short Arm Support Tubes (pr)	Standard
ZARMTL	Long Arm Support Tubes (pr)	N/C
ZARMSL	Short Arm Support Pad <sup>1</sup> 10"L – Left	Standard
ZARMSR	Short Arm Support <sup>1</sup> Pad 10"L – Right	Standard
ZARMLL	Long Arm Support Pad 13"L – Left	N/C
ZARMLR	Long Arm Support Pad 13"L – Right	N/C
ZALK	Arm Support Locks (pr) for Fixed Back (ZFB1818)	\$75
Z2PLT	2-Point Arm Supports <sup>2</sup> (pr) with 13"L pads (std with reclining backs)	\$345
Z2PST	2-Point Arm Supports <sup>2</sup> (pr) with 10"L pads (std with reclining backs)	\$345
ZFAPLL	Long Molded Flat Armpad <sup>3</sup> – Left (14"L x 3.25"W – fits ZARM13)	\$90
ZFAPLR	Long Molded Flat Armpad <sup>3</sup> – Right (14"L x 3.25"W – fits ZARM13)	\$90
ZFAPSL	Short Molded Flat Armpad <sup>3</sup> – Left (10"L x 3"W – fits ZARM10)	\$90
ZFAPSR	Short Molded Flat Armpad <sup>3</sup> – Right (10"L x 3"W – fits ZARM10)	\$90
ZTRL	Molded Arm Trough <sup>3</sup> – Left (18.5"L x 4.5"W)	\$110
ZTRR	Molded Arm Trough <sup>3</sup> – Right (18.5"L x 4.5"W)	\$110
ZPAML	Pivoting Arm Mount – Left (for ZFAPLL/ZTRL)	\$165
ZPAMR	Pivoting Arm Mount – Right (for ZFAPLR / ZTRR)	\$165
ZASK05	Arm Support Spacer Kit – 1/2" per side for all Adjusta-Backs (pr)	\$60
ZASK10	Arm Support Spacer Kit – 1" per side for all Adjusta-Backs (pr)	\$115
ZPHGL	Pivoting Hand Grips (pr) for 13"L arms <sup>4</sup>	\$395
ZPHGS	Pivoting Hand Grips (pr) for 10"L arms <sup>4</sup>	\$395
ZAPSB_ _ *	Anterior Posture Support Bar <sup>5</sup> * Enter frame width to complete code (14" / 16" / 18" / 20")	\$225
ZAPSB_ _ *	Anterior Posture Support Bar <sup>5,6</sup> * Enter frame width to complete code (22" / 24")	\$245
<b>Reclining Back Arm Support Accessories</b>		
Z2PLTR	2-Point Flip-up Arm Supports (Long) (pr)	N/C
Z2PSTR	2-Point Flip-up Arm Supports (Short) (pr)	N/C
ZASKR05	Arm Support Spacers 1/2" per side (pr) for Reclining Back Support	\$60
ZASKR10	Arm Support Spacers 1" per side (pr) for Reclining Back Support	\$115
<sup>1</sup> Not compatible with ZFAPL / ZFAPS / ZTR <sup>2</sup> For chairs with same back width and frame width only <sup>3</sup> May require ZPAM / ZASK05 / ZASK10 for Arm Support to flip up, depending on medial / lateral position of the Flat Arm Pad or Arm Trough <sup>4</sup> ZPHGL ZPHGS combined with 2 Point Arms requires a custom up charge <sup>5</sup> Not compatible with ZFAPL / ZFAPS / ZTR <sup>6</sup> Special length required with ZASK05 (+1") and ZASK (+2")		

## 9. Foot and Leg Supports

ZFSL	Foot Support Left – 16"-20.5"	N/C
ZFSR	Foot Support Right – 16"-20.5"	N/C
ZFET2L	Footrest Extension Tube <sup>1</sup> Left – 2" longer	\$25
ZFET2R	Footrest Extension Tube <sup>1</sup> Right – 2" longer	\$25
ZMX8CL	MX8 Foot Support Left – 8"- 12"	N/C
ZMX8CR	MX8 Foot Support Right – 8"- 12"	N/C
ZMX12CL	MX12 Foot Support Left – 12"-16.5"	Standard
ZMX12CR	MX12 Foot Support Right – 12"-16.5"	Standard
ZADABS18L	Angle/Depth Adj ABS Footplate Left For 18"W frame	\$98
ZADABS18R	Angle/Depth Adj ABS Footplate Right For 18"W frame	\$98
ZADML	Adj Angle/Depth Molded Footplate <sup>2</sup> - Left	\$60
ZADMR	Adj Angle/Depth Molded Footplate <sup>2</sup> - Right	\$60
ZFET2AL	Footrest Extension Tube Left for ZADML – 2" longer	\$25
ZFET2AR	Footrest Extension Tube Right for ZADMR – 2" longer	\$25
ZVS_ _ *	V-Style Foot Support (1-piece) 16"-20.5" * Enter frame width to complete code (16" / 18")	\$295
ZVS20	V-Style Foot Support (1-piece) 16"-20.5" for 20"W frame	\$345
ZFBT12	MX12 Flip-Back Foot Support (1-piece) 12"-16.5" with calf pads	\$455
ZFBT	Flip-Back Foot Support (1-piece) 16.5"-20" with calf pads	\$455
ZELRL	Elevating Leg Support <sup>3</sup> – Left (includes ZADML)	\$495
ZELRR	Elevating Leg Support <sup>3</sup> – Right (includes ZADMR)	\$495
ZRLSL	Universal Residual Limb Support – Left	\$195
ZRLSR	Universal Residual Limb Support – Right	\$195
<b>Foot Support Options for 16"W Frames</b>		
ZFS16L	Foot Support for 16"W frame Left – (16"-20.5")	N/C
ZFS16R	Foot Support for 16"W frame Right – (16"-20.5")	N/C
ZADABS16L	Angle/Depth Adj ABS Footplate Left for 16"W frame	\$98
ZADABS16R	Angle/Depth Adj ABS Footplate Right For 16"W frame	\$98
ZMX816L	MX8 Foot Support for 16"W Frame Left – 8"- 12"	\$20
ZMX816R	MX8 Foot Supports for 16"W Frame Right – 8"- 12"	\$20
ZMX1216L	MX12 Foot Support for 16"W Frame Left – 12" - 16.5"	\$35
ZMX1216R	MX12 Foot Support for 16"W Frame Right – 12" - 16.5"	\$35
<sup>1</sup> Not compatible with ZMX12/ZADM/ZVS/ZFBT		
<sup>2</sup> Not compatible with Z30016		
<sup>3</sup> ZELRL/R compatible with frame widths ≥ 18" and adjusted seat depths ≤ 18" or ≤ 21" w/ 3" longer seat frames (Z3LNGSF3)		

**10. Foot and Leg Support Accessories**

ZFSLK	Foot Support Lock (pr)	\$95
ZCS_*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter frame width to complete code (16" / 18" / 20")	\$105
ZCSIC_*	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code (16" / 18" / 20")	\$95
ZCPL	Calf Panel – Bodypoint Aeromesh <sup>1</sup> Fits frame widths (16" / 18" / 20")	\$225
ZHS1820	H-Strap – for 18" / 20"W frames – hook and loop attachment	\$90
ZAHL	Adjustable Heel Loops (pr) (not available for ZVS / ZFBT)	\$60
ZNC	Neoprene Cover for Standard Footplate (pr)	\$90
ZNCA	Neoprene Cover for Angle/Depth Adjustable Footplates	\$90
ZNC16	Neoprene Covers for Narrow Molded Footplates – 16"W frame (pr)	\$90
ZNCABS16	Neoprene Cover for Narrow Molded ABS Angle/Depth Adjustable Footplates	\$90
ZNCVS	Neoprene Cover for V-Style Foot Support	\$90
ZNCFB	Neoprene Cover for Flip-Back Foot Support	\$90
ZLOR	Lateral Offset Foot Support Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 ____	\$95
ZMOR	Medial Offset Foot Support Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 ____	\$95
ZFOR	Forward Offset Foot Support Receiver – (ea) Moves foot support forward 2.5" to accomodate deeper seats and reduced knee flexion range. Specify 1 or 2 ____	\$95
<sup>1</sup> Not compatible with ZMFX12 / ZFBT / ZVS		

**11. Positioning Pads – Hip Adductors are not compatible with Reclining Backrest unless ordered on-chair.**

<b>Head Supports</b>			
ZHSP1	Head Support Pad – Small (10" x 4")		Standard
ZHSP2	Head Support Pad – Large (12" x 6")		\$50
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment		\$75
<b>Lateral Thoracic Supports</b> – Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For up to 4.5" medial offset, specify ZLAEP extension plate.			
ZLTS1L	Lateral Thoracic Support w/ Small Pad – Left		\$225
ZLTS1R	Lateral Thoracic Support w/ Small Pad – Right		\$225
ZLTS2L	Lateral Thoracic Support w/ Large Pad – Left		\$225
ZLTS2R	Lateral Thoracic Support w/ Large Pad – Right		\$225
<b>Leg Adductors</b>			
ZLGT1L	Leg Adductor w/ Small Pad – Left		\$225
ZLGT1R	Leg Adductor w/ Small Pad – Right		\$225
ZLGT2L	Leg Adductor w/ Large Pad – Left		\$225
ZLGT2R	Leg Adductor w/ Large Pad – Right		\$225
<b>Hip Adductors for 4" narrower backs</b>			
ZHP4C1L	Hip Adductor for 4" narrower back w/ Small Pad– Left		\$225
ZHP4C1R	Hip Adductor for 4" narrower back w/ Small Pad– Right		\$225
ZHP4C2L	Hip Adductor for 4" narrower back Large Pad– Left		\$225
ZHP4C2R	Hip Adductor for 4" narrower back w/ Large Pad– Right		\$225
<b>Hip Adductors for backs that are same width as seat frame or 2" narrower</b>			
ZHP02C1L	Hip Adductor for 0-2" narrower back w/ Small Pad – Left		\$225
ZHP02C1R	Hip Adductor for 0-2" narrower back w/ Small Pad – Right		\$225
ZHP02C2L	Hip Adductor for 0-2" narrower back w/ Large Pad – Left		\$225
ZHP02C2R	Hip Adductor for 0-2" narrower back w/ Large Pad – Right		\$225
<b>Hip Adductors For Reclining Backs</b>			
ZHPRC1L	Hip Adductor for Reclining Back w/ Small Pad – Left		\$225
ZHPRC1R	Hip Adductor for Reclining Back w/ Small Pad – Right		\$225
ZHPRC2L	Hip Adductor for Reclining Back w/ Large Pad – Left		\$225
ZHPRC2R	Hip Adductor for Reclining Back w/ Large Pad – Right		\$225
<b>Extension Plates for Lateral Thoracic Supports and Hip / Leg Adductors</b>			
These plates allow for an additional 1.5" medially and can be used on Lateral Thoracic Support, Hip and Leg Adductors			
ZLAEPL	Lateral/Adductor Extension Plate – Left		\$28
ZLAEPR	Lateral/Adductor Extension Plate – Right		\$28
Lateral Thoracic Support	Left	Right	
Leg Adductor	Left	Right	
Hip Adductor	Left	Right	
Total Quantity ____			
<b>Medial Thigh Supports (Pommel)</b>			
ZMTST1	Medial Thigh Support w/ Small Pads		\$495
ZMTST2	Medial Thigh Support w/ Large Pads		\$495

**12. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!**

Pelvic Belts				Chest Belts			
ZPBS	Pelvic Belt – Small (frame widths 14-16", 2-piece w/ side-release buckle)		\$50	ZCBS	Chest Belt <sup>2</sup> – Small (back widths 14"-16", 1-piece w/ hook and loop closure)		\$45
ZPBM	Pelvic Belt – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)		\$50	ZCBM	Chest Belt <sup>2</sup> – Medium (back widths 16"-18", 1-piece w/ hook and loop closure)		\$45
ZPBL	Pelvic Belt – Large (frame widths 18"-20", 2-piece w/ side-release buckle)		\$50	ZCBL	Chest Belt <sup>2</sup> – Large (back widths 18"-20", 1-piece w/ hook and loop closure)		\$45
ZPBICS	Pelvic Belt <sup>1</sup> – Infection Control, Small 2-piece w/ frame clamps, mounting plates, side-release buckle (frame widths 14"-16")		\$75	ZCBICS	Chest Belt <sup>1,2</sup> – Infection Control, Small 1-piece w/ side-release buckle (back widths 14-16")		\$75
ZPBICM	Pelvic Belt <sup>1</sup> – Infection Control, Medium 2-piece w/ frame clamps, mounting plates, side-release buckle (frame widths 16"-18")		\$75	ZCBICM	Chest Belt <sup>1,2</sup> – Infection Control, Medium 1-piece w/ side-release buckle (back widths 16"-18")		\$75
ZPBICL	Pelvic Belt <sup>1</sup> – Infection Control, Large 2-piece w/ frame clamps, mounting plates, side-release buckle (frame widths 18"-20")		\$85	ZCBICL	Chest Belt <sup>1,2</sup> – Infection Control, Large 1-piece w/ side-release buckle (back widths 18"-20")		\$85
ZPBBPM	Pelvic Belt – Bodypoint Aeromesh, Med		\$215	ZCBBPM	Chest Belt – Bodypoint Aeromesh, Med		\$215
ZPBBPL	Pelvic Belt – Bodypoint Aeromesh, Large		\$215	ZCBBPL	Chest Belt – Bodypoint Aeromesh, Large		\$215
<b>Recommended sizes for Bodypoint belts</b>	Pelvic belt	ZPBBPS	N/A	ZPBBPM	16"W base frame	ZPBBPL	18"/20"W base frame
	Chest belt	ZCBBPS	14"W back frame	ZCBBPM	16"W back frame	ZCBBPL	18"/20"W back frame
<sup>1</sup> Non-absorbent belts – recommended for infection-control environments <sup>2</sup> Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL/ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.							

**13. Urine Management**

ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$67
ZSS	SplaShield – Molded urine deflector (must be held in place while using)	\$125
ZSC17	Splash Curtain <sup>1</sup> for 17"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	\$325
ZSCC	Custom Splash Curtain for Custom Seats – Raises seat 1/2"	\$345
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$345
ZBUDBH	Bard Urine Drainage Bag Holder	\$155
<sup>1</sup> Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for order codes and pricing.		



**14. Other**

ZDCLT	Directional Caster Lock – fits on one rear caster (facilitates steering)	\$165	ZPSK1	Commode Pan Spacer Kit <sup>1</sup> – 1" lower	\$200
ZATRR	Anti-Tippers – Rear (pr)	\$345	ZPSK2	Commode Pan Spacer Kit <sup>1</sup> – 2" lower	\$215
ZATFR	Anti-Tippers – Front (pr)	\$345	ZPSK3	Commode Pan Spacer Kit <sup>1</sup> – 3" lower	\$230
ZCUT4	4" Caster Upgrade (4x)	\$120	ZPSK4	Commode Pan Spacer Kit <sup>1</sup> – 4" lower	\$245
ZCUT6	6" Caster Upgrade (4x)	\$350	ZSPKC22	CAT Self-Propel Kit with 22" wheels – incl anti-tippers (tilt limited to 25°)	\$1,425
ZTLS	Tilt Limit Spacers – 1/2" spacers snap onto gas struts to limit tilt (each reduces tilt range approx. 5°. Can be stacked for additional tilt limiting. Quantity ____	N/C	ZSPKC24	CAT Self-Propel Kit with 24" wheels – incl anti-tippers (tilt limited to 25°)	\$1,425

<sup>1</sup> Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Increases height required to clear toilet.

**15. Credits and Omissions**

Z905	Omit Commode Pan and Mounting Brackets	–\$20	Z911	Omit Molded Head Support Pad	–\$20
Z907	Omit Foot Supports	–\$45	Z912	Omit Complete Head Support Assembly	–\$50

**Some assembly of Raz products is required.**

**16. Notes****Specifications****Raz-CAT (Compact Attendant Tilt)**

Weight Capacity	300 lbs (136 kg)
Overall Width	Seat width plus 4.5" (with Standard Arm Support Pads)
Overall Length (with MFX12 foot supports and headrest adjusted to middle positions)	
<ul style="list-style-type: none"> <li>• at 0° tilt</li> <li>• at 40°</li> </ul>	36.75" (93.5cm) (rear casters trailing) 44" (112cm)
Seat Depth	11"-16" (28-41cm) depending on seat/cushion length
Clearance Height (under seat)	
<ul style="list-style-type: none"> <li>• with pan brackets</li> <li>• without pan brackets</li> </ul>	17.25"-21.25" (44-54cm) (1" increments) 17.75"-21.75" (45-55cm) (1" increments)
Clearance Width (to fit around toilet)	15.75"(40cm) for standard 18"W (46cm) seat frame. For other widths, seat frame width minus 2.25" (5.75cm)
Arm Support Height - above molded seat (1/2" increments) (1.25cm)	
<ul style="list-style-type: none"> <li>• 15"H Adjusta-Back</li> <li>• 18"H Adjusta-Back</li> <li>• 18"H Reclining Adjusta-Back</li> </ul>	5.5"-8" (14-20cm) 6.5"-9" (16.5-23cm) 6.5"-9" (16.5-23cm)
Seat-to-Floor Height (with molded seat)	21"-25" (53.5-63.5cm)
Product Weight	55 lbs (25 kg)
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)