


Catalyst Half Folding Tension Back Upholstery



| Straight Back | | | | 8° Bend Back | | | |
|---------------------------|--------|-----------------|------------------|---------------------------|--------|-----------------|------------------|
| Backrest Range and Height | Wings | Total Strap Qty | Total D-Ring Qty | Backrest Range and Height | Wings | Total Strap Qty | Total D-Ring Qty |
| Medium 12.5" - 13.5" | 003094 | 4 | 2 | Medium 13.5" - 14.5" | 003094 | 5 | 4 |
| | 003095 | | | | 003095 | | |
| Medium 14.5" - 15.5" | 003094 | 5 | 4 | Medium 15.5" | 003094 | 6 | 6 |
| | 003095 | | | | 003095 | | |
| Tall 16.5" - 17.5" | 003096 | 6 | 2 | Tall 16.5" - 17.5" | 003096 | 6 | 2 |
| | 003097 | | | | 003097 | | |
| Tall 18.5" - 19.5" | 003096 | 7 | 4 | Tall 18.5" - 19.5" | 003096 | 7 | 4 |
| | 003097 | | | | 003097 | | |

Attaching hardware is not included with the back upholstery assemblies

| Item Number | Part Number | Description | Quantity |
|-------------------|-------------|---|----------|
| 1a, 2a, 3a, 4, 5a | 506510 | Back Upholstery Assembly Tension Adjustable with Standard Cover 14W 12.5-13.5H, Half Folding Back Straight | 1 |
| 1b, 2a, 3a, 4, 5a | 506511 | Back Upholstery Assembly Tension Adjustable with Standard Cover 14W 14.5-15.5H, Half Folding Back Straight | 1 |
| 1b, 2a, 3a, 4, 5a | 506512 | Back Upholstery Assembly Tension Adjustable with Standard Cover 14W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 1b, 2a, 3a, 4, 5a | 506513 | Back Upholstery Assembly Tension Adjustable with Standard Cover 14W 15.5H, Half Folding Back 8 Degree Bend | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------------|-------------|---|----------|
| 1c, 2b, 3b, 4, 5a | 506514 | Back Upholstery Assembly Tension Adjustable with Standard Cover 14W 16.5-17.5, Half Folding Back | 1 |
| 1d, 2b, 3b, 4, 5a | 506515 | Back Upholstery Assembly Tension Adjustable with Standard Cover 14W 18.5-19.5, Half Folding Back | 1 |
| 1e, 2a, 3a, 4, 5b | 506516 | Back Upholstery Assembly Tension Adjustable with Standard Cover 15W 12.5-13.5H, Half Folding Back Straight | 1 |
| 1f, 2a, 3a, 4, 5b | 506517 | Back Upholstery Assembly Tension Adjustable with Standard Cover 15W 14.5-15.5H, Half Folding Back Straight | 1 |
| 1f, 2a, 3a, 4, 5b | 506518 | Back Upholstery Assembly Tension Adjustable with Standard Cover 15W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 1f, 2a, 3a, 4, 5b | 506519 | Back Upholstery Assembly Tension Adjustable with Standard Cover 15W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 1g, 2b, 3b, 4, 5b | 506520 | Back Upholstery Assembly Tension Adjustable with Standard Cover 15W 16.5-17.5, Half Folding Back | 1 |
| 1h, 2b, 3b, 4, 5b | 506521 | Back Upholstery Assembly Tension Adjustable with Standard Cover 15W 18.5-19.5, Half Folding Back | 1 |
| 1i, 2a, 3a, 4, 5c | 506522 | Back Upholstery Assembly Tension Adjustable with Standard Cover 16W 12.5-13.5H, Half Folding Back Straight | 1 |
| 1j, 2a, 3a, 4, 5c | 506523 | Back Upholstery Assembly Tension Adjustable with Standard Cover 16W 14.5-15.5H, Half Folding Back Straight | 1 |
| 1j, 2a, 3a, 4, 5c | 506524 | Back Upholstery Assembly Tension Adjustable with Standard Cover 16W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 1j, 2a, 3a, 4, 5c | 506525 | Back Upholstery Assembly Tension Adjustable with Standard Cover 16W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 1k, 2b, 3b, 4, 5c | 506526 | Back Upholstery Assembly Tension Adjustable with Standard Cover 16W 16.5-17.5, Half Folding Back | 1 |
| 1l, 2b, 3b, 4, 5c | 506527 | Back Upholstery Assembly Tension Adjustable with Standard Cover 16W 18.5-19.5, Half Folding Back | 1 |
| 1m, 2a, 3a, 4, 5d | 506528 | Back Upholstery Assembly Tension Adjustable with Standard Cover 17W 12.5-13.5H, Half Folding Back Straight | 1 |
| 1n, 2a, 3a, 4, 5d | 506529 | Back Upholstery Assembly Tension Adjustable with Standard Cover 17W 14.5-15.5H, Half Folding Back Straight | 1 |
| 1n, 2a, 3a, 4, 5d | 506530 | Back Upholstery Assembly Tension Adjustable with Standard Cover 17W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 1n, 2a, 3a, 4, 5d | 506531 | Back Upholstery Assembly Tension Adjustable with Standard Cover 17W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 1o, 2b, 3b, 4, 5d | 506532 | Back Upholstery Assembly Tension Adjustable with Standard Cover 17W 16.5-17.5, Half Folding Back | 1 |

| Item Number | Part Number | Description | Quantity |
|--------------------|-------------|---|----------|
| 1p, 2b, 3b, 4, 5d | 506533 | Back Upholstery Assembly Tension Adjustable with Standard Cover 17W 18.5-19.5, Half Folding Back | 1 |
| 1q, 2a, 3a, 4, 5e | 506534 | Back Upholstery Assembly Tension Adjustable with Standard Cover 18W 12.5-13.5H, Half Folding Back Straight | 1 |
| 1r, 2a, 3a, 4, 5e | 506535 | Back Upholstery Assembly Tension Adjustable with Standard Cover 18W 14.5-15.5H, Half Folding Back Straight | 1 |
| 1r, 2a, 3a, 4, 5e | 506536 | Back Upholstery Assembly Tension Adjustable with Standard Cover 18W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 1r, 2a, 3a, 4, 5e | 506537 | Back Upholstery Assembly Tension Adjustable with Standard Cover 18W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 1s, 2b, 3b, 4, 5e | 506538 | Back Upholstery Assembly Tension Adjustable with Standard Cover 18W 16.5-17.5, Half Folding Back | 1 |
| 1t, 2b, 3b, 4, 5e | 506539 | Back Upholstery Assembly Tension Adjustable with Standard Cover 18W 18.5-19.5, Half Folding Back | 1 |
| 1u, 2a, 3a, 4, 5f | 506540 | Back Upholstery Assembly Tension Adjustable with Standard Cover 19W 12.5-13.5H, Half Folding Back Straight | 1 |
| 1v, 2a, 3a, 4, 5f | 506541 | Back Upholstery Assembly Tension Adjustable with Standard Cover 19W 14.5-15.5H, Half Folding Back Straight | 1 |
| 1v, 2a, 3a, 4, 5f | 506542 | Back Upholstery Assembly Tension Adjustable with Standard Cover 19W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 1v, 2a, 3a, 4, 5f | 506543 | Back Upholstery Assembly Tension Adjustable with Standard Cover 19W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 1w, 2b, 3b, 4, 5f | 506544 | Back Upholstery Assembly Tension Adjustable with Standard Cover 19W 16.5-17.5, Half Folding Back | 1 |
| 1x, 2b, 3b, 4, 5f | 506545 | Back Upholstery Assembly Tension Adjustable with Standard Cover 19W 18.5-19.5, Half Folding Back | 1 |
| 1y, 2a, 3a, 4, 5g | 506546 | Back Upholstery Assembly Tension Adjustable with Standard Cover 20W 12.5-13.5H, Half Folding Back Straight | 1 |
| 1z, 2a, 3a, 4, 5g | 506547 | Back Upholstery Assembly Tension Adjustable with Standard Cover 20W 14.5-15.5H, Half Folding Back Straight | 1 |
| 1z, 2a, 3a, 4, 5g | 506548 | Back Upholstery Assembly Tension Adjustable with Standard Cover 20W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 1z, 2a, 3a, 4, 5g | 506549 | Back Upholstery Assembly Tension Adjustable with Standard Cover 20W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 1aa, 2b, 3b, 4, 5g | 506550 | Back Upholstery Assembly Tension Adjustable with Standard Cover 20W 16.5-17.5, Half Folding Back | 1 |
| 1ab, 2b, 3b, 4, 5g | 506551 | Back Upholstery Assembly Tension Adjustable with Standard Cover 20W 18.5-19.5, Half Folding Back | 1 |

| Item Number | Part Number | Description | Quantity |
|--------------------|-------------|---|----------|
| 1ac, 2a, 3a, 4, 5h | 506552 | Back Upholstery Assembly Tension Adjustable with Standard Cover 21W 12.5-13.5H, Half Folding Back Straight | 1 |
| 1ad, 2a, 3a, 4, 5h | 506553 | Back Upholstery Assembly Tension Adjustable with Standard Cover 21W 14.5-15.5H, Half Folding Back Straight | 1 |
| 1ad, 2a, 3a, 4, 5h | 506554 | Back Upholstery Assembly Tension Adjustable with Standard Cover 21W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 1ad, 2a, 3a, 4, 5h | 506555 | Back Upholstery Assembly Tension Adjustable with Standard Cover 21W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 1ae, 2b, 3b, 4, 5h | 506556 | Back Upholstery Assembly Tension Adjustable with Standard Cover 21W 16.5-17.5, Half Folding Back | 1 |
| 1af, 2b, 3b, 4, 5h | 506557 | Back Upholstery Assembly Tension Adjustable with Standard Cover 21W 18.5-19.5, Half Folding Back | 1 |
| 1ag, 2a, 3a, 4, 5h | 506558 | Back Upholstery Assembly Tension Adjustable with Standard Cover 22W 14.5-15.5H, Half Folding Back Straight | 1 |
| 1ag, 2a, 3a, 4, 5h | 506559 | Back Upholstery Assembly Tension Adjustable with Standard Cover 22W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 1ag, 2a, 3a, 4, 5h | 506560 | Back Upholstery Assembly Tension Adjustable with Standard Cover 22W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 1ah, 2b, 3b, 4, 5h | 506561 | Back Upholstery Assembly Tension Adjustable with Standard Cover 22W 16.5-17.5, Half Folding Back | 1 |
| 1ai, 2b, 3b, 4, 5h | 506562 | Back Upholstery Assembly Tension Adjustable with Standard Cover 22W 18.5-19.5, Half Folding Back | 1 |
| 2a, 3a, 4, 5a, 7a | 107950 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 14W 12.5-13.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5a, 7b | 107951 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 14W 14.5-15.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5a, 7b | 107952 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 14W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 2a, 3a, 4, 5a, 7b | 107953 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 14W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 2b, 3b, 4, 5a, 7c | 107954 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 14W 16.5-17.5, Half Folding Back | 1 |
| 2b, 3b, 4, 5a, 7d | 107955 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 14W 18.5-19.5, Half Folding Back | 1 |
| 2a, 3a, 4, 5b, 7e | 107956 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 15W 12.5-13.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5b, 7f | 107957 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 15W 14.5-15.5H, Half Folding Back Straight | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------------|-------------|---|----------|
| 2a, 3a, 4, 5b, 7f | 107958 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 15W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 2a, 3a, 4, 5b, 7f | 107959 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 15W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 2b, 3b, 4, 5b, 7g | 107960 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 15W 16.5-17.5, Half Folding Back | 1 |
| 2b, 3b, 4, 5b, 7h | 107961 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 15W 18.5-19.5, Half Folding Back | 1 |
| 2a, 3a, 4, 5c, 7i | 107962 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 16W 12.5-13.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5c, 7j | 107963 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 16W 14.5-15.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5c, 7j | 107964 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 16W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 2a, 3a, 4, 5c, 7j | 107965 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 16W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 2b, 3b, 4, 5c, 7k | 107966 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 16W 16.5-17.5, Half Folding Back | 1 |
| 2b, 3b, 4, 5c, 7l | 107967 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 16W 18.5-19.5, Half Folding Back | 1 |
| 2a, 3a, 4, 5d, 7m | 107968 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 17W 12.5-13.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5d, 7n | 107969 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 17W 14.5-15.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5d, 7n | 107970 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 17W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 2a, 3a, 4, 5d, 7n | 107971 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 17W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 2b, 3b, 4, 5d, 7o | 107972 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 17W 16.5-17.5, Half Folding Back | 1 |
| 2b, 3b, 4, 5d, 7p | 107973 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 17W 18.5-19.5, Half Folding Back | 1 |
| 2a, 3a, 4, 5e, 7q | 107974 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 18W 12.5-13.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5e, 7r | 107975 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 18W 14.5-15.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5e, 7r | 107976 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 18W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |

| Item Number | Part Number | Description | Quantity |
|--------------------|-------------|---|----------|
| 2a, 3a, 4, 5e, 7r | 107977 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 18W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 2b, 3b, 4, 5e, 7s | 107978 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 18W 16.5-17.5, Half Folding Back | 1 |
| 2b, 3b, 4, 5e, 7t | 107979 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 18W 18.5-19.5, Half Folding Back | 1 |
| 2a, 3a, 4, 5f, 7u | 107980 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 19W 12.5-13.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5f, 7v | 107981 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 19W 14.5-15.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5f, 7v | 107982 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 19W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 2a, 3a, 4, 5f, 7v | 107983 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 19W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 2b, 3b, 4, 5f, 7w | 107984 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 19W 16.5-17.5, Half Folding Back | 1 |
| 2b, 3b, 4, 5f, 7x | 107985 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 19W 18.5-19.5, Half Folding Back | 1 |
| 2a, 3a, 4, 5g, 7y | 107986 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 20W 12.5-13.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5g, 7z | 107987 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 20W 14.5-15.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5g, 7z | 107988 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 20W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 2a, 3a, 4, 5g, 7z | 107989 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 20W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 2b, 3b, 4, 5g, 7aa | 107990 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 20W 16.5-17.5, Half Folding Back | 1 |
| 2b, 3b, 4, 5g, 7ab | 107991 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 20W 18.5-19.5, Half Folding Back | 1 |
| 2a, 3a, 4, 5h, 7ac | 107992 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 21W 12.5-13.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5h, 7ad | 107993 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 21W 14.5-15.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5h, 7ad | 107994 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 21W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 2a, 3a, 4, 5h, 7ad | 107995 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 21W 15.5H, Half Folding Back 8 Degree Bend | 1 |

| Item Number | Part Number | Description | Quantity |
|--------------------|-------------|---|----------|
| 2b, 3b, 4, 5h, 7ae | 107996 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 21W 16.5-17.5, Half Folding Back | 1 |
| 2b, 3b, 4, 5h, 7af | 107997 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 21W 18.5-19.5, Half Folding Back | 1 |
| 2a, 3a, 4, 5h, 7ag | 107998 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 22W 14.5-15.5H, Half Folding Back Straight | 1 |
| 2a, 3a, 4, 5h, 7ag | 107999 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 22W 13.5-14.5H, Half Folding Back 8 Degree Bend | 1 |
| 2a, 3a, 4, 5h, 7ag | 108000 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 22W 15.5H, Half Folding Back 8 Degree Bend | 1 |
| 2b, 3b, 4, 5h, 7ah | 108001 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 22W 16.5-17.5, Half Folding Back | 1 |
| 2b, 3b, 4, 5h, 7ai | 108002 | Back Upholstery Assembly Tension Adjustable with Air Mesh Cover 22W 18.5-19.5, Half Folding Back | 1 |
| 1a | 002498 | Backrest Cover 13" Rigid / 14" Cat x 11-14" | 1 |
| 1b | 002499 | Backrest Cover 13" Rigid / 14" Cat x 13-16" | 1 |
| 1c | 002500 | Backrest Cover 13" Rigid - 14" CAT/Liberty x 15-18" | 1 |
| 1d | 003043 | Backrest Cover 13" Rigid - 14" CAT/Liberty x 17-20" | 1 |
| 1e | 002502 | Backrest Cover 14" Rigid / 15" Cat x 11-14" | 1 |
| 1f | 002503 | Backrest Cover 14" Rigid / 15" Cat x 13-16" | 1 |
| 1g | 002504 | Backrest Cover 14" Rigid - 15" CAT/Liberty x 15-18" | 1 |
| 1h | 003044 | Backrest Cover 14" Rigid - 15" CAT/Liberty x 17-20" | 1 |
| 1i | 002506 | Backrest Cover 15" Rigid / 16" Cat x 11-14" | 1 |
| 1j | 002507 | Backrest Cover 15" Rigid / 16" Cat x 13-16" | 1 |
| 1k | 002508 | Backrest Cover 15" Rigid - 16" CAT/Liberty x 15-18" | 1 |
| 1l | 003045 | Backrest Cover 15" Rigid - 16" CAT/Liberty x 17-20" | 1 |
| 1m | 002510 | Backrest Cover 16" Rigid / 17" Cat x 11-14" | 1 |
| 1n | 002511 | Backrest Cover 16" Rigid / 17" Cat x 13-16" | 1 |
| 1o | 002512 | Backrest Cover 16" Rigid - 17" CAT/Liberty x 15-18" | 1 |
| 1p | 003046 | Backrest Cover 16" Rigid - 17" CAT/Liberty x 17-20" | 1 |
| 1q | 003049 | Backrest Cover 17" Rigid / 18" Cat x 11-14" | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 1r | 003050 | Backrest Cover 17" Rigid / 18" Cat x 13-16" | 1 |
| 1s | 003051 | Backrest Cover 17" Rigid - 18" CAT/Liberty x 15-18" | 1 |
| 1t | 003052 | Backrest Cover 17" Rigid - 18" CAT/Liberty x 17-20" | 1 |
| 1u | 003055 | Backrest Cover 18" Rigid / 19" Cat x 11-14" | 1 |
| 1v | 003056 | Backrest Cover 18" Rigid / 19" Cat x 13-16" | 1 |
| 1w | 003057 | Backrest Cover 18" Rigid - 19" CAT/Liberty x 15-18" | 1 |
| 1x | 003058 | Backrest Cover 18" Rigid - 19" CAT/Liberty x 17-20" | 1 |
| 1y | 003061 | Backrest Cover 19" Rigid / 20" Cat x 11-14" | 1 |
| 1z | 003062 | Backrest Cover 19" Rigid / 20" Cat x 13-16" | 1 |
| 1aa | 003063 | Backrest Cover 19" Rigid - 20" CAT/Liberty x 15-18" | 1 |
| 1ab | 003064 | Backrest Cover 19" Rigid - 20" CAT/Liberty x 17-20" | 1 |
| 1ac | 003067 | Backrest Cover 20" Rigid / 21" Cat x 11-14" | 1 |
| 1ad | 003068 | Backrest Cover 20" Rigid / 21" Cat x 13-16" | 1 |
| 1ae | 003069 | Backrest Cover 20" Rigid - 21" CAT/Liberty x 15-18" | 1 |
| 1af | 003070 | Backrest Cover 20" Rigid - 21" CAT/Liberty x 17-20" | 1 |
| 1ag | 003072 | Backrest Cover 21" Rigid / 22" Cat x 13-16" | 1 |
| 1ah | 003073 | Backrest Cover 21" Rigid - 22" CAT/Liberty x 15-18" | 1 |
| 1ai | 003074 | Backrest Cover 21" Rigid - 22" CAT/Liberty x 17-20" | 1 |
| 2a | 003095 | Tension Adjustable Upholstery Wing 5.5" Height Left | 1 |
| 2b | 003097 | Tension Adjustable Upholstery Wing 9.5" Height Left | 1 |
| 3a | 003094 | Tension Adjustable Upholstery Wing 5.5" Height Right | 1 |
| 3b | 003096 | Tension Adjustable Upholstery Wing 9.5" Height Right | 1 |
| 4 | 100973 | Back Upholstery D-Ring | 4 |
| 5a | 003487 | Back Upholstery Strap 15.3" <i>Used on 14" chair widths</i> | 2 |
| 5b | 002532 | Back Upholstery Strap 17.3" <i>Used on 15" chair widths</i> | 2 |
| 5c | 002533 | Back Upholstery Strap 19.3" | 2 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| | | <i>Used on 16" chair widths</i> | |
| 5d | 002534 | Back Upholstery Strap 21.3" <i>Used on 17" chair widths</i> | 2 |
| 5e | 002535 | Back Upholstery Strap 23.3" <i>Used on 18" chair widths</i> | 2 |
| 5f | 002536 | Back Upholstery Strap 25.3" <i>Used on 19" chair widths</i> | 2 |
| 5g | 002537 | Back Upholstery Strap 27.3" <i>Used on 20" chair widths</i> | 2 |
| 5h | 002538 | Back Upholstery Strap 29.3" <i>Used on 21" chair widths</i> | 3 |
| 5i | 002539 | Back Upholstery Strap 31.3" <i>Used on 22" chair widths</i> | 3 |
| 6 | 100734 | M5 x 16 PHMS ZC | 2 |
| 7a | 003897 | Air Mesh Backrest Cover 14C 13R Width x 11-14" Height | 1 |
| 7b | 003898 | Air Mesh Backrest Cover 14C 13R Width x 13-16" Height | 1 |
| 7c | 003899 | Air Mesh Backrest Cover 14C 13R Width x 15-18" Height | 1 |
| 7d | 003913 | Air Mesh Backrest Cover 14C 13R Width x 17-20" Height | 1 |
| 7e | 003901 | Air Mesh Backrest Cover 15C 14R Width x 11-14" Height | 1 |
| 7f | 003902 | Air Mesh Backrest Cover 15C 14R Width x 13-16" Height | 1 |
| 7g | 003903 | Air Mesh Backrest Cover 15C 14R Width x 15-18" Height | 1 |
| 7h | 003914 | Air Mesh Backrest Cover 15C 14R Width x 17- 20" Height | 1 |
| 7i | 003905 | Air Mesh Backrest Cover 16C 15R Width x 11-14" Height | 1 |
| 7j | 003906 | Air Mesh Backrest Cover 16C 15R Width x 13-16" Height | 1 |
| 7k | 003907 | Air Mesh Backrest Cover 16C 15R Width x 15-18" Height | 1 |
| 7l | 003915 | Air Mesh Backrest Cover 16C 15R Width x 17-20" Height | 1 |
| 7m | 003910 | Air Mesh Backrest Cover 17C 16R Width x 11-14" Height | 1 |
| 7n | 003911 | Air Mesh Backrest Cover 17C 16R Width x 13-16" Height | 1 |
| 7o | 003912 | Air Mesh Backrest Cover 17C 16R Width x 15-18" Height | 1 |
| 7p | 003916 | Air Mesh Backrest Cover 17C 16R Width x 17-20" Height | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| 7q | 003919 | Air Mesh Backrest Cover 18C 17R Width x 11-14" Height | 1 |
| 7r | 003920 | Air Mesh Backrest Cover 18C 17R Width x 13-16" Height | 1 |
| 7s | 003921 | Air Mesh Backrest Cover 18C 17R Width x 15-18" Height | 1 |
| 7t | 003922 | Air Mesh Backrest Cover 18C 17R Width x 17-20" Height | 1 |
| 7u | 003925 | Air Mesh Backrest Cover 19C 18R Width x 11-14" Height | 1 |
| 7v | 003926 | Air Mesh Backrest Cover 19C 18R Width x 13-16" Height | 1 |
| 7w | 003927 | Air Mesh Backrest Cover 19C 18R Width x 15-18" Height | 1 |
| 7x | 003928 | Air Mesh Backrest Cover 19C 18R Width x 17-20" Height | 1 |
| 7y | 003931 | Air Mesh Backrest Cover 20C 19R Width x 11-14" Height | 1 |
| 7z | 003932 | Air Mesh Backrest Cover 20C 19R Width x 13-16" Height | 1 |
| 7aa | 003933 | Air Mesh Backrest Cover 20C 19R Width x 15-18" Height | 1 |
| 7ab | 003934 | Air Mesh Backrest Cover 20C 19R Width x 17-20" Height | 1 |
| 7ac | 003937 | Air Mesh Backrest Cover 21C 20R Width x 11-14" Height | 1 |
| 7ad | 003938 | Air Mesh Backrest Cover 21C 20R Width x 13-16" Height | 1 |
| 7ae | 003939 | Air Mesh Backrest Cover 21C 20R Width x 15-18" Height | 1 |
| 7af | 003940 | Air Mesh Backrest Cover 21C 20R Width x 17-20" Height | 1 |
| 7ag | 003942 | Air Mesh Backrest Cover 22C 21R Width x 13-16" Height | 1 |
| 7ah | 003943 | Air Mesh Backrest Cover 22C 21R Width x 15-18" Height | 1 |
| 7ai | 003944 | Air Mesh Backrest Cover 22C 21R Width x 17-20" Height | 1 |