

AQUATEC OCEAN / OCEAN XL



OCEAN basic model

- ... Stable, stainless steel frame
- ... Individually adjustable seat height (simple to adjust, no tools required)
- ... Washable backrest with quick-release buckles – can be adjusted to the user's needs
- ... Comfortable due to the ergonomically designed seat
- ... Removable seat – easy to clean
- ... Height-adjustable, swivel-mounted and detachable footrests for easy transfer and adjustment
- ... Structured fold-away footrests with adjustable heel loops provide a safe hold
- ... Four braked castors for reliable stopping and fixation
- ... Easily manoeuvrable due to smooth castors
- ... Ergonomic fold-away armrests
- ... Capable for pool-use
- ... Knock-down concept allows easy transport and use on travels
- ... Space-saving packing size
- ... Comprehensive accessories for specific requirements
- ... Smooth surfaces and flowing contours for easy cleaning
- ... Additionally available as XL version with 150 kg loading capacity and more space between the armrests
- ... Can also be retrofitted to the XL version



OCEAN with 24" self-propelling wheels




design award
winner
2004




Armrest width extension OCEAN XL

AQUATEC OCEAN / OCEAN XL


AQUATEC OCEAN 24"-wheels




AQUATEC OCEAN XL 24"-wheels



AQUATEC OCEAN



AQUATEC OCEAN XL



Order numbers

AQUATEC OCEAN

AQUATEC OCEAN
24"-WHEELS

AQUATEC OCEAN XL

AQUATEC OCEAN XL
24"-WHEELS

7.10.900

7.12.900

7.13.900

7.14.900

Dimensions

Height to armrests (A)	680-810 mm	680-810 mm	680-810 mm	680-810 mm
Seat height (B)	477-600 mm	477-600 mm	477-600 mm	477-600 mm
Seat depth (C)	450 mm	450 mm	450 mm	450 mm
Overall width (D)	560 mm	675 mm	640 mm	675 mm
Overall length (E)	934 mm	1060 mm	934 mm	1060 mm
Length without footrests (F)	730 mm	856 mm	730 mm	856 mm
Footrest depth (G)	154 mm	154 mm	154 mm	154 mm
Overall height (H)	964-1094 mm	964-1094 mm	964-1094 mm	964-1094 mm
Height seat to lockable footrest (I)	280-530 mm	280-530 mm	280-530 mm	280-530 mm
Height ground to lower seat edge:				
with holder for sanitary pan (J)	390-520 mm	390-520 mm	390-520 mm	390-520 mm
without holder for sanitary pan (J)	410-540 mm	410-540 mm	410-540 mm	410-540 mm
Width between armrests (K)	459 mm	459 mm	539 mm	539 mm
Width between 24"-wheels	-	530 mm	-	530 mm
Height adjustment steps	25 mm	25 mm	25 mm	25 mm
Wheel diameter	127 mm	127 mm	127 mm	127 mm
Wheel diameter of 24" self-propelling wheels	-	610 mm	-	610 mm

Other

Weight	14.5 kg	19.5 kg	14.5 kg	19.5 kg
Load capacity	130 kg	130 kg	150 kg	150 kg
Colour	white	white	white	white

Heel loops are included in delivery of **AQUATEC OCEAN** and **OCEAN XL**.



Soft seat
with hygiene recess.
16341 for OCEAN / OCEAN XL



Insert
for soft seat (16341).
16343 for OCEAN / OCEAN XL



Soft seat
with small oval recess.
16342 for OCEAN / OCEAN XL



Splash guard
easily fitted. For soft seat with hygiene recess (16341).
16511 for OCEAN / OCEAN XL

NEW!



Safety bar
for additional security.
10167-10 for OCEAN
length: 54 cm



Armrest lock
locks the armrest in lower position.
16597 for OCEAN / OCEAN XL (1 piece)

NEW!



Armrest width extension
widens the distance between the armrests by 8 cm (4 cm on each side). Using the armrest width extension the load capacity increases to 150 kg.
16515 Upgrade kit OCEAN XL



24"-self-propelling wheels
set including brakes and quick release axle.
16379 for OCEAN / OCEAN XL



Amputation support
suits left or right side. Adjustable in height and angle. Upholstered pad.
16616 for OCEAN / OCEAN XL

NEW!

ACCESSORIES



Calf support

suits left or right side. Adjustable in height and angle. Upholstered pad.

16615 for OCEAN / OCEAN XL

NEW!



Sanitary pan

with handle and lid.

10230 for OCEAN / OCEAN XL
length: 46 cm

16335 Sanitary pan holder
for OCEAN / OCEAN XL (not ill.)



Transport bag

with castors. Suitable for all OCEAN-models.

16521 for OCEAN / OCEAN XL

NEW!