


Liberty O2 Holder



Rigid Frame Hardware

Item Number	Part Number	Description	Quantity
1-9	506200	O2 Holder Assembly, Liberty	1
4, 5, 10, 11	118062	O2 Holder Rigid Hardware Kit, Liberty <i>Used with Rigid Liberty Frame Only</i>	1
1	001416	O2 Holder	1
2	003338	Bracket, O2 Holder, Lower	1
3	003339	Bracket, O2 Holder, Upper	1
4	100658	M6 Nylock Nut Blk Zn	3
5	100746	M6 Flat Washer Blk Zn	5

Item Number	Part Number	Description	Quantity
6	101594	M6 x 14 CLS10.9 DIN7991 FSHCS BLKZC w/Patch	1
7	102327	M6 x 16 HHCS Blk Zn	2
8	103026	M5 x 8 HHCS, Blk Zn, w/Patch	1
9	103027	M6 x 80 HHCS, Partial Thread, Blk Zn	1
10	116440	Spacer, .59 OD x .251 ID x .250 L	1
11	116854	HHCS, M6x1, 0x30, CLS10.9, PT, BLZ, DIN931	1