


Arc Seat Pan



Pan (Items 1 and 2) Part Numbers					
	Frame Depth	Frame Width			
		Small (10-12")	Medium (12-14")	Large (14-16")	Extra Large (16-18")
Item 1	Extra Short (12-15")	113269	113270	113271	113272
	Short (14-17")	113273	113274	113275	113276
	Medium (16-19")	113277	113278	113279	113280
	Long (18-20")	113281	113282	113283	113284
Item 2	All	111430	111431	111432	111433

Item Number	Part Number	Description	Quantity
1a, 2a, 3-8	113253	Seat Pan Assembly, Extra Short Depth, 10-12" Wide, Arc	1
1b, 2b, 3-8	113254	Seat Pan Assembly, Extra Short Depth, 12-14" Wide, Arc	1
1c, 2c, 3-8	113255	Seat Pan Assembly, Extra Short Depth, 14-16" Wide, Arc	1
1d, 2d, 3-8	113256	Seat Pan Assembly, Extra Short Depth, 16-18" Wide, Arc	1

Item Number	Part Number	Description	Quantity
1e, 2a, 3-8	113257	Seat Pan Assembly, Short Depth, 10-12" Wide, Arc	1
1f, 2b, 3-8	113258	Seat Pan Assembly, Short Depth, 12-14" Wide, Arc	1
1g, 2c, 3-8	113259	Seat Pan Assembly, Short Depth, 14-16" Wide, Arc	1
1h, 2d, 3-8	113260	Seat Pan Assembly, Short Depth, 16-18" Wide, Arc	1
1i, 2a, 3-8	113261	Seat Pan Assembly, Medium Depth, 10-12" Wide, Arc	1
1j, 2b, 3-8	113262	Seat Pan Assembly, Medium Depth, 12-14" Wide, Arc	1
1k, 2c, 3-8	113263	Seat Pan Assembly, Medium Depth, 14-16" Wide, Arc	1
1l, 2d, 3-8	113264	Seat Pan Assembly, Medium Depth, 16-18" Wide, Arc	1
1m, 2a, 3-8	113265	Seat Pan Assembly, Long Depth, 10-12" Wide, Arc	1
1n, 2b, 3-8	113266	Seat Pan Assembly, Long Depth, 12-14" Wide, Arc	1
1o, 2c, 3-8	113267	Seat Pan Assembly, Long Depth, 14-16" Wide, Arc	1
1p, 2d, 3-8	113268	Seat Pan Assembly, Long Depth, 16-18" Wide, Arc	1
1a, 4	113269	Seat Pan with AB Tape, Front, Extra Short Depth, 10-12" Wide, Arc	1
1b, 4	113270	Seat Pan with AB Tape, Front, Extra Short Depth, 12-14" Wide, Arc	1
1c, 4	113271	Seat Pan with AB Tape, Front, Extra Short Depth, 14-16" Wide, Arc	1
1d, 4	113272	Seat Pan with AB Tape, Front, Extra Short Depth, 16-18" Wide, Arc	1
1e, 4	113273	Seat Pan with AB Tape, Front, Short Depth, 10-12" Wide, Arc	1
1f, 4	113274	Seat Pan with AB Tape, Front, Short Depth, 12-14" Wide, Arc	1
1g, 4	113275	Seat Pan with AB Tape, Front, Short Depth, 14-16" Wide, Arc	1
1h, 4	113276	Seat Pan with AB Tape, Front, Short Depth, 16-18" Wide, Arc	1
1i, 4	113277	Seat Pan with AB Tape, Front, Medium Depth, 10-12" Wide, Arc	1
1j, 4	113278	Seat Pan with AB Tape, Front, Medium Depth, 12-14" Wide, Arc	1
1k, 4	113279	Seat Pan with AB Tape, Front, Medium Depth, 14-16" Wide, Arc	1
1l, 4	113280	Seat Pan with AB Tape, Front, Medium Depth, 16-18" Wide, Arc	1
1m, 4	113281	Seat Pan with AB Tape, Front, Long Depth, 10-12" Wide, Arc	1
1n, 4	113282	Seat Pan with AB Tape, Front, Long Depth, 12-14" Wide, Arc	1
1o, 4	113283	Seat Pan with AB Tape, Front, Long Depth, 14-16" Wide, Arc	1

Item Number	Part Number	Description	Quantity
1p, 4	113284	Seat Pan with AB Tape, Front, Long Depth, 16-18" Wide, Arc	1
2a	111430	Seat Pan, Rear, 10-12 Wide	1
2b	111431	Seat Pan, Rear, 12-14 Wide	1
2c	111432	Seat Pan, Rear, 14-16 Wide	1
2d	111433	Seat Pan, Rear, 16-18 Wide	1
3	111407	Seat Pan Mount, Snap-on	4
4	001206	AB Tape 3x3 Die Cut Tilt	2
5	102856	Seat Pan Stop Collar	2
6	102634	M6 X 1.0 X 16 BHCS BLK ZC	8
7	100746	M6 Flat Washer Blk Zn	4
8	100658	M6 Nylock Nut Blk Zn	8