


Clik Camber Tube - Growth



Black Camber Tube kits based on seat width. Color Anodized Camber Tube kits used on all chair widths. <p>

New camber tube required when changing the degree of camber by more than 2 degrees. <p>

If Wheelchair Drive Ready Option is present, replacement parts need to be processed through Alber USA.

Item Number	Part Number	Description	Quantity	Retail
1a	503238	Camber Tube Rigid 8" Kit Black	1	\$172.94
1b	503239	Camber Tube Rigid 9" Kit Black	1	\$172.94
1c	500720	Camber Tube Rigid 10" Kit Black	1	\$176.47
1d	500721	Camber Tube Rigid 11" Kit Black	1	\$176.47
1e	003369	One Arm Drive Camber Tube 11" Kit, Clik <i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i>	1	\$176.47
1f	500722	Camber Tube Rigid 12" Kit Black	1	\$176.47
1g	003370	One Arm Drive Camber Tube 12" Kit, Clik <i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i>	1	\$176.47

Item Number	Part Number	Description	Quantity	Retail
1h	500723	Camber Tube Rigid 13" Kit Black	1	\$176.47
1i	003371	One Arm Drive Camber Tube 13" Kit, Klik <i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i>	1	\$176.47
1j	500724	Camber Tube Rigid 14" Kit Black	1	\$176.47
1k	003372	One Arm Drive Camber Tube 14" Kit, Klik <i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i>	1	\$176.47
1l	500725	Camber Tube Rigid 15" Kit Black	1	\$176.47
1m	003373	One Arm Drive Camber Tube 15" Kit, Klik <i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i>	1	\$176.47
1n	500726	Camber Tube Rigid 16" Kit Black	1	\$176.47
1o	003374	One Arm Drive Camber Tube 16" Kit, Klik <i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i>	1	\$176.47
1p	500727	Camber Tube Rigid 17" Kit Black	1	\$176.47
1q	003708	One Arm Drive Camber Tube 17" Kit, Klik <i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i>	1	\$176.47
1r	500728	Camber Tube Rigid 18" Kit Black	1	\$176.47
1s	003709	One Arm Drive Camber Tube 18" Kit, Klik <i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i>	1	\$176.47
1t	501749	Camber Tube Rigid Blue Anodized	1	\$190.93
1u	501750	Camber Tube Rigid Bright Green Anodized	1	\$190.93
1v	501751	Camber Tube Rigid Orange Anodized	1	\$190.93
1w	501752	Camber Tube Rigid Pink Anodized	1	\$190.93
1x	501753	Camber Tube Rigid Pewter Anodized	1	\$190.93
1y	501754	Camber Tube Rigid Purple Anodized	1	\$190.93
1z	501755	Camber Tube Rigid Red Anodized	1	\$190.93
1aa	501756	Camber Tube Rigid Gold Anodized	1	\$190.93
2	100667	M6 x 12 BTN SHCS BLKZC w/Patch	4	\$1.58
3	000066	M6 x 27 Threaded Barrel	2	\$4.46

Color Anodizing Samples

